

Spirit and Heine Sphygmomanometers

Spirit Aneroid Portable Sphygmomanometer

- 2 Tube Style Sphygmomanometer with separate bulb and aneroid dial.
- Latex free
- Accurately calibrated gauge with measurement from 0 - 300mmHG.
- Includes inflation system and zipper vinyl bag with adult cuff bladder



Product Code	Description
SP110	SPHYG - SPIRIT ANEROID PORTABLE (EACH)

Spirit Aneroid Sphygmomanometer Wall and Mobile Units

- Professional Model
- Latex free and specifically designed for use in hospitals and clinical environment.
- No pin stop calibrated scale measures from 0 - 300mmHg
- Enlarged blue coloured dial with white numerical measures for easy reading
- Wall mounted unit includes: Inflation system with adult size cuff and bladder, extendable rubber coil tubing, wall mounting hardware and cuff storage basket (not pictured).
- Mobile unit includes: Mobile adjustable stand with 42cm diameter base with 5 wheels, adjustable height from 102 cm - 135cm, inflation storage basket, sphygmomanometer gauge and inflation system including adult cuff / bladder, coil tubing and bulb.



Product Code	Description
SPCK141A	SPHYG - SPIRIT ANEROID WALL MOUNTED W/ BASKET (SET)
SPCK142	SPHYG - SPIRIT ANEROID ON MOBILE STAND (SET)

SPCK142

Heine Gamma G5 Sphygmomanometer w/Adult Cuff

- Integrated Dial and Bulb to allow for single hand use with 5 year warranty*
- Single tube latex free adult cuff and bladder with velcro fastener. Range of sizes available as option from Infant to Thigh.
- Shockproof pressure distribution system and Pressure proof manometer to ensure durability + reliability
- Large 56mm diameter easy to read scale and Adjustable spoon with bulb allowing use in left or right hand



HM00009230

Product Code	Description
HM00009230	SPHYG - HEINE ANEROID GAMMA G5 W/ ADULT CUFF (EACH)

Heine Gamma XXL Sphygmomanometer

- Top Quality Professional Model designed for hospital use with 10 year warranty*
- Wall, desk and adjustable mobile stand mounting options available.
- Durable polycarbonate construction with oversized 130mm diameter blue on white scale.
- High quality pressure proof manometer with precision air release valve, precisely calibrated + durable aneroid dial.
- 130mm oversize dial that can be pivoted for clear easy reading for ease of use by medical personnel.
- Supplied with Adult size cuff and bladder (latex free). Range sizes available from Infant to Thigh.

* against defects in design, materials and workmanship for a period of 10 years.



HM00009326

Product Code	Description
HM0009323	SPHYG - HEINE ANEROID GAMMA XXL W/ WALL MT (EACH)
HM0009326	SPHYG - HEINE ANEROID GAMMA XXL W/ WHEELED STAND (EACH)
HM0009322	SPHYG - HEINE ANEROID GAMMA XXL W/ DESK MT (EACH)



W: www.dbhealth.com.au
 E: customerservice@denyer.com.au
 P: (02) 9638 2100
 ABN: 37 000 459 174