

Customisable Options for Emery Trolleys

Upgrade Options	Product Code
Shelf Rails Bottom	HQAO-SRB
Braking Caster	HQAO-CAST
Bucket holder LHS Mounted	HQAO-BUCKET HOLDER
Chart Holder***	HQAO-CH
Drawer Dividers / Trays	HQAO-PDD001
Drawer Dividers / Trays	HQAO-PDD001-1
Drawer Dividers / Trays	HQAO-PDD002
Drawer Dividers / Trays	HQAO-PDD005
Glove Box Holder*** (3 x Ansell Boxes)	HQAO-GLOVE BH
IV Pole Mount on corner	HQAO-IV-TROLLEY
IV Pole Flat Top Trolley Mount on corner	HQAO-IV TROLLEY FT
Keyed Lock - RHS Mounted**	HQAO-LC****
Magnifying Lamp	HQAO-MAG LIGHT TROLLEY
Mesh Drawers Bottom	HQAO-MESH DRAWER
Mesh Base SS Catheter Holder***	HQAO-CH-150-800
Push Handle	HQAO-PH
Drop down Shelf RHS Mount	HQAO-DDS
Size C Gas Bottle Holder	HQAO-GAS BOTTLE HOLDER
Size C Gas Bottle Holder Portable***	HQAO-GBHP

*HQAO-PDD001

Drawer Dividers to fit Standard 500 x 500 Unit
3mm Clear Polycarbonate to fit One drawer
2 of 400 x 105 high dividers with 5 slots
2 of 365 x 105 high dividers with 5 slots
Creates 9 Compartments

*HQAO-PDD001-1

Drawer Dividers to fit Standard 500 x 500 Unit
3mm Clear Polycarbonate to fit One drawer
1 of 400 x 105 high dividers with 5 slots
1 of 365 x 105 high dividers with 5 slots
Creates 4 Compartments

*HQAO-PDD002

3mm Clear Drawer Dividers – T Divider
Fit drawers on Standard 500 x 500 unit
1 of 390 x 115 high T dividers
Run front to back only
Creates 2 Compartments

*HQAO-PDD005

Drawer Dividers to fit Standard 750 x 500 Unit
3mm Clear Polycarbonate to fit One drawer
2 of 665 x 105 high dividers with 5 slots
4 of 365 x 105 high dividers with 5 slots
Creates 15 Compartments



HQAO-PDD001



*HQAO-PDD002

All accessories are added during manufacture

* Please select required Drawer Dividers

** Fitted to Right hand side unless another accessory is on the right hand side

*** For Trolleys with Side Rails only

**** Locking System can be for all drawers or a combination of drawers. Please state how many drawers require to have the locking system when placing your order.



HQAO-SRB - Shelf Rails Bottom



HQAO-CAST - Braking



HQAO-BUCKET HOLDER - Bucket holder LHS Mounted



HQAO-CH - Chart Holder***



HQAO-GLOVE BH - Glove Box Holder*** (3 x Ansell Boxes)



HQAO-IV-TROLLEY - IV Pole Mount on corner

HQAO-IV TROLLEY FT - IV Pole Flat Top Trolley Mount on corner



HQAO-LC**** - Keyed Lock - RHS Mounted**



HQAO-MAG LIGHT TROLLEY - Magnifying Lamp



HQAO-MESH DRAWER -
Mesh Drawers Bottom



HQAO-CH-150-800 -
Mesh Base SS Cathe-



HQAO-PH - Push Handle



HQAO-DDS - Drop down Shelf
RHS Mount



HQAO-GAS BOTTLE HOLDER
- Size C Gas Bottle Holder



HQAO-GBHP - Size C Gas
Bottle Holder Portable***

Emery Powdercoat Colours

PEARL WHITE



BIRCH WHITE



LEMON YELLOW



PINK BLUSH



FLAME RED



VIPER GREEN



POMMEL BLUE



FRENCH BLUE



DEEP VIOLET

