

# Spirit Stethoscopes

## Spirit Deluxe Dual Head Stethoscopes

Available in Adult 47mm, Paediatric 38mm and Infant 32mm (Head ø)  
Stainless Steel Dual Head Design with Floating Diaphragm and deep cone shaped bell  
Ultraflex™ binaural PVC tubing available in 20 colours.  
Includes spare diaphragm, retaining ring and soft ear tips.

Product Code	Description
STS601BLACK	STETHOSCOPE - SPIRIT DELUXE ADULT D/H BLACK (EACH)
STS606BLACK	STETHOSCOPE - SPIRIT DELUXE PAED D/H BLACK (EACH)
STS607BLACK	STETHOSCOPE - SPIRIT DELUXE INFANT D/H BLACK (EACH)



STS601PINK

## Stethoscope Colours Available

(bold colours are commonly available).  
Product codes will include colour name or shortened version of name—such as ST099LBLU—Light Blue Spirit Majestic Nurse Stethoscope.

Grey	<u>Black</u>	Dark Green	<u>Pink</u>
<u>Royal Blue</u>	<u>Yellow</u>	Dark Blue	Orange
Light Green	Burgundy	White	Teal
<u>Light Blue</u>	Magenta	Lavender	<u>Red</u>
<u>Purple</u>	Aqua	Neon Orange	Navy Blue

## Spirit Deluxe Dual Head Teaching Stethoscopes

Same features as the Spirit Deluxe D/H Adult but with 2 sets of head sets.

Product Code	Description
STS621BLACK	STETHOSCOPE - SPIRIT DELUXE TEACHING D/H BLACK (EACH)



STS621DBLU

## Spirit Majestic Nurse Stethoscopes

Aluminium Dual Head Adult design (47mm head ø) with screw on retaining ring for diaphragm.  
Simplex™ binaural PVC tubing available in 20 colours.  
Includes stethoscope with hard screw ear tips.

Product Code	Description
ST099	STETHOSCOPE - SPIRIT MAJESTIC LIGHT D/H BLACK (EACH)
ST099LBLU	STETHOSCOPE - SPIRIT MAJESTIC LIGHT D/H LIGHT BLUE (EACH)
ST099PINK	STETHOSCOPE - SPIRIT MAJESTIC LIGHT D/H PINK (EACH)



ST099

## Spirit Cardiology Dual Head Stethoscopes

Stainless Steel Dual Head Design with Floating Diaphragm and deep cone shaped bell  
Ultraflex™ binaural PVC tubing available in 8 colours.  
Includes spare diaphragm, retaining ring and soft ear tips.

Product Code	Description
STS747	STETHOSCOPE - SPIRIT CARDIOLOGY ADULT D/H BLACK (EACH)



STS747

## Spirit Deluxe Series III Cardiology Single Head Stethoscopes

Stainless Steel Single Head Design with Multi Frequency Ring-O-Raised Diaphragm  
Ultraflex™ binaural PVC tubing available in 8 colours.  
Includes spare diaphragm, retaining ring and soft ear tips.

Product Code	Description
STS748	STETHOSCOPE - SPIRIT CARDIOLOGY III ADULT S/H BLACK (EACH)



STS748



W: [www.dbhealth.com.au](http://www.dbhealth.com.au)  
E: [customerservice@denyer.com.au](mailto:customerservice@denyer.com.au)  
P: (02) 9638 2100  
ABN: 37 000 459 174

Suppliers of quality products to the Healthcare Industry est 1889. Australian Owned & Operated.