


Warwick Sterilisation Products

Warwick Re-usable Holloware Product Range

Manufactured from a high quality polypropylene material that is approved for use in medical devices.

- Registered as CE Marked Class 1 Medical Devices.
- Low Cost, Robust and hard wearing—Guaranteed for 1000 cycles
- 50% lighter than Stainless Steel & easier to clean blood protein
- Can be sterilised using Autoclave, ETO or Hydrogen Plasma
- Smooth corners & surfaces reduce the risk of piercing wraps.



WAPIT3025

Warwick Mesh & Perforated Instrument Trays

- Perforated or mesh base for improved water penetration and steam circulation,
- Select trays also have perforations in the side walls
- Strengthened sides
- Tag holes for tray tags.




WAMIT1515

WAMIT2709



WAMIT2715



WAMIT2015



WADINT2525


Product Code	Description	Carton Qty
WAMIT1515	INSTRUMENT TRAY - MESH 200 X 180 X 30MM (EACH)	50
WAMIT2015	INSTRUMENT TRAY - MESH 200 X 150 X 51MM (EACH)	60
WAMIT2709	INSTRUMENT TRAY - MESH 275 X 110 X 30MM (EACH)	60
WAMIT2715	INSTRUMENT TRAY - MESH 270 X 150 X 30MM (EACH)	60
WAMIT2817	INSTRUMENT TRAY - MESH 280 X 170 X 55MM (EACH)	25
WAMIT3025	INSTRUMENT TRAY - MESH 300 X 250 X 52MM (EACH)	24
WAMIT4030	INSTRUMENT TRAY - MESH 424 X 305 X 75MM (EACH)	12
WAPIT3025	INSTRUMENT TRAY - PERFORATED 300 X 250 X 52MM (EACH)	25
WAPIT4030	INSTRUMENT TRAY - PERFORATED 424 X 305 X 75MM (EACH)	12
WADINT2525	DENTAL INSTRUMENT TRAY - 255 X 255X 50MM (EACH)	20
WADINT4825	DENTAL INSTRUMENT TRAY - 480 X 255 X 50MM (EACH)	25




WAMIT3025



WAMIT2817



WAMIT4030



WADINT4825



W: www.dbhealth.com.au
 E: customerservice@denyer.com.au
 P: (02) 9638 2100
 ABN: 37 000 459 174

DBH1604


Suppliers of quality products to the Healthcare Industry est 1889. Australian Owned & Operated.

Warwick Sterilisation Products


Warwick Silicone Protection Pads and Rolls

- Designed to protect instruments during washing, sterilisation and transportation.
- Pliable pad with perforations for maximum water penetration and steam circulation with long fingers for instrument positioning.

Product Code	Description
WASP2313	SILICONE PAD - PINK #230 X 130MM (B/6)
WASP2823	SILICONE PAD - PINK #280 X 230MM (B/6)
WASP2719	SILICONE PAD - PINK #270 X190MM (B/6)
WASR10	SILICONE ROLL - PINK 1M #230 X1000MM (EACH)



WASP2313



WASR10

Warwick Tray Tags with S/S Rings

- Designed to enable tracking of trays .
- Tray Tags are autoclavable with area for barcode or label.

Product Code	Description
WATR30	ATTACHMENT RING FOR TRAY TAG S/S (EACH)
WATAG070	TRAY TAG - WARWICK 105 X 40MM BLACK (EACH)



WATAG***




WATR30

Tag Code	Colour	Tag Code	Colour
010	WHITE	050	GREEN
015	YELLOW	055	GREY
020	BLUE	060	LAVENDAR
022	LIGHT GREEN	065	TAN
025	ORANGE	070	BLACK
030	RED		

Warwick Container for Long Instruments and Scopes

- Designed to process up to 3 long instruments at a time.
- Each Container can also be stacked on each other
- Mesh base, sides & lid allow water penetration and steam circulation.
- Moveable silicone grips to hold instruments securely in place.
- Rapid water drainage with raised floor and internal curves.
- Positioning lock rib to centralise lid.
- Handle area designed to help lid removal.



WASMT5510

Product Code	Description
WASMT5510	DISINFECTION SET - CONTAINER + LID FOR RIGID SCOPES / LONG INSTRUMENTS (EACH)



W: www.dbhealth.com.au
 E: customerservice@denyer.com.au
 P: (02) 9638 2100
 ABN: 37 000 459 174

DBH1604