


Sentry LabelBands Wristbands

Sentry LabelBand® Wristbands








- Available in 4 Models - Adult, Adult/Paediatric, Infant and Mother-Infant Set of 2
- Non-transferable, secure, plastic snap provides maximum security
- Media support text, linear and 2D bar codes, photos and graphics
- Clear adhesive shield provides optimal protection from solvents for accurate and repeatable barcode scanning
- Large label area accommodates both laser and thermal labels
- Soft, smooth, non-abrasive material and latex free
- Water-resistant, durable and strong, yet lightweight and comfortable

Wristband Type	Wristband Size (L x W):	Bands/Box:	Max Label:
LabelBand Adult PD5080-**	298mm x 31mm	500	25mm x 76mm
LabelBand Adult/Paed PD5090-**	273mm x 31mm	500	25mm x 51mm
LabelBand Infant PD5095-11	176mm x 19mm	500	16mm x 36mm
LabelBand Mother-Infant PD5096-11 Mother Band	298mm x 31mm	250 sets of 2 - 1 Adult and 1 Infant	25mm x 76mm
LabelBand Mother-Infant PD5096-11 Infant Band	176mm x 19mm		16mm x 36mm



From top:
 PD5080-11
 PD5080-16
 PD5090-11
 PD5090-16
 PD5095-11

Colour Range for LabelBand

	11—White
	13—Blue
	14—Yellow
	16—Red
	17—Orange
	22—Kelly Green
	23—Purple

PD5080 and PD5090 colour options available, ** substitutes for colour code

Product Code	Description
PD5080-**	ID BAND SENTRY LABELBAND ADULT (B/500)
PD5090-**	ID BAND SENTRY LABELBAND ADULT/PAED (B/500)
PD5095-11	ID BAND CONFI-ID-ENT SHIELD INFANT WHITE (B/500)
PD5096-11	ID BAND CONFI-ID-ENT SHIELD DUO-BAND WHITE (B/250)
DBSD9163-11	LABEL SHEET - A4 14-UP PATIENT CHART (PKT/250)

Laser A4 Label Sheets for Medical Charts and LabelBands

Currently 1 custom A4 label sheet is available with 14 labels which has 12 labels for Medical Charts and 2 Labels for usage on the 5080 LabelBand. 12 Labels - 99 x 38mm and 2 Labels 75 x 25mm. With Printed Arrow to show how sheet should be orientated in printer. Designed for use with Laser Printers.



W: www.dbhealth.com.au
 E: customerservice@denyer.com.au
 P: (02) 9638 2100
 ABN: 37 000 459 174

DBH1612