

# Sentry Speedi-Print and SuperBand Write On Wristbands

## Sentry® SuperBand® Adult Wristband

- Write on white area
- Lightweight non-stretch polyester material and Latex Free.
- SecurSnap closure to prevent transfer
- ClearImage® surface protects handwritten and embossed information
- Water-resistant

### PD5040 Series

Information Area: 90 x 25mm

Wristband Size: 292 x 27mm


## Sentry® SuperBand® Write-On Adult/Paed Wristband

- Full colour wristband
- Lightweight non-stretch polyester material and Latex Free
- SecurSnap closure to prevent transfer
- Hospital Marker (PD1400) can be used to write-on this wristband.
- Water-resistant

### PD5070 Series

Information Area: 66 x 25mm

Wristband Size: 262 x 27mm



From top:  
PD5040-23  
PD5070-14  
PD510-11  
PD220S-12

Colour Range for Superband & ShortStay



## Speedi-Print® Infant Wristband

- Soft vinyl material and Latex Free
- SecurSnap closure to prevent transfer
- ClearImage® surface protects handwritten and embossed information

### PD510-\*\* Series

Information Area: 57 x 19mm

Wristband Size: 169 x 19mm

## Speedi-Print® Mother-Infant Wristband

- Soft vinyl material and Latex Free
- SecurSnap closure to prevent transfer
- ClearImage® surface protects handwritten and embossed information
- Serialized with alpha-numeric characters
- DuoBand - Set of 2 Wristbands - 1 Adult/Paed and 1 Infant with the same alpha-numerial serial number.

### PD220S-\*\* Series

Adult/Paed Information Area: 60 x 19mm

Adult/Paed Wristband Size: 237 x 25mm

Infant Information Area: 47 x 17mm

Infant Wristband Size: 169 x 19mm

## Short Stay® Tyvek Wristband

- Tyvek® material with tamper-evident with adhesive closure
- Ideal for day only admissions and procedure
- Write-on wristband with permanent marker or pen
- Anti-microbial Additive.


### PD3000-\*\*

Information Area: 101 x 25mm

Wristband Size: 255 x 25mm



Colour Range for Speedi-Print



Product Code	Description
PD5040-**	ID BAND SENTRY SUPER ADULT WRITE-ON (B/500)
PD5070-**	ID BAND SENTRY SUPER AD/PAED WRITE-ON (B/500)
PD510-**	ID BAND SPEEDI-PRINT INFANT IMPRINT (B/250)
PD220S-**	ID BAND SPEEDI-PRINT DUO-BAND WRITE-ON (B/150 SETS)
PD3000-**	ID BAND SECURLINE SHORT STAY WRITE-ON (B/1000)



W: [www.dbhealth.com.au](http://www.dbhealth.com.au)

E: [customerservice@denyer.com.au](mailto:customerservice@denyer.com.au)

P: (02) 9638 2100

ABN: 37 000 459 174

DBH1613