

Warwick Ward Products

Warwick Hollowware Jugs


Long Life Re-usable Polypropylene Jugs are designed for up to 1000 washes and sterilisations.

- Clear, easy to read graduated measurements.
- Lightweight, strong, durable and easy to handle.
- Can also be washed in automatic washer disinfectant or by hand using detergents and disinfectants.
- Can be sterilised by Autoclave, Ethylene Oxide or Hydrogen Plasma and Autoclavable up to 135C.
- Class 1 Medical Device (CE Marked)

Discription	Carton Qty	Product code
JUG NATURAL 500ML (EACH)	50	WAJG500
JUG NATURAL 1000ML (EACH)	45	WAJG1000
JUG NATURAL 2000ML (EACH)	40	WAJG2000



WAJG1000



WAJG2000

Warwick Holloware Spa Range

- Spa Products are designed for a lengthy service - guaranteed for 1000 washes and sterilizations.
- Scouring and abrasive cleaners will reduce the performance of products.
- Can be cleaned in automatic washer disinfectors or suitable for hand washing using detergents and disinfectants.

Warwick Spa Bedpans

Long Life Re-usable Polypropylene Bed Pan, Ideal for patients confined to bed due to illness or after surgery

- Lightweight, strong and durable - with a weight tolerance up to 250kg (WASPN has a max user limit of 160kg).
- Can be cleaned in automatic washer disinfectors or hand washed.
- Smooth material surfaces ensure ease of use and cleaning
- Class 1 Medical Device (CE Marked) and Autoclaveable.

Description	Carton Qty	Product code
BEDPAN - WARWICK BUDGET 370 X 310 X 90MM (EACH)	30	WABPN
BEDPAN - VECTOR W/ LID 335 X 308 X 108MM (EACH)	12	WAVPN
BEDPAN - SELINA SUPER-SLIPPER 335 X 245MM (EACH)	25	WASPN



WAVPN



WASPN

Warwick Spa Male and Female Urinals

Re-usable Polyethylene Female and Male Urinals

Comfortable and easy to clean.

Manual washing (maximum 90C) and processing in Washer Disinfectant Machines.

Description	Carton Qty	Product code
FEMALE URINAL - CYGNET BOTTLE W/ HANDLE 1000ML (EACH)	30	WAUF1000
URINAL MALE - TOPPER BOTTLE W/ HANDLE 1200ML (EACH)	45	WAUMC1000



WAUF1000



WAUMC1000

Warwick Spa Vomit Bowls and Denture Cup

Re-usable Polypropylene Vomit Bowl and Denture Cups

Manual washing (maximum 90C) and processing in Washer Disinfectant Machines.

Description	Carton Qty	Product code
VOMIT BOWL W/ HANDLE 300 X 65MM 1500ML (EACH)	25	WAVB300
DENTURE CUP W/ HINGED LID & LIFT OUT STAINER (EACH)	75	WADNT



WAVB300



WADNT



W: www.dbhealth.com.au
 E: customerservice@denyer.com.au
 P: (02) 9638 2100
 ABN: 37 000 459 174

DBH1618