

Warwick Gallipots, Bowls & Kidney Dishes

Long Life Re-usable Polypropylene Hollowware is Designed for up to 1000 washes and sterilisations - scouring and abrasion will reduce the performance of products.

Features of the Warwickware Hollowware range are:

- Lightweight, strong and durable.
- Improved removal of blood residue compared to Stainless Steel.
- Grip feature moulded below the rim to improve holding.
- Smooth edges and surfaces eliminate paper rips during packing and sterilization.
- Can be sterilised by Autoclave, Ethylene Oxide or Hydrogen Plasma and Autoclavable up to 135C
- Class 1 Medical Device (CE Marked)

Warwick Gallipots and Lotions Bowls

Non drip with strengthened design for secure grip and anti vacuum pips on the bases.

Description	Carton Qty	Product Code
GALLIPOT 35ML #40 X 30MM (EACH)	100	WAGP40
GALLIPOT 95ML #60 X 40MM (EACH)	300	WAGP60
GALLIPOT 280ML #80 X 60MM (EACH)	80	WAGP80
LOTION BOWL 280ML #100 X 45MM (EACH)	200	WALB100
LOTION BOWL 900ML #150 X 70MM (EACH)	60	WALB150
LOTION BOWL 2L #200 X 90MM (EACH)	30	WALB200
LOTION BOWL 3L #250 X 105MM (EACH)	20	WALB250



WALB100



WAGP60

Warwick Kidney Dishes

Non drip with strengthened design for secure grip and graduations

Description	Carton Qty	Product Code
KIDNEY DISH 300ML 150 X 75 X 35MM (EACH)	210	WAKD150
KIDNEY DISH 500ML 200 X 98 X 45MM (EACH)	150	WAKD200
KIDNEY DISH 750ML 250 X 130 X 55MM (EACH)	60	WAKD250
KIDNEY DISH 1.5L 300 X 144 X 65MM (EACH)	25	WAKD300



WAKD150

Warwick Wash Bowls

Non drip with strengthened design for secure grip and graduations

Description	Carton Qty	Product Code
WASH BOWL 3L 310 X 118MM (EACH)	25	WAWB3112
WASH BOWL 5L 350 X 120MM (EACH)	25	WAWB3513



WAWB3513