

Grip N Assist Mobility Assistance Belts

Grip N Assist Professional Mobility Assistance Belts

Designed to provide the patient and nursing staff with Mobility Assistance Belt that provides comfort and safety.

All models have the following features:

- Ergonomic Handles—Reinforced and easy to hold handles
- Reinforced belt designed to reduce slippage and be comfortable for the wearer.
- Disinfect with Hospital Grade Solution
- Latex Free and Waterproof
- Adjustable to suitable individual patients with a secure and reliable velcro strap and locking buckle clip
- Available in 2 ranges from Standard to Lift Capable Ranges

Usage

- Can be worn backwards or forwards by both nursing staff or patient
- Ideally for transfers, moving from seating to standing and supporting movement.

Standard Range

Grip N Assist Lite

The GNA LITE is designed specifically for use with small adults and young adults.

Fits waist size: 61cm to 95cm
Handles: 2 Handles



LCLPGNAW006

Grip N Assist Pro White and Blue

The GNA PRO is designed for adults and is the most popular of all Grip N Assist Mobility Assistance Belts.

PRO WHITE

Fits waist size: 68.5cm to 137cm
Handles: 4 Handles



LCPCNAW001

PRO BLUE

Fits waist size: 68.5cm to 137cm
Handles: 5 Handles



LCPCNAW001-B

Product Description	Product Code
TRANSFER BELT - GRIP N ASSIST LITE #WHITE #61-96.5CM # 2 HANDLES (EACH)	LCLGNAW006
TRANSFER BELT - GRIP N ASSIST PROFESSIONAL #WHITE #68.5—137CM # 4 HANDLES (EACH)	LCPCNAW001
TRANSFER BELT - GRIP N ASSIST PROFESSIONAL #BLUE # 68.5—137CM #5 HANDLES(EACH)	LCPCNAW001-B

Grip N Assist Pro Blue worn backwards



Grip N Assist Mobility Assistance Belts

Lift Capable Range

These mobility assistance belt include the following features to provide additional safety and security

- Centre D-Ring which can be used with Lifting Equipment.
- Harness which can be employed as leg support or upper body support

Grip N Assist Pro Petite

The GNA PRO PETITE is designed specifically for use with small adults and young adults.

Fits waist size: 53cm to 81cm
Handles: 3 Handles



LCPGNA600

Grip N Assist Pro Lift

The GNA PRO LIFT is designed for adults with the included benefit of upper body or leg supports.

Fits waist size: 68.5cm to 137cm
Handles: 5 Handles



LCPGNA200

Grip N Assist Pro Plus

The GNA PRO PLUS is designed for large adults with the included benefit of upper body or leg supports.

Fits waist size: 101cm to 183cm
Handles: 6 Handles



LCPGNA300

Note: Please ensure for larger adults that proper technique is used and it is not recommended to use on patients that can't assist the nursing staff.

Grip N Assist Pro Lift showing different positions on wearer.

A—Worn with front of Belt on the Back of the wearer with Upper Support—View from Back of Wearer

B—Same as A—View from Front of Wearer

C—Worn with front of Belt on the front of the wearer with Leg Support—View from Back of Wearer

D—Same as C—View from Front of Wearer



Description	Product Code
TRANSFER BELT - GRIP N ASSIST PRO PETITE #TEAL #53-81CM WAIST (EACH)	LCPGNA600
TRANSFER BELT - GRIP N ASSIST PRO PLUS #ORANGE #101-183CM WAIST # (EACH)	LCPGNA300
TRANSFER BELT - GRIP N ASSIST PRO LIFT #YELLOW #68.5-137CM WAIST #5 HANDLES (EACH)	LCPGNA200



W: www.dbhealth.com.au
E: customerservice@denyer.com.au
P: (02) 9638 2100
ABN: 37 000 459 174

Suppliers of quality products to the Healthcare Industry est. 1889. Australian Owned & Operated.