

Examination Lights

Heine EL10 LED Examination Light

The Heine EL 10 Examination Light distinguishes itself through an absolutely homogenous field of light. The light intensity can be regulated continuously between 6,500 and 45,000 Lux. The spot size can be adjusted in a diameter between 8.5 and 14cm to suit any situation at a working distance of 30cm.

- Controls with ergonomic grip allow for intuitive and fast adjustment of brightness and spot size. The compact illumination head (Ø approx. 60mm) allows almost coaxial illumination.
- Bright and homogeneous. Perfect illumination in every examination situation.
- 10 Watt high performance LED with authentic colour rendering
- Easy and efficient cleaning and disinfection. Easy and intuitive adjustments via coiled goose neck.
- An electrical control system monitors all functions and guarantees safety.

Description	Product code
EXAM LIGHT - HEINE EL10 WHITE LED W/ WALL MT (EACH)	HJ0827001
EXAM LIGHT - HEINE EL10 WHITE LED W/ CLAMP MT (EACH)	HJ0827002
EXAM LIGHT - HEINE EL10 WHITE LED W/ MOBILE STAND (EACH)	HJ0827003



HJ0827003



Adjusting spot size on EL10

Heine EL3 LED Examination Light

- The EL3 LED provides bright shadow free illumination precisely at the point of interest. At a working distance of 30cm > 42,000 lux / >350 lumens.
- The LED illumination provides optimal light intensity, homogeneity and colour rendering for the most accurate diagnosis.
- Compact Illumination Head (diameter 6cm) allows for convenient, precise positioning, particularly in difficult application situations.
- Cool Touch closed head design makes disinfection safe and easy.
- The lightweight head stable gooseneck allow easy and intuitive adjustment to keep the light in the required position.

Description	Product code
EXAM LIGHT - HEINE EL3 WHITE LED W/ WALL MT (EACH)	HJ0827011
EXAM LIGHT - HEINE EL3 WHITE LED W/ CLAMP MT (EACH)	HJ0827013
EXAM LIGHT - HEINE EL3 WHITE LED W/ MOBILE STAND (EACH)	HJ0827014



HJ0827011

FlexLED 15 Examination Light

The FlexLED is a LED Examination Light perfect for use in all medical examinations.

- Seven powerful long-life LED's with focused beam
- LED's provide 50,000 Lux at 50cm
- 1000mm Long flexible gooseneck arm
- Variable intensity and colour control
- Low heat and no infrared radiation with High energy efficiency
- Low maintenance. No globes to change.

Description	Product code
EXAM LIGHT - FLEXLED 15 W/ WALL MT (EACH)	DA1500
EXAM LIGHT - FLEXLED 15 W/ MOBILE STAND (EACH)	DA1501

DA1501



W: www.dbhealth.com.au
 E: customerservice@denyer.com.au
 P: (02) 9638 2100
 ABN: 37 000 459 174

Suppliers of quality products to the Healthcare Industry est. 1889. Australian Owned & Operated.

Examination Lights

Maggylamp Halogen Examination Light

Our Maggylamp exam lights are an entry level range of lamps that offers a cost effective solution for surgery examination lighting.

- 700mm gooseneck.
- Sure grip handle with on/off switch.
- Warm white colour 2300K colour.
- 2 years warranty.
- Wall or mobile mounted (View all options below).
- 4,000 hour bulb life.
- 45,000 lux at 30cm.
- 35 watt, 24° beam.



RBML308

Maggylamp LED Examination Light

The Maggylamp LED Examination Light with its flexible gooseneck arm and specifically selected colour corrected LED's offers flexibility, correct light output and colour required for surgery examinations.

- 40,000 hour lamp life.
- Cool white light, 5500K colour.
- Low energy consumption, <5 watt.
- 2 years warranty.
- Low heat output.
- 3 x 1 watt LED giving 16,000 lux at 30cm.
- Concise beam angle of 30°.
- Wall or mobile mounted.
- 700mm gooseneck.



RBML308LED

Description	Product code
EXAM LIGHT - MAGGYLAMP HALOGEN (EACH)	RBML308
EXAM LIGHT - MAGGYLAMP LED (EACH)	RBML308LED
OPTION—WALL BRACKET RIGHT ANGLE FOR MAGGYVUE & MAGGYLAMP(EACH)	RBML307
OPTION—MOBILE STAND FOR MAGGYVUE & MAGGYLAMP #95CM (EACH)	RBML130

Maggyvue Opal Magnifier Lamp

The Maggyvue Magnifier Lamp is a quality German designed lamp. It is very popular with the medical, electronics and inspection industries where sharp definition is required.

- Quality German optics.
- Circular lens, Ø 120mm, 2 x magnification.
- Optional clip on lens to increase magnification #ML151.
- Spring balanced joints, easy to adjust.
- Supplied with bench clamp and diffuser.
- 22 watt daylight colour fluorescent tube.



RBML150

Description	Product code
MAGNIFIER LAMP—MAGGYVUE OPAL (EACH)	RBML150
REPLACEMENT TUBE—22W FOR ML150 (EACH)	RBML235
OPTION—BENCH CLAMP FOR MAGGYVUE #0-65MM (EACH)	RBML110
OPTION—WALL BRACKET RIGHT ANGLE FOR MAGGYVUE & MAGGYLAMP (EACH)	RBML307
OPTION—MOBILE STAND FOR MAGGYVUE AND MAGGYLAMP #95CM (EACH)	RBML130
OPTION—CLIP ON LENS FOR MAGGYVUE #3X MAGNIFICATION (EACH)	RBML151