

Examination Couches

Fixed Height Examination Couch

The Fixed Height Exam Couch are cost effective exam couches essential for consulting and treatment rooms. Combining comfort, stability, strength and aesthetic design ensure that these exam tables remain an essential part of your surgery.

Features of the DA1000 2 Section Exam Couch are:

- Sturdy rectangular steel frame and easy to assemble—flat packed (4 bolts)
- Adjustable backrest with Vinyl covered top available in a Grey or Black.
- Available in Grey Powder Coat
- Dimensions: 185 x 61 x 80cm (High) and Max User Weight 175kg



DA1000

Features of the DA1004 2 Section Heavy Duty Exam Couch are:

- Sturdy rectangular steel frame and easy to assemble—flat packed (boltless)
- Adjustable backrest with Vinyl covered top available in a Grey or Black.
- Available in Grey Powder Coat
- Dimensions: 183 x 61 x 78cm (High) and Max User Weight 250kg



DA1004

Product	Product Code
EXAM COUCH - HEAVY DUTY P/C #250KG SWL (EACH)	DA1004
EXAM COUCH - STANDARD P/C #175KG SWL (EACH)	DA1000

High-Low Examination Couch

The High-Low Examination Couches provide a cost effective alternative for medical professionals on a budget.

- Safe working load - 200kg
- Sturdy and reliable heavy duty construction.
- Quiet electronic motor controlled by foot-switch operation for effortless positioning.
- Retractable castors available on all models.
- Height adjustment from 46 to 91cm.
- Dimensions (approx.): 195.5 x 71cm.

Product	Product Code
EXAM COUCH - ELECTRIC HI-LO P/C #2 SECTION W/ MANUAL HEADREST (EACH)	DA1242M
EXAM COUCH - ELECTRIC HI-LO P/C #2 SECTION W/ ELEC HEADREST (EACH)	DA1242E
EXAM COUCH - ELECTRIC HI-LO P/C #3 SECTION (EACH)	DA1243M
EXAM COUCH - ELECTRIC HI-LO P/C #3 SECTION #3 MOTORS (EACH)	DA1243E



Retractable Castors



DA1242E



DA1243M



DA1243E

Examination Couches

Access Fixed Height Exam Table

The Access Exam Table is a fixed height exam table without electrically functions. The exam table is very stable and features an inbuilt step to safely allow patients to get on or off this couch. The backrest is adjusted via a hydraulic adjustment and includes stirrups for gynaecological examinations..



MIACCESS

Product	Product Code
EXAM TABLE - BREWER ACCESS W/MANUAL BACKREST (EACH)	MIACCESS

Access High-Low Exam Table

The Access High-Low Exam Table has been designed to transform from a true chair position into the exam table position allowing patient to be easily moved from consultation to examination.

- Abundant storage with pass-through drawers.
- Ergonomical integrated stirrups.
- Comfortable easy to clean seamless upholstery.
- Safety grab bars to enable wheelchair transfers with ease.
- Height adjustment from 94cm down to a low of 46cm.
- Electrical adjustment of height and backrest.
- Pelvic tilt and drawer heater.



mi6501

Product	Product Code
EXAM TABLE - BREWER ACCESS HIGH-LOW (EACH)	MI6501

Assist Procedure Exam Table

The Assist Procedure Exam Table has resulted in a Exam Table that is the lowest, safest, and easiest to use power procedure table on the market.

Features of this Exam Table are:

- The Brewer Assist was the first to lower to 48cm, making wheelchair transfers and patient accessibility effortless.
- Brewer Assist also safely supports 204kg.
- Convenient leg rest that extends an additional 12.7cm.
- Additionally this exam table features a wider range of adjustment including headrest and positions from the standard position to specialised positions for lithotomy, proctology and gynaecology.



DA1000

Product	Product Code
EXAM TABLE - BREWER ASSIST POWER PROCEDURE (EACH)	MI7000

Examination Couches

Cupboard and MultiBank Melamine Exam Couch

The Cupboard and MultiBank Melamine Exam Couch have been designed and manufactured to enable easy of use, comfort and allow for storage for healthcare professionals.

- Constructed from impact and abrasion resistant melamine the surface of the exam table is sealed and is moisture resistant.
- The vinyl tops have a high level of heat / flame retardant and can be cleaned with mild detergent and antiseptics.
- Fully adjustable dual ratchet backrest
- Choice of vinyl and melamine colours.
- Dimensions on all couches without options—178.5 x 56.5 x 85.5cm


Cupboard Exam Couch - 2 Sliding Doors (full length), 1 Instrument Drawer and Retractable armslide

MultiBank Exam Couch - 2 Hinged Doors, 2 banks of drawers at either end of the couch and Retractable armslide


Available options:

- Gynaecological End with footrest,
- Pelvic Tilt,
- Face Hole in Backrest,
- Conceal Drawers behind sliding doors
- Lockable Doors,
- Liftup Extension
- Large range of Melamine and Vinyl Colour available.

Product	Product Code
COUCH EXAM - WOODEN CUPBOARD (EACH)	BBCUPBOARD
COUCH EXAM - WOODEN MULTIBANK (EACH)	BBMBANK



BBCUPBOARD



BBMBANK