

X-Ray Viewers

Dalcross Slimline LCD X-ray Viewers

- Surface Mounted Single, Double and Triple Bay Slimline LCD X-Ray Viewers
- Stylish design with 2 key-hole slots in rear for wall mounting
- Supplied with plug lead that exits at bottom right of light box
- High brightness and whiteness with even, flicker-free illumination
- Aluminium housing with good heat dissipation
- Positive grip holds X-ray film secure
- Maintenance free. No globes to change

Description	Product Code
X-RAY VIEWER - SINGLE SLIM LINE LCD #47 X 54.5X 2.4CM (EACH)	DA1600SL
X-RAY VIEWER - DOUBLE SLIM LINE LCD #85 X54.5 X2.4CM (EACH)	DA1601SL
X-RAY VIEWER - TRIPLE SLIM LINE LCD #123 X54.5 X2.4CM (EACH)	DA1602SL



DA1600SL



DA1601SL



DA1602SL

Dalcross Halogen X-ray Viewer

- Surface Mounted Single, Double and Triple Bay Halogen X-ray Viewers
 - 2 Key-hole slots in rear for wall mounting
 - Single and double can be upright on base, for use on desk or bench. Triple Bay viewer is Wall Mounted only
 - Supplied with plug lead that exits at bottom right of light box
 - Excellent light diffusion and uniform illumination is achieved from: circular fluorescent tube, continuous opal perspex viewing panel and internal white reflector surface
 - Stainless steel strip contains grip rollers to hold X-rays and Each individually switched bay is internally divided
- White powdercoated steel case is easy to clean and Made in Australia

Description	Product Code
X-RAY VIEWER - SINGLE #37.7X50X10.5CM (EACH)	DA1600
X-RAY VIEWER - DOUBLE #75X50X10.5CM (EACH)	DA1601
X-RAY VIEWER - TRIPLE #112.5X50X10.5CM (EACH)	DA1602

Also available in Flush Mount and on a mobile stand.

Available in 4 in a line, 6 in a line, 8 in a line, and spinal X-Ray. Also available in 2 tiers—2 x 2, 3x 3 and 4 x4.



DA1600



DA1601



DA1602