

# Fiab Eco Gel Ultrasound Gel

Gel ECO is a good quality echo-transmission gel for Ultrasound and IPL.

## Features

- Gel ECO is push-vacuum-produced and the gel is therefore free from air microbubbles making the resultant Ultrasound results more accurate.
- Particular structure assists the user in movement of the Ultrasound Probe thus making the users work easier.
- ECO Gel doesn't stain, is not greasy, does not ionize and does not oxidize when in normal use for Ultrasound Scanning.
- Does not contain any active ingredients that might damage the probe transducer.
- The preservatives used are non-toxic and strictly in accordance with the quality standards.
- Salt-free and without string effect.
- Does not contain formaldehyde,
- Irritant Free,
- Water Soluble,
- Drying-Resistant,
- Odourless,
- Easily removable,
- Transparent and Clear Appearance,
- Alcohol Free

## Storage Requirements

- Storage Temperature: +10°C - +40°C -
- Storage Relative humidity: 30% - 75%
- Expiry Length—36 months from production date
- Packaging Material — for 260ml in PE nontoxic bottle, open-close cap

ARTG Number 206370

MSDS available please contact Customer Service today for your copy.



G008A

Product Code	Description	Carton Quantity
G008A	ULTRASOUND GEL - GEL ECO CLEAR 260ML (EACH)	Ctn/40
G0084	ULTRASOUND GEL - GEL ECO CLEAR 5L (EACH)	Ctn/4



**W:** [www.dbhealth.com.au](http://www.dbhealth.com.au)  
**E:** [customerservice@denyer.com.au](mailto:customerservice@denyer.com.au)  
**P:** (02) 9638 2100  
**ABN:** 37 000 459 174

Suppliers of quality products to the Healthcare Industry est 1889. Australian Owned & Operated.