

Medication Carts

Dosage Medication Carts

The 30 and 50 Dosage Medication Trolley are designed to store medication in secure fashion. Each of the small 10cm include aluminium compartments (maximum of 12 compartments per drawer) to separate patient medication.

Features of this medication trolley are:

- Central locking
- 25cm deep bottom drawer
- Rail on 3 sides
- Drawer using full extension roller bearing runners that gives access to entire drawer
- Wheels: 12.5cm standard with 2 locking
- Available in Stainless Steel or Powder Coated in a range of colours.

30 Dosage - 3 x 10cm Drawers and 1 x 25cm Drawer

50 Dosage - 5 x 10cm Drawers and 1 x 25cm Drawer



OX18-4660



OX55-4560

Description	Code
MEDICATION TROLLEY - OXFORD 30 DOSAGE S/S W/ 4 DRAW #60 X 50 X 100CM (EACH)	OX18-4560
MEDICATION TROLLEY - OXFORD 50 DOSAGE S/S W/ 6 DRAW #60 X 50 X 100CM (EACH)	OX18-4660
MEDICATION TROLLEY - OXFORD 30 DOSAGE P/C W/ 4 DRAW #60 X 50 X 100CM (EACH)	OX55-4560
MEDICATION TROLLEY - OXFORD 50 DOSAGE P/C W/ 6 DRAW #60 X 50 X 100CM (EACH)	OX55-4660

Unit Medication Carts

The 9-18, 12-24, 15-30 and 20-40 Unit Medication Trolley are designed to store medication in secure fashion. Each of the small 10cm drawers include aluminium compartments (maximum of 12 compartments per drawer for 60cm wide and 15 compartments for 75cm wide Carts) to separate patient medication. Additionally includes Stainless Steel top with insert tray, document basket and lift-up side shelf

- Central locking
- Rail on 2 sides
- Drawers using full extension roller bearing runners that gives access to entire drawer (40kg capacity),
- Wheels: 12.5cm standard with 2 locking.
- Available in Powder Coat from a range of colours.

9-18 Unit - 3 x 10cm Drawer, 1 x 20cm and 1 x 25cm Drawer at bottom

12-24 Unit - 3 x 10cm Drawer, 1 x 20cm and 1 x 25cm Drawer at bottom

15-30 Unit - 5 x 10cm Drawer and 1 x 25cm Drawer at bottom

20-40 Unit - 5 x 10cm Drawer and 1 x 25cm Drawer at bottom



OX55-1224



OX55-1530

Description	Code
MEDICATION TROLLEY - OXFORD 9-18 UNIT P/C #60 X 50 X 100CM (EACH)	OX55-0918
MEDICATION TROLLEY - OXFORD 12-24 UNIT P/C #75 X 50 X 100CM (EACH)	OX55-1224
MEDICATION TROLLEY - OXFORD 15-30 UNIT P/C #60 X 50 X 100CM (EACH)	OX55-1530
MEDICATION TROLLEY - OXFORD 20-40 UNIT P/C #75 X 50 X 100CM (EACH)	OX55-2040

Medication Carts

Half Size Rounds Carts

The ½ Size Rounds Medication Trolleys are designed for very small group of patient requiring medication.

- Lockable sloping front compartment
- Drawer 20cm deep keyed alike as front sloping compartment
- Lift up side shelf 30cm long
- Push handle on one side.
- Available in Stainless Steel or Powder Coated in variety of colours.

Description	Code
ROUNDS TROLLEY - OXFORD 1/2 SIZE S/S W/ DRAWER #50 X 50 X 119CM (EACH)	OX18-9670
ROUNDS TROLLEY - OXFORD 1/2 SIZE P/C W/ DRAWER #50 X 50 X 119CM (EACH)	OX55-9670



OX55-4560

Kerry Webster Medication Carts

The Kerry Flexicart Webster 2 Drawers and 3 Drawers Medication Trolley are designed to store medication packaged in Webster drug system in secure fashion.

- 2 or 3 drawers which are all 26cm high which allows the webster packs to be stored upright in each drawer.
- Central locking, Push handle
- Full wrap-around bumper
- Wheels: 12.5cm heavy duty swivel
- Choice of ABS top (plain or pill cup extrusion).
- Available with a variety of options - below are 3 common options that are selected with these Carts
- Available in Epoxy Coated Only from a range of colours.

Description	Code
MEDICATION TROLLEY - KERRY WEBSTER FLEXICART 3 DRAWER P/C #63.5 X 48.5 X 104.5CM (EACH)	KH220WEB
MEDICATION TROLLEY - KERRY WEBSTER FLEXICART 2 DRAWER P/C #63.5 X 48.5 X 101.5CM (EACH)	KH222WEB
OPTION - KERRY BASKET (EACH)	KH22B
OPTION - KERRY FOLDING SHELF W/ FLAT TOP (EACH)	KH22FSF
OPTION - KERRY WASTE BIN (EACH)	KH22W



KH220WEB with options folding shelf and waste bin




KH222WEB with options folding shelf and waste bin

Colour Options for Medical Carts

Oxford Powdercoat Colours

Below are the colours for Oxford Products along with white gloss. Please colours are approximate representations of the final powdercoating colour.

When selecting colour please advise drawer colour and body colour.



Standard Powdercoat Colours



Non-Standard Powdercoat Colours (additional cost may apply)

Flexicart Powdercoat Colours

Below are the standard colours for Flexicart Medication Trolleys along with White and Signal Red.

