

# Procedure Carts

## Emergency Crash Cart

The Emergency Crash Cart is a durable and fully featured Emergency Cart ready for a range of possibilities.

- Rails: 3 sides standard
- Drawer: 3 x 10cm, 1 x 15cm and 1 x 25cm deep that use full extension roller bearing runners that gives access to entire drawer.
- Top 3 drawer have removable drawer dividers.
- Lift-Up Side Shelf: 30 (L) x 50cm (W).
- Wheels: 12.5cm standard single with 2 locking.
- Central Locking - available with key lock or key lock and tamper seal.
- Adjustable IV Pole, Oxygen Cylinder Holder for C-size and Backboard.
- Powder Coated with Red Body and White Drawers



OX55-4471

Description	Code
EMERGENCY TROLLEY - OXFORD CRASH CART P/C #71 X 50 X 94CM (EACH)	OX55-4471
OPTION - OXFORD ROTATING TRAY FOR DEFIB (EACH)	OX60-0320
TAMPER SEAL PADLOCK - OXFORD #DISPOSABLE (B/100)	OX60-5160

## Paediatric NICU Cart

Strong, fully welded construction with S/Steel recessed top surface.

- Rail on rear side standard and removable double maxi-bin rails
- Drawers: 1x 100, 9x 75mm deep with 40kg capacity on full extension slides
- Keyed central locking drawers
- Top drawer with 3x full width divider storage trays, 2nd drawer without divider storage trays and 3rd -10th drawers: with 2x full width divider storage trays
- Standard colour: Broselow Hinkle trolley and drawers
- 350 x 550mm fold down side shelf
- 2x rear IV sockets with 760mm IV ramshead pole
- Castor features: 125mm single rubber with 2x braking with bumpers



OX55-4301

Description	Code
ICU TROLLEY - OXFORD PAEDIATRIC EMERGENCY P/C # 80 x 50 x 104.5CM (EACH)	OX55-4301

## Endoscopy Carts

The range of Endoscopy Carts is designed to provide safe storage of your valuable endoscopy equipment.

- Rails: 3 sides standard.
- Drawers: Heavy duty that use full extension roller bearing runners that gives access to entire drawer with 40kg capacity.
- Wheels: 12.5cm standard single with 2 locking.
- Available in Stainless Steel only



OX18-5980



OX18-5710



OX18-4210

Description	Code
ENDOSCOPY TROLLEY - OXFORD S/S W/6 DRAWS #90 X 50 X 90CM (EACH)	OX18-6290
ENDOSCOPY TROLLEY - OXFORD S/S W/ 4 DRAWS #100 X 50 X 90CM(EACH)	OX18-4210
ENDOSCOPY TROLLEY - OXFORD S/S W/ 4 DRAWERS #100 X 50 X 90CM (EACH)	OX18-5980
ENDOSCOPY TROLLEY - OXFORD ST LUKES S/S W/ 5 DRAWERS #100 X 50 X 90CM(EACH)	OX18-5710



OX18-6290

# Procedure Carts

## General Purpose Equipment Cart

The General Purpose Equipment Cart is designed to provide safe storage of your valuable medical equipment.

- Rails: 3 sides standard.
- Drawers: 5 x 10cm and 1 x 20cm deep that use full extension roller bearing runners
- Wheels: 12.5cm standard single with 2 locking.
- Central Locking. Corner Bumpers.
- Available in Stainless Steel or Powder Coated from a range of colours.

Description	Code
EQUIPMENT TROLLEY - OXFORD GENERAL PURPOSE P/C #105 X 50 X 91(EACH)	OX55-4300
EQUIPMENT TROLLEY - OXFORD GENERAL PURPOSE S/S #105 X 50 X 91(EACH)	OX18-4300



## ICU Recovery Cart

The ICU Recovery Cart is designed for use in ICU to be used as a mobile table especially when required at night.

- Drawers: 5 x 10cm and 1 x 25cm deep.
- Tilt adjusted writing platform.
- Wheels: 12.5cm standard single with 2 locking.
- Available in Stainless Steel and Powder Coated

Description	Code
ICU RECOVERY TROLLEY - OXFORD S/S W/ 6 DRAWS #60 X 50 X 100CM (EACH)	OX25-5100
ICU RECOVERY TROLLEY - OXFORD P/C W/ 6 DRAWS #60 X 50 X 100CM (EACH)	OX55-5100



## ICU Workstation - St Vincent's

The ICU Workstation is designed for use in ICU specifically to hold essential supplies.

- Drawers Left: 2 x 12.5cm, 1 x 16.5cm and 1 x 25cm deep.
- Drawers Right: 2 x 12.5cm
- Shelves Right: 1 x 15cm and 1 x 24cm clearance
- Double Maxi Bin Rail (Bins not included)
- Push Handle and Sharp container bracket
- Wheels: 12.5cm standard single with 2 locking.
- Available in Powder Coat Only from a range of colours.

Description	Code
ICU WORKSTATION - OXFORD ST VINCENTS P/C #110 X 60 X 90CM (EACH)	OX55-5200



## ICU Basket Trolleys

- Stainless Steel Construction design to hold nylon coated baskets for storage in ICU Wards.
- Available in Single, Double and Double Enclosed styles.
- Basket Sizes: 43 x 53 x 8.5cm and 43 x 53 x 18.5cm
- Push Handle on Double Enclosed Basket Trolley only.
- Wheels: 12.5cm standard single with 2 locking.
- Available in Powder Coat Only from range of colours.



Description	Code
ICU BASKET TROLLEY - OXFORD S/S SINGLE #49 X 54 X 100CM (EACH)	OX25-1000
ICU BASKET TROLLEY - OXFORD S/S DOUBLE #96 X 54 X 100CM (EACH)	OX25-2000
ICU BASKET TROLLEY - OXFORD S/S DOUBLE ENCLOSED #96 X 50 X 90CM (EACH)	OX25-3000



# Procedure Carts

## Plaster Cart

The Plaster Cart is designed for use in busy orthopaedic area and with some unique features can provide significant benefits to staff working in this area.

- Drawer: 5 x 15cm and 1 x 30cm deep with the drawer configured to make best use of the space available.
- Lift up side shelf and Waste Bin.
- Recessed top to easily contain the mess.
- Large bowl attached to one side of the trolley for easy hydration of the plaster rolls.
- Wheels: 12.5cm standard single with 2 locking
- Available in Stainless Steel



OX18-4890

Description	Code
PLASTER TROLLEY - OXFORD S/S W/ 6 DRAWERS #90 X 60 X 84CM (EACH)	OX18-4890

## Resuscitation Carts

The Range of Resuscitation Carts are designed for use in storage of critical supplies for resuscitation as well as equipment.

- Rails: 3 sides standard.
- Bottom shelf
- Pullup side shelf
- Wheels: 12.5cm standard with 2 locking
- Available in Stainless Steel only



OX18-5360



OX18-5080



OX18-6480

Description	Code
RESUSCITATION TROLLEY - OXFORD 2 LVL S/S W/3 DRAW & SIDE SHELF#60 X 50 X 90CM (EACH)	OX18-5360
RESUSCITATION TROLLEY - OXFORD 2 LVL S/S W/5 DRAW & SIDE SHELF #80 X 50 X 90CM (EACH)	OX18-5080
RESUSCITATION TROLLEY - OXFORD 2 LVL S/S W/5 DRAW, DOUBLE MAXIBIN RAIL & SIDE SHELF #80 X 50 X 90CM (EACH)	OX18-6480

## Urology Cart

- Strong, fully welded construction
- Rails on 3 sides standard
- Drawer features: lightweight with 40kg capacity full extension slides
- Drawers: 2x 100, 2x 200 and 1x 300mm deep
- Manufactured from 304 stainless steel or from mild steel with zinc coating
- Standard colour: blue weave trolley and pearl white drawers
- Castor features: 125mm single rubber with 2x braking with bumpers



OX55-4410


Description	Code
UROLOGY TROLLEY - OXFORD S/S W/ 5 DRAWS #105 X 55 X 112CM (EACH)	OX18-4410
UROLOGY TROLLEY - OXFORD P/C W/ 5 DRAWS #105 X 55 X 112CM (EACH)	OX55-4410

# Colour Options for Medical Carts

## Oxford Powdercoat Colours

Below are the colours for Oxford Products along with white gloss. Please colours are approximate representations of the final powdercoating colour.

When selecting colour please advise drawer colour and body colour.



Standard Powdercoat Colours



Non-Standard Powdercoat Colours (additional cost may apply)

## Flexicart Powdercoat Colours

Below are the standard colours for Flexicart Medication Trolleys along with White and Signal Red.

