

Operating Theatre Carts

Anaesthetic Carts

The Range of Anaesthetic Carts are designed to meet the requirements for Anaesthetics.

- Manufactured from 304 grade stainless steel to a satin finish
- Central locking drawer with full extension roller bearing runners that gives access to entire drawer - 40kg capacity
- Rails: 3 sides standard
- Wheels: 12.5cm standard single with 2 locking
- Available in Stainless Steel or Powder coated from A range of colours



OX55-4060



OX18-4065



OX18-4090



OX55-4070

Description	Code	Description	Code
ANAESTHETIC TROLLEY - OXFORD S/S W/ 6 DRAW & SHELVE #90 X 50 X 90CM (EACH)	OX18-4090	ANAESTHETIC TROLLEY - OXFORD P/C W/ 6 DRAW & SHELVE #90 X 50 X 90CM (EACH)	OX55-4090
ANAESTHETIC TROLLEY - OXFORD S/S W/ 6 DRAW & TOP SHELVE #60 X 50 X 90CM (EACH)	OX18-4070	ANAESTHETIC TROLLEY - OXFORD P/C W/ 6 DRAW & TOP SHELVE #60 X 50 X 90CM (EACH)	OX55-4070
ANAESTHETIC TROLLEY - OXFORD S/S W/ 5	OX18-4060	ANAESTHETIC TROLLEY - OXFORD P/C W/ 5 DRAW #60 X 50 X 90CM (EACH)	OX55-4060
ANAESTHETIC TROLLEY - OXFORD S/S W/ 4 DRAW #60 X 50 X 90CM(EACH)	OX18-4065	ANAESTHETIC TROLLEY - OXFORD P/C W/ 4 DRAW #60 X 50 X 90CM(EACH)	OX55-4065

Mayo Instrument Tables

The Mayo Instrument Tables are designed to be height adjustable, table top rotatable or non-rotatable and mobile. Available in 2 styles with 2 sizes small and large they are designed depending on the requirements of the user.

- Height Adjustment: 90 - 130cm.
- Tray Dimensions - Small: 40 (L) x 55 (W) x 2cm (H).
- Tray Dimensions - Large: 70 (L) x 50 (W) x 3cm (H).
- Manufactured from 304 grade stainless steel to a satin finish.
- Rotating Top - Wheels: 3 castors - 7cm standard.
- Non-Rotating Top - Wheels: 4 castors - 7cm standard.



OX10-3612



OX10-3615

Description	Code
MAYO INSTRUMENT TABLE - OXFORD SML W/ROT. TOP #60 X 67 X 90-130CM (EACH)	OX10-3600
MAYO INSTRUMENT TABLE - OXFORD LGE W/ROT. TOP #60 X 95 X 90-130CM (EACH)	OX10-3612
MAYO INSTRUMENT TABLE - OXFORD SML W/NON-ROT #51 X 60 X 90-130CM (EACH)	OX10-3602
MAYO INSTRUMENT TABLE - OXFORD LGE W/NON-ROT #51 X 70 X 90-130CM (EACH)	OX10-3615

Operating Theatre Carts

Instrument Tables, Arm Table, Instrument Over-Table and Gynae Table

The Nest of Instrument Tables comprise of 3 stainless steel instrument tables. These tables are available as either individually or as a nest of 3 tables.

- The features of these tables are:
- Wheels: 7cm standard.
- Can be stacked to minimise space usage.

The Instrument Over-Table, Gynae Table and Arm Table are height adjustable Stainless Steel tables.

- The Arm Table includes radiolucent top.
- The Gynae Table includes one curved side.
- Wheels: 12.5cm standard with 2 locking
- Height adjustment with hand tightened knobs for each leg.



OX19-2790



OX19-2675



OX18-2500



OX19-2480

Description	Code
INSTRUMENT TABLE - OXFORD 1 LVL NEST S/S (SET)	OX18-2500
INSTRUMENT TABLE - OXFORD 1 LVL S/S LARGE #105 X 52 X 92CM (EACH)	OX18-2510
INSTRUMENT TABLE - OXFORD 1 LVL S/S MEDIUM #95 X 50 X 88CM (EACH)	OX18-2595
INSTRUMENT TABLE - OXFORD 1 LVL S/S SMALL #85 X 48 X 84CM (EACH)	OX18-2585
ARM TABLE - OXFORD RADIOLUCENT TOP #75 X 40 X 75-119CM (EACH)	OX19-2675
INSTRUMENT OVERTABLE - OXFORD ADJ. S/S #90 X 50 X 75-125CM (EACH)	OX19-2790
GYNAE TABLE - OXFORD S/S #80 X 45 X 70-125CM (EACH)	OX19-2480

Bowl, Bucket and Kick Bucket Stands

The Bowl Stands are designed to hold either 1 or 2 stainless steel bowls (30 or 35cm diameter bowls).

- These stands are mounted on a 5 leg base with 5 wheels.
- Only available in stainless steel.

The Bucket Stand is designed to provide a stable stand for a stainless steel bucket.

- The stand base has 5 legs and only being 60cm high is very stable.

The Kick Bucket Stand is designed to provide a stable stand for a range of round buckets.

- The stand has a diameter of 30cm and a height of 37.5cm.



OX10-3936



OX10-3930



OX10-3820



OX10-3840

Description	Code
BOWL STAND - OXFORD SINGLE 30CM S/S (EACH)	OX10-3930
BOWL STAND - OXFORD DOUBLE 30CM S/S (EACH)	OX10-3931
BOWL STAND - OXFORD SINGLE 35CM S/S (EACH)	OX10-3935
BOWL STAND - OXFORD DOUBLE 35CM S/S (EACH)	OX10-3936
BUCKET STAND - OXFORD STD S/S W/ ALLOY BASE (EACH)	OX10-3830
BUCKET STAND - OXFORD STD S/S W/ NYLON BASE (EACH)	OX10-3840
BUCKET STAND - OXFORD KICK BASE S/S (EACH)	OX10-3820

Operating Theatre Carts

Surgeon's Step Stool and Wall Mounted Shoe Racks

The Surgeon's Step Stool is made from stainless steel tubular frame with non-slip aluminium tread surface. Dimensions are customisable.

The Wall Mounted Shoe Racks are made from stainless steel and available in 12 or 18 pair sizes.



Description	Code
FOOT STOOL - OXFORD SURGEON'S SINGLE #40 X 30 X 20 (EACH)	OX10-3850
FOOT STOOL - OXFORD SURGEON'S SINGLE #70 X 30 X 20CM (EACH)	OX10-3860
FOOT STOOL - OXFORD SURGEON'S DOUBLE #70 X 30 X 40CM (EACH)	OX10-3870
SHOE RACK - OXFORD 12 PAIR S/S #69 X 106CM (EACH)	OX40-1090
SHOE RACK - OXFORD 18 PAIR S/S #69 X 142CM (EACH)	OX40-1100



Db Health has a extensive range of Surgical products and supplies. Our product range cover Sterilisable Hollowware including Kidney Dishes to Theatre bowls, AMS Operating Surgical Equipment (scalpels, forceps, pincers, scissors) and Diagnostic Instruments from the leading manufacturers.

Skin and Laboratory Markers



Disposable and Reusable instrument



Surgical Headlights



WarwickWare Hollowware

