

# CSSD Tables and Carts

## CSSD Tables

CSSD Tables are manufactured from 304 grade stainless steel to a satin finish, each leg has a nylon cap and one is height adjustable.

Options available include:  
Wheels



## CSSD Prep Table

CSSD Prep Tables are manufactured from 304 grade stainless steel to a satin finish, each leg has a nylon cap and one is height adjustable, bottom shelf and top section with top shelf and rail for reels of sterilisation pouches.

Options available include:  
Wheels



Description	Code
CSSD TABLE - OXFORD S/S #110 X 60 X 86CM (EACH)	OX18-2211
CSSD TABLE - OXFORD S/S #170 X 60 X 86CM (EACH)	OX18-2217
CSSD TABLE - OXFORD S/S #230 X 60 X 86CM (EACH)	OX18-2223
CSSD TABLE - OXFORD S/S #110 X 75 X 86CM (EACH)	OX18-2311
CSSD TABLE - OXFORD S/S #170 X 75 X 86CM (EACH)	OX18-2317
CSSD TABLE - OXFORD S/S #230 X 75 X 86CM (EACH)	OX18-2323
CSSD PREP TABLE - OXFORD S/S #110 X 60 X 86CM (EACH)	OX19-5316

## Premium Case Cart

- Strong, fully welded construction
- 3x slide out wire shelves
- 2x 270° fold back doors with hold back and closed stay latches
- Twin 270° fold back doors
- Trolley finishes: stainless steel with stain finish
- Manufactured from 304 stainless steel
- Castors standard: 150mm stainless steel hi-temp, 2x plain and 2x fixed



## Bedpan and Bottle Rack

The Bedpan and Bottle rack is a Stainless Steel Rack designed to fit bedpans and urine bottles. Can be wall mounted or seat on a desktop.

Description	Code
BEDPAN & BOTTLE RACK - OXFORD S/S 4 UNIT #54 X 26 X 52CM (EACH)	OX30-1020
BEDPAN & BOTTLE RACK - OXFORD S/S 6 UNIT #78 X 26 X 52CM (EACH)	OX30-1030
CASE CART - PREMIUM #87 X 680 X 1180CM (EACH)	OX20-0014

