

# Maternity Bassinet Trolleys and Cots

## Baby Bassinet Trolley

The Baby Bassinet Trolley is available in 2 sizes Single or Twins.

- Mattress base height: 87.5 cm.
- Wheels: 7.5cm standard with 2 locking.
- Clear acrylic bassinet tub and 2 storage trays.
- Trendelenberg and Reverse Trendelenberg positioning.
- 5cm Mediflex covered mattress.
- Available in Powder Coated and Stainless Steel



Description	Code
BASSINET TRLY - OXFORD SINGLE S/S #80 X 42 X 106.5CM (EACH)	OX20-3815
BASSINET TRLY - OXFORD SINGLE P/C #80 X 42 X 106.5CM (EACH)	OX20-3812
BASSINET TRLY - OXFORD TWINS S/S #80 X 64 X 106.5CM (EACH)	OX20-3842
BASSINET TRLY - OXFORD TWINS P/C #80 X 42 X 106.5CM (EACH)	OX20-3840

Description	Code
BASSINETT MATTRESS - OXFORD SINGLE MEDIFLEX #66 X 32 X 5CM (EACH)	OX20-4100
BASSINET ACRYLIC TUB - OXFORD SINGLE #71 X 36 X 19-29CM (EACH)	OX20-4000
BASSINET TRAY - OXFORD WHITE PLASTIC #36 X 32 X 8CM (EACH)	OX20-4200
BASSINET ACRYLIC TUB - OXFORD TWINS (EACH)	OX20-4001
BASSINETT MATTRESS - OXFORD TWINS MEDIFLEX #66 X 50 X 5CM (EACH)	OX20-4101

## Hi-Lo Baby Bath

The Hi-Lo Baby Bath Trolley is designed specifically for the purpose of giving baby baths without the need to empty heavy bathtubs and to be able to adjust the bath to a comfortable height. Electronically adjustable height from 63cm to 93cm

- Wheels: 10 cm standard with 2 locking.
- Clear acrylic bassinet single tub with drainage and plug.
- Controlled via Foot pedal control.
- Available in Powder Coated Only.



## Single Bassinette Bath Trolley

The Single Bassinette Bath Trolleys provides a cost effective alternative to Baby Baths without height adjustment.

- Bassinette features: clear acrylic twin tub with plug/waste drain hose and valve
- Tray feature: removable 2 x 390 x 350 x 80mm deep, white acrylic
- Tub height: 820mm high
- Manufactured from 304 stainless steel tube
- Castors standard: 75mm single rubber with 2x braking



Description	Code
BABY BATH - OXFORD HI-LO (EACH)	OX20-3860
BASSINET BATH TRLY - OXFORD SINGLE S/S W/ TUBING FOR BATH (EACH)	OX20-3850

# Maternity Bassinet Trolleys and Cots

## Mobile Cots

The range of Mobile Cot are designed and manufactured to meet AS2130.

- Light weight side rails with squeeze grip handle for easy single handed operation.
- On each corner Socket with adjustable knob for either IV Poles / Detachable C Size Oxygen Bottle Holder can placed as required (IV Poles and Oxygen Bottle Holder is optional).
- Bottom shelf for additional storage.
- 10cm wall bumpers on each corners.
- Trendelenberg Angle Gauge
- Wheels: 12.5cm standard with 2 locking (diagonally positioned - expect for King Cot).
- Mediflex covered mattress—7.5cm high

	Trendelenberg Adjust	Mattress Height	End Panels	Extras
Neo-Natal Cot	Both Ends Manual	86cm Height	Fixed end panels	
Premier Cot	Both Ends Manual	86cm Height	Fixed end panels	
King Cot	Both Ends Electrical	95.5cm Height	Fixed end panels	Central Locking Castors
Quattro Cot	Both Ends Manual	86cm Height	Removable end panels	
King Quattro Cot	Both Ends Electrical	95.5cm Height	Removable end panels	Central Locking Castors



OX20-3802



OX20-3811



OX20-3800



OX20-3804

## Cot Accessories

Include Telescopic IV Poles 76cm long with Ramshead and Detachable C Size Oxygen Bottle Holders.



OX13-3012


Description	Code
MOBILE COT - OXFORD NEO-NATAL P/C #102 X 68.5 X 145CM (EACH)	OX20-3802
MOBILE COT - OXFORD PREMIER P/C #143 X 79 X 145CM (EACH)	OX20-3800
MOBILE COT - OXFORD KING ELECTRIC P/C #143 X 79 X 152.5CM (EACH)	OX20-3804
MOBILE COT - OXFORD QUATTRO P/C #143 X 79 X 145CM (EACH)	OX20-3811
MOBILE COT - OXFORD KING ELECTRIC QUATTRO P/C #143 X 79 X 152.5CM (EACH)	OX20-3818

# Colour Options for Medical Carts

## Oxford Powdercoat Colours

Below are the colours for Oxford Products along with white gloss. Please colours are approximate representations of the final powdercoating colour.

When selecting colour please advise drawer colour and body colour.



Standard Powdercoat Colours



Non-Standard Powdercoat Colours (additional cost may apply)

## Flexicart Powdercoat Colours

Below are the standard colours for Flexicart Medication Trolleys along with White and Signal Red.

