

# Proctological Devices

## Heine Instrument Head for UniSpec Tubes

- Heine Instrument Head for UniSpec disposable tubes provide annular fiber optic illumination.
- Bright, reflex-free illumination.
- Magnifying lens 1.5X which produces a high resolution image.
- Hinged viewing window with insufflation port.
- Stainless steel construction.
- Can be powered by Lamp handle with plug-in transformer, Illumination adaptor with 3.5V Beta 4 handle or Projector with F.O. cable.



Product	Product Code
Instrument Head for Disposable UniSpec Tubes #viewing window and swivel lens w/ 1.5X Magnification (Each)	HE0318100
Filter for Heine UniSpec Proctology Instrument Head (Each)	HEFILTER
Lamp Handle 6V w/ switch (Each)	HX0499607
Plug-in Transformer E7 (Each)	HX9516104
Heine Illumination Adaptor #use with 2.5V Handle and Instrument Head (Each)	HX0298002
Heine Beta 4 NT Rechargeable Handle with Li-ion Battery (Each)	HX0799396
NT4 Table Recharger #2 Handles (Each)	HX0299494

## Heine UniSpec Disposable Tubes

Heine UniSpec disposable tubes are designed to easily used with the Instrument Head found in the previous listing.

- Tubes are disposable meaning no cross infection.
- No cleaning required.
- Moulded in patient friendly thermoplastic.
- Light transmission through tube wall with distular light exit.
- Graduated in cm on the inside of the tube.



HE0319925



HE0319825



HE0318825

Description	Product Code
UniSpec Disposable Anoscopes 85 (L) x 20mm (W) (B/25)	HE0319925
UniSpec Disposable Proctoscopes 130 (L) x 20mm (W) (B/25)	HE0319825
UniSpec Disposable Sigmoidoscopes 250 (L) x 20mm (W) (B/25)	HE0318825

## Heine Proctology Accessories

Heine provides a number of accessories that are designed for use in the proctology area. These include:

- Rubber Insufflation Bulb - improved vision
- Sponge Holder - 400mm long and constructed from Stainless Steel making this item easy to clean, disinfect and sterilise.
- Suction Tube - 400mm long and constructed from Stainless Steel making this item easy to clean, disinfect and sterilise.



Product Description	Product Code
Heine Rubber Inflation Bulb w/ Hygiene Filter (Each)	HE0018105
Heine Stainless Steel Sponge Holder 400mm (Each)	HE0018906
Heine Stainless Steel Suction Tube 400mm (Each)	HE0018907