

Dermatoscopes

Heine Delta 20 T Dermatoscope

For dermatoscopic examination for early diagnostic of malignant melanoma and examination of non-melanocytic lesions, basal cell carcinoma and dermatofibroma.

- Precision Heine Optics. High quality, achromatic optical system provides unmatched sharpness and resolution. Focusing system provides for 10x to 16x magnification for recognition of even the smallest details and structures.
- Adjustable eyepiece, Individual focusing, range of correction approx. -6 to +6.
- LED in HQ technology. Red is Red, blue is blue. Colour tem. 5000 K, CRI >80 on maximum scale of 100.
- Optional Small contact plate available - Special small contact plate for difficult to access lesions
- Digital documentation capabilities connect to major brand SLR digital cameras with Heine digital lens system,. Specially designed lenses integrate the camera and Delta 20 T optics together to provide the highest quality imaging possible.
- Conveniently switch between polarised and non-polarised views at the press of a button this allows recognition of "Blink Sign".
- A range of accessories are available for the Delta 20 T include SLR Photo adaptor, contact plates with and without scale, carry case and Dermatoscopy oil.



Beta 4 USB Handle with Plug-in Power Supply

Product Code	Description
HK26210118	DERMATOSCOPE SET - HEINE DELTA 20T 2.5V W/ BETA HANDLE (SET)
HK26228388	DERMATOSCOPE SET - HEINE DELTA 20T 3.5V W/ BETA 4 USB HANDLE AND PLUG-IN POWER SUPPLY (SET)
HK26224420	DERMATOSCOPE SET - HEINE DELTA 20T 3.5V W/ BETA 4 NT HANDLE AND NT 4 TABLE CHARGER (SET)

Sets include Delta 20 T Dermatoscope with contact plate 23mm diameter with scale, 10ml dermatoscopy oil, dermatoscopy compendium and hard case. Dermatoscope Heads are available individually as well.



NT 4 Table Charger with Beta 4 NT Handle



Delta 20 T Instrument Head and Handle

Dermatoscopes



Photo Accessory Set

Photo Accessory Sets are available that allow SLR Digital Cameras to take high quality images. These sets are available for Canon, Nikon and Olympus SLR Cameras. These sets include SLR Photo Adapter, Adaptor Code, Distance Ring and Beta Belt Clip.



Small Contact Plate



Dermatology Oil

Product Code	Description
HK0034216	ACCESSORY - HEINE DERMATOSCOPY CONTACT PLATE 23MM DIA W/ SCALE (EACH)
HK0034217	ACCESSORY - HEINE DERMATOSCOPY CONTACT PLATE 23MM DIA W/ SCALE (EACH)
HK0034205	ACCESSORY - HEINE DERMATOSCOPY SMALL CONTACT PLATE 8MM DIA (EACH)
HK0034005	ACCESSORY - HEINE DERMATOSCOPY OIL 10ML (B/6)
HX0099231	ACCESSORY - HEINE ADAPTOR CORD 1.3M FOR INSTRUMENT AND HANDLE (EACH)
HX0099008	ACCESSORY—HEINE BELT CLIP FOR BETA HANDLE (EACH)

Heine NC1 Dermatoscope

Speed and diagnostic accuracy for the daily demands of a busy practice.

- Polarised illumination for detailed views of the crystalline vascular structures without any need for contact liquids
- Precision Heine optics with 6x magnification - High quality aspherical optical system provides unmatched sharpness and resolution.
- LED in HQ technology, red is red blue is blue. Colour temperature : 5000K colour rendering index (CRI) > 80, special index for red colours (R9) >35 on a maximum scale of 100
- Optional magnification plate with additional magnification. Adding the optional contact plate transforms the NC1 to a traditional contact dermatoscope with increased magnification and scale.



NC1 Dermatoscope with Mini 3000 Handle and Contact Plate.

Product Code	Description
HD0878132	DERMATOSCOPE SET - HEINE NC1 W/ MINI 3000 BATTERY HANDLE (EACH)
HD0878134	DERMATOSCOPE SET - HEINE NC1 W/ MINI 3000 BATTERY HANDLE AND CONTACT PLATE (EACH)
HD0078115	ACCESSORY—HEINE NC1 CONTACT PLATE W/ SCALE

Sets selected are for standard handles which use 2 x AA batteries. Also available are rechargeable sets and instrument head.