

Silent Knight Pill Crusher and Pill Cutter

Silent Knight Pill Crushers

Offers healthcare facilities many helpful features, including the use of OxiTitan™, an antimicrobial coating, to help fight against touch contamination infections.

- Is distinguished for its smooth, quiet operation assuring that residents will not be disturbed, day or night.
- The Silent Knight Pill Crushers effectively crushes multiple tablets at once.
- Patented single use no mess pouches designed to eliminate cross contamination.
- Portable & sturdy, designed to fit on medication trolleys or in our Carstens wall units
- Simple and easy to use, no special training required.
- Used in most top hospitals and over 500 aged care and rehabilitation institutions across Australia.

NS999a - Silent Knight



Silent Knight Pill Crusher Pouches

- Purpose designed pouches for tablet during crushing for use with Silent Knight Pill Crushers.
- Pouch ensure no-mess crushing.
- Eliminates cross-contamination.
- Tear resistant.
- You can directly add fluids to the crushed tablets inside the pouch for oral or tube feeding, eliminating touch contamination spillage and loss



NS9997

Ezy Dose Original Pill Cutter

Simple to operate mechanical pill cutter - easily cuts pills in half for easier swallowing.

- Only to be used on 1 tablet at a time.
- Must be washed between uses to avoid cross contamination.
- Not to be used with Enteric-coated tablets.



AMC67015

Product Code	Description	Individual Unit Qty	Carton Qty
NS9999a	PILL CRUSHER - SILENT KNIGHT	EACH	P.O.A.
NS9997	PILL CRUSHER POUCH - DELUXE #ROUND BASE	B/1000	8 boxes of 1000
AMC67015	ORIGINAL TABLET CUTTER EZY DOSE	EACH	72 units