

# Aspire Wheelchair—Assist 2 Deluxe Vinyl

The Aspire ASSIST 2 Deluxe Vinyl wheelchair has all the features of the ASSIST 2, but with durable deluxe vinyl, suitable for the hospital & clinic environment. The Aspire ASSIST 2 Deluxe Vinyl has been carefully designed with fully welded seams and concealed mounting of upholstery to ensure that rigorous infection control standards can be met.

## FOLDING DESIGN

The centre cross-brace design allows for the Aspire ASSIST 2 Deluxe Vinyl to be folded up for unoccupied transport or storage.

## IMPRESSIVE SAFE WORKING LOAD (S.W.L)

140kg: Seat widths 350-500mm  
180kg (HD): Seat widths 570 - 610mm

## DURABLE TRANSPORT SOLUTION

The Aspire ASSIST 2 Deluxe Vinyl is a versatile transport wheelchair. 24 inch quick release wheels with solid, maintenance-free tyres provide on-call reliability and are suitable across common areas, bedroom and bathroom.

## ROBUST & RELIABLE

With a strong steel frame, removable swing-away legrests and quick release self-propelling rear wheels, the Aspire ASSIST 2 Deluxe Vinyl is a practical yet robust patient transport solution.



Assist 2 is registered in the ARTG, AUST R ARTG 291670



Product Code	Description
AD449340	Wheelchair—Aspire Assist 2 Deluxe Vinyl 35cm Wide 140kg (Each)
AD449350	Wheelchair—Aspire Assist 2 Deluxe Vinyl 40cm Wide 140kg (Each)
AD449360	Wheelchair—Aspire Assist 2 Deluxe Vinyl 45cm Wide 140kg (Each)
AD449370	Wheelchair—Aspire Assist 2 Deluxe Vinyl 50cm Wide 140kg (Each)
AD449380	Wheelchair—Aspire Assist 2 Deluxe Vinyl 57cm Wide 180kg #HD (Each)
AD449390	Wheelchair—Aspire Assist 2 Deluxe Vinyl 61cm Wide 180kg #HD (Each)



W: [www.dbhealth.com.au](http://www.dbhealth.com.au)  
E: [customerservice@denyer.com.au](mailto:customerservice@denyer.com.au)  
P: (02) 9638 2100  
ABN: 37 000 459 174

DBH1802

Suppliers of quality products to the Healthcare Industry est 1889. Australian Owned & Operated.