

ASPIRE CLASSIC BEDSIDE COMMODE












MODELS

1	Removable Backrest – Enables easy storage or transportation
2	Raised Padded Armrests – Facilitate easier transfers and offer improved comfort
3	Easy clean vinyl surfaces
4	Height adjustable legs with easy push button adjustment and non-slip rubber tips for stability

5	Removable bowl with lid and carry handle
6	Removable, Padded Seat – with lift-up toilet seat and fixed bowl carrier beneath
7	Lightweight Design – Lightweight yet sturdy steel design for easy manoeuvrability (5.5kg)
8	Three-sided Shroud – Discreetly conceals bowl
9	Powder Coated Steel Frame – Strong, durable and easy to clean

Slate Vinyl

1
Year
Warranty

SEAT DEPTH	OVERALL DEPTH	SEAT WIDTH	OVERALL WIDTH	SEAT HEIGHT	BACKREST HEIGHT	OVERALL HEIGHT	UNIT WEIGHT	S.W.L
								
440mm	600mm	450mm	560mm	475-615mm	450mm	900-1040mm	5.5kg	125kg

AD045055 Classic Bedside Commode - Champagne Vinyl (Each)



W: www.dbhealth.com.au
E: customerservice@denyer.com.au
P: (02) 9638 2100
ABN: 37 000 459 174

Suppliers of quality products to the Healthcare Industry est 1889. Australian Owned & Operated.

ASPIRE DELUXE BEDSIDE COMMODE



Champagne Vinyl

1	Generously padded seat lifts and tilts neatly out of the way to reveal a lift-up toilet seat and fixed bowl carrier beneath
2	Padded armrests are raised to provide an easily reachable support point for lowering onto the seat
3	Removable bowl with lid and carry handle
4	Contoured, deluxe padded and slightly reclined backrest helps maintain posture while supporting an open and comfortable hip position

5	Height adjustable legs with easy push button adjustment and non-slip tips for stability
6	Easy clean vinyl surfaces
7	Three-Sided Shroud – Discreetly conceals bowl
8	Powder Coated Steel Frame – Strong, durable and easy to clean

1
Year
Warranty

SEAT DEPTH	OVERALL DEPTH	SEAT WIDTH	OVERALL WIDTH	SEAT HEIGHT	BACKREST HEIGHT	OVERALL HEIGHT	UNIT WEIGHT	S.W.W.L
470mm	600mm	460mm	600mm	410-565mm	440mm	865-1015mm	12kg	160kg

AD045115 Deluxe Bedside Commode - Champagne Vinyl (Each)



W: www.dbhealth.com.au
E: customerservice@denyer.com.au
P: (02) 9638 2100
ABN: 37 000 459 174

Suppliers of quality products to the Healthcare Industry est 1889. Australian Owned & Operated.