

Aspire Shower Commodes

Aspire Shower Commodes—Weight Bearing Footrest

Aspire Attendant Propelled 460mm and 530mm Weight Bearing Platform Shower Commode with 200kg max user weight.

- Assists both users and carer with toileting and showering
- Swing back arms for easier transfer
- Deeper access trough for easier hygiene
- Easy slide weight bearing platform
- Safety arms as standard
- Soft PU seat with 70mm of travel (forward & back) - **Seat sold separately**
- 4 x lockable stainless steel casters as standard
- Warranty - Corrosion: 5 Years and Component: 1 Year



Aspire Shower Commodes—Heavy Duty Weight Bearing Footrest

Aspire Attendant Propelled 600mm Weight Bearing Platform Shower Commode with 400kg Max user weight

- Assists users and carers alike for toileting and showering
- Removable Safety Arms as Standard
- Deeper Access Trough for Easier Hygiene
- Easy-Slide 400kg Weight Bearing Platform
- Height Adjustable Armrests
- Soft PU Seat with 70mm of Travel (Fwd & Back) - **Seat sold separately**
- 4 x Stainless Steel Lockable Castors as Standard
- Warranty - Corrosion: 5 Years and Component: 1 Year



Code	Overall Height	Overall Width	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height
AD066000	995mm	535mm	460mm	305mm	225mm	590mm	560mm
AD066050	995mm	605mm	530mm	305mm	225mm	590mm	560mm
AD061400	995mm	725mm	600mm	300mm	290mm	590mm	560mm

Product Code	Description
AD066000	SHOWER COMMUNE - ASPIRE S/S 46CM #200KG W/ WB PLATFORM (EACH)
AD066050	SHOWER COMMUNE - ASPIRE S/S 53CM #200KG W/ WB PLATFORM (EACH)
AD061400	SHOWER COMMUNE - HEAVY DUTY ASPIRE S/S 60CM #400KG W/ WB PLATFORM (EACH)

Aspire Shower Commodes

Aspire Shower Commodes—Swing Away Footrest

Aspire Attendant Propelled 460mm, 530mm and 600mm Swingaway Footrest Shower Commode with 200kg max user weight

- Assists both users and carer with toileting and showering
- Swing back arms for easier transfer
- Deeper access trough for easier hygiene
- Swing Away Footrests
- Safety arms as standard
- Soft PU seat with 70mm of travel (forward & back) - **Seat sold separately**
- 4 x lockable stainless steel casters as standard
- Warranty - Corrosion: 5 Years and Component: 1 Year



Code	Overall Height	Overall Width	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height
AD065900	995mm	535mm	460mm	305mm	225mm	590mm	560mm
AD066060	995mm	605mm	530mm	305mm	225mm	590mm	560mm
AD061300	995mm	675mm	600mm	305mm	225mm	590mm	560mm

Aspire Shower Commodes—Tilt in Space

Aspire Attendant Propelled 460mm Tilt in Space Swingaway Footrest Shower Commode with 250kg max user weight.

- Backrest may be configured as full length or 3/4 length with the top 1/4 section removable
- Open back straps provide greater access for bathing needs
- 30° of backrest angle adjustment in 10° increments for effortless repositioning and maximum user comfort
- Safety arms prevent falls and swing back out of the way for transfers
- 35° of gas-assisted tilt-in-space functionality minimises manual handling efforts during the showering process
- Tool-free height adjustment to 4 settings allow the commode to be easily adapted to varying needs
- Height adjustable, swing away footrests for optimal positioning and ease of transfers
- 4 lockable stainless steel castors maximise safety
- Warranty - Corrosion: 5 Years and Component: 1 Year



Code	Overall Height	Overall Width	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height
AD066030	1315-1485mm	560mm	460mm	900mm	275mm	460-630mm	460-620mm

Product Code	Description
AD065900	SHOWER COMMUNE - ASPIRE S/S 46CM #200KG W/ S/AYAW FOOTREST (EACH)
AD066060	SHOWER COMMUNE - ASPIRE S/S 53CM #200KG W/ S/AYAW FOOTREST (EACH)
AD061300	SHOWER COMMUNE - ASPIRE S/S 60CM #200KG W/ S/AYAW FOOTREST (EACH)
AD066030	SHOWER COMMUNE - ASPIRE TILT IN SPACE S/S 46CM #250KG W/ S/AYAW FOOTRESTS (EACH)

Aspire Shower Commodes

Aspire Shower Commodes—Seats


PU seamless seats offer excellent durability and easy cleaning. Assists both users and carer with toileting and showering.

Aspire Shower Commodes—Commode and Bedpan Holders

Compatible with Aspire Shower Commodes.

Carries a pan below the commode seat

Carrier slides easily into the shower commode seat



Aspire Shower Commodes—Options



Oxygen Bottle Holders, Headrests, Elevating Legrests, Stump Supports, IV Poles are also available.

Product Code	Description
AD066600	OPTION - OPEN SEAT FOR ASPIRE SHOWER COMMODE 46CM (EACH)
AD066700	OPTION - CLOSEDSEAT FOR ASPIRE SHOWER COMMODE 46CM (EACH)
AD066750	OPTION - OPEN SEAT FOR ASPIRE SHOWER COMMODE 53CM (EACH)
AD066800	OPTION - OPEN SEAT FOR ASPIRE SHOWER COMMODE 60CM (EACH)
AD066400	OPTION – ASPIRE BOWL & PAN CARRIER 46/53/60CM WIDE (EACH)
AD081600	COMMODOE BOWL – ASPIRE W/HDL & LID (EACH)