

---

## Lab and Skin Markers

---

### Securline Lab Marker

Securline Laboratory Markers are designed for marking on plastic, metal, and porcelain labware. The unique ink has been formulated for writing on labware even when the surface is damp. These markers are available with either Black, Blue, Green or Red ink.



Product Code	Description
PD1400	MARKER - SECURLINE LAB MARKER BLACK FINE (PKT/10)
PD1400-13	MARKER - SECURLINE LAB MARKER BLUE FINE (PKT/10)
PD1400-15	MARKER - SECURLINE LAB MARKER GREEN FINE (PKT/10)
PD1400-16	MARKER - SECURLINE LAB MARKER RED FINE (PKT/10)

### Securline Surgical Skin Markers

Designed with a taper tip that allows for both broad and fine lines with the use of a single marker that does not scratch the skin. This marker uses a gentian non-toxic violet ink that is resistant to most skin preparation solutions. These markers are supplied in sealed sterilised packets which include a 15cm rule.



Product Code	Description
PD1000	MARKER- SURGICAL SKIN STERILE W/ RULER #VIOLET (PKT/10)

### Viscot Dual Tip Sterile Markers

Designed with a permanent Ink marker tip and a surgical skin marker tip in the one marker.

Product Code	Description
PDSTER-DPR	MARKER - PDC DUAL TIP SKIN + PERMANENT MARKER STERILE (B/100)



---

## Custom Sterile Label Packets with Markers

---

Custom Sterile Label Packs are available to order with the choice of sterile labels you wish to have included . Sterile Markers can also be included in the same pack.