


## BetterShield Face Shield

BetterShield is a Full Face Shield that is designed to provide improved protection even on the sides and bottom of the face. The BetterShield provides optimum coverage against airborne pathogens, and provides the widest range of coverage for the healthcare worker protecting them against high level diseases.

The BetterShield can be used together with an N95 mask and glasses or goggles if required.

True Anti-fog film, distortion-free optics, secure fit with one size fits all, easy adjustability, airflow & ventilation, with a 3D curved design.



### Better how?

- **Better Price:** Premium performance shouldn't come with a premium price tag. That's why with BetterShield, your purchasing department will see it's better too.
- **Better Shielding:** Innovative and patent-pending, the unique 3D CURVED DESIGN provides unmatched fluid protection, even on sides and bottom of the face.
- **Better Fit:** One-size-fits-all, extensive design and testing has yielded a premium face shield that securely fits the widest range of users.
- **Better Liked:** Multi-site, multi-ward hospital testing revealed that health care workers prefer BetterShield over the nearest competitor.
- **Better Storage:** BetterShield is designed for ultra-compact stacking.

### Even Better Now since the start of the pandemic

- Improved thickness – moved from a 1.0mm thickness to a 1.2mm thickness - still light and flexible
- The application of the antifog is part of the material of the shield – so the shields are VERY clear, and fog much less!
- Moved from glue tabs for the elastic, to stapled elastic for reinforced strength
- Meets the following international standards : EN 166:2001 & ANSI/ISEA Z87.1
- Lastly, although the BetterShield has been sold as single use disposable, with all these improvements, we find the useful life of the BetterShield to be longer and reusable in some instances.

Product Code	Description
MXFS01	BETTERSHIELD DISPOSABLE FACE SHIELD (CTN/100)