

## SureSafe® Impervious Isolation Gown

SureSafe Impervious Isolation Gowns are non-sterile, latex-free gowns made from laminated impervious fabric that is suitable for heavy fluid applications. SureSafe Impervious Isolation Gowns feature knit cuffs, tie-neck closure and tie waste.

### Features

- Laminated impervious fabric
- Suitable for heavy fluid applications
- Tie neck closure
- Knit cuffs
- Tie waist
- Non sterile
- Latex-free
- Size - 125 x 150cm
- Material: 35gsm Polypropylene laminated with Polyethylene
- Recyclable



Product Code	Description
SA64002Y	GOWNS - SURESAFE IMPERVIOUS WATERPROOF ISOLATION GOWNS UNIVERSAL YELLOW W/ KNITTED CUFF (C/50)

Other PPE available from DB Health includes: Sterile and non sterile gloves, Better-Shield FaceShield, hand wash and sanitisers, surface wipes and disinfectants, and other PPE and cleaning products.



Suppliers of quality products to the Healthcare Industry est 1889.  
Australian Owned & Operated.

Web [dbhealth.com.au](http://dbhealth.com.au)

Email [customerservice@denyer.com.au](mailto:customerservice@denyer.com.au)

Freecall: 1800 451 446

ABN: 37 000 459 174

DBH2201