

Surgical Instruments, Electro and Cryosurgical Devices

Product range for Surgical Instruments, Electrosurgical and Cryosurgical Devices includes: High Frequency Desiccators, Electrosurgical Generator, Disposable Electrodes, Electrosurgical pencils, Nylon Coated Vaginal Speculum, Bovie Smoke Shark Evacuator, Cautery Pen, Cryospray and Cryopen, Disposable and Reusable Surgical Instruments



Web: dbhealth.com.au

Email: customerservice@denyer.com.au

Freecall: 1800 451 446

ABN: 37 000 459 174

Suppliers of quality products to the Healthcare Industry Est 1889.
Australian Owned & Operated.

Bovie Derm 102 High Frequency Desiccator

The Derm 102 is an economical, 10 watt high frequency desiccator, thoughtfully designed to give practitioners the ability to perform in-office, minor skin procedures such as the removal of skin tags, benign lesions and premalignant lesions. This Bovie unit allow practitioners to provide better patient care in the office, at a sensible price, saving the patient time and money, while increasing their practice revenue.

- Monopolar mode
- Output power adjusts by 1/10 watt increments up to 10 watts
- Ground reference unit means no grounding pad is required
- Large illuminated digital displays
- Made in USA.
- 4 Year warranty

Standard Inclusions:

- Wall mount kit
- 1 Power control pencil
- 2 Hand piece sheaths
- 21 assorted dermal tips:
- Users guide

Description	Product Code
ELETROSURGICAL GENERATOR—BOVIE DERM101 (EACH)	DAD101
OPTION—DESKTOP STAND FOR BOVIE E/S GENERATORS S/S (EACH)	DAA813
OPTION—MOBILE STAND FOR BOVIE E/S GENERATORS P/C STEEL BASE #ADJ 100-150CM (EACH)	DA1580ARN



Bovie A942 High Frequency Desiccators

Designed for years of trouble-free use, the Bovie Derm 942 is engineered to be the most reliable and durable desiccator available today. The unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat. Post-procedure and long-term cleanup are simplified by its non-absorbing, non-staining, high impact plastic case.

- Automatic high and low power modes
- Output power adjusts by 1/10 watt increments up to 10 watts, then by 1 watt increments up to 40 watts
- Superior quality power control pencils
- Electrodes are firmly gripped, with no irritating swivel
- Lower operating costs. Aaron replacement parts and accessories are on average 30% less expensive than competitors
- No patient plate required
- Made in the USA
- 4 Year warranty

Standard Inclusions:

- Model A940 HF Desiccator unit
- 1 Power control pencil (A901) + 1 extra **BONUS** pencil
- 50 sharp and 50 blunt disposable dermal tips
- Wall mount and power cable
- Kit manual

Description	Product Code
ELETROSURGICAL GENERATOR—BOVIE A940 (EACH)	DAA942
OPTION—DESKTOP STAND FOR BOVIE E/S GENERATORS S/S (EACH)	DAA813
OPTION—MOBILE STAND FOR BOVIE E/S GENERATORS P/C STEEL BASE #ADJ 100-150CM (EACH)	DA1580ARN
ACCESSORY—BOVIE FOOT SWITCH #BIPOLAR (EACH)	DAA803



Aaron Specialist Pro Electrosurgical Generator

The Bovie Specialist Pro is the electrosurgery solution for in demand physician offices and specialty facilities. Designed to operate like an operating room generator in a compact lower watt package. The multipurpose unit is ideal for monopolar and bipolar office surgical procedures.

- Five output modes - cut, blend, coag, fulguration & bipolar
- Large illuminated digital displays
- Neutral electrode monitoring safety system
- Discreet outputs (FCFS™ - First Come First Served)
- Digital error detection
- With international fittings

Standard Inclusions:

- Unit & power cord
- 1 Sterile disposable pencil (ESP1)
- 3 Sterile disposable electrodes (1x ES01 blade, 1x ES02 needle, 1x ES20 ball)
- 2 Reusable grounding cords (A1252C)
- Five disposable split ground pads (ESRE-1)
- User's guide

Description	Product Code
ELETROSURGICAL GENERATOR—BOVIE SPECIALIST PRO #A1250S (EACH)	DAA1250S
OPTION—DESKTOP STAND FOR BOVIE E/S GENERATORS S/S (EACH)	DAA813
OPTION—MOBILE STAND FOR BOVIE E/S GENERATORS P/C STEEL BASE #ADJ 100-150CM (EACH)	DA1580ARN
ACCESSORY—BOVIE FOOT SWITCH #SUITS A1250S (EACH)	DAB1253B
ACCESSORY—ADAPTOR BOVIE FOOT SWITCH #SUITS A1250S (EACH)	DAA1255A
SPLIT RETURN ELECTRODE—ADULT (B/50)	DAWR021



Disposable Dermal Tip Electrodes

All disposable electrodes utilize the standard 2.3 mm stainless steel shaft. All molded electrodes feature safety grip insulators combining patient and user safety with easy insertion into and removal from the surgical pencil. All disposable electrodes are manufactured to the highest standards. They come individually packaged sterile or non-sterile box packed.

Description	Product Code
DISPOSABLE ELECTRODES - AARON N/S SHARP DERMAL TIP (B/100)	DAA804
DISPOSABLE ELECTRODES - AARON N/S BLUNT DERMAL TIP (B/100)	DAA806
DISPOSABLE ELECTRODES - AARON STERILE SHARP DERMAL TIP (B/100)	DAA805
DISPOSABLE ELECTRODES - AARON STERILE BLUNT DERMAL TIP (B/100)	DAA807

Disposable Non-Stick Electrodes

The build-up of carbonised material (eschar) on the electrosurgery tip during surgery is radically reduced with the Prima non-stick p.t.f.e. coated electrodes. The eschar layer prohibits the electrosurgery current being passed through the electrode, necessitating higher output settings which can cause damage to the peripheral tissue. Nonstick electrodes are designed to focus the energy on the target tissue, making surgical time more productive.

Benefits of using Non-Stick Electrodes:

- They are frustration free, allowing you to focus on the surgery
- Quicker procedure times
- Use lower power settings
- Reduction of peripheral electro surgery damage
- Less smoke
- Easier cleaning of the tip

Disposable Prima Non-Stick Electrodes

Shielded tip electrodes offer a greater level of safety when working in cavities and areas of restricted vision.

Description	Product Code
DISPOSABLE ELECTRODES - PRIMA NONSTICK 63MM BLADE STERILE (B/10)	DAPMS2920
DISPOSABLE ELECTRODES - PRIMA NONSTICK 63MM BLADE STERILE #SHIELDED TIP (B/10)	DAPMS2950
DISPOSABLE ELECTRODES - PRIMA NONSTICK 70MM NEEDLE STERILE (B/10)	DAPMS2922
DISPOSABLE ELECTRODES - PRIMA NONSTICK 70MM NEEDLE STERILE (B/10)	DAPMS2972
DISPOSABLE ELECTRODES - PRIMA NONSTICK 5MM BALL X 130MM STERILE (B/10)	DAPMS2906L
DISPOSABLE ELECTRODES - PRIMA NONSTICK 4MM BALL X 50MM STERILE (B/10)	DAPMS290



Disposable Lletz Loop Electrodes

- Ultra-fine rigid tungsten wire for minimum eschar buildup
- Coloured for easy identification of shape and size
- Packed in protective packaging
- Box of 5
- Sterile Packed

Description	Product Code
DISPOSABLE LLETZ LOOP ELECTRODE—AQUA 10X10MM STERILE (B/5)	DAPMS1742
DISPOSABLE LLETZ LOOP ELECTRODE—BLUE 15X15MM STERILE (B/5)	DAPMS1743
DISPOSABLE LLETZ LOOP ELECTRODE—YELLOW 20X20MM STERILE (B/5)	DAPMS1744
DISPOSABLE LLETZ LOOP ELECTRODE—RED 20X15MM STERILE (B/5)	DAPMS1745
DISPOSABLE LLETZ LOOP ELECTRODE—ORANGE 15X8MM STERILE (B/5)	DAPMS1746
DISPOSABLE LLETZ LOOP ELECTRODE—PINK 15X12MM STERILE (B/5)	DAPMS1747
DISPOSABLE LLETZ LOOP ELECTRODE—BLUE 20X10MM STERILE (B/5)	DAPMS1748
DISPOSABLE LLETZ LOOP ELECTRODE—GREY 20X12MM STERILE (B/5)	DAPMS1751



Handpieces

Description	Product Code
ELECTROSURGICAL PENCIL - AARON POWER CONTROL REUSEABLE #3 BUTTON (EACH)	DAA902
ELECTROSURGICAL PENCIL - AARON POWER CONTROL REUSEABLE #2 BUTTON (EACH)	DAESPR
PENCIL SHEATH AARON DISPOSABLE (B/100)	DAA910
PENCIL SHEATH AARON DISPOSABLE #X-LONG (B/25)	DAA910ST

DAA902—Reusable 3 button power pencil for use with DERM101, A942 and A1250S—Autoclaveable at 121C.
DAESPR—Reusable 2 button and 3 prong connection power pencil for use with A1250S.

Pencil Sheath designed to keep your pencils clean without having to sterilise.

Bipolar Disposable Instruments and Cables

- No more cleaning, no more sterilising, no more fault checking, no more reinsulating or repairs, no more mismatched cables, no more wasted time!
- TGA compliant 2-pin bipolar plugs, with adaptors available for all electrosurgical units
- No more uncertainty. Every time you open a new sterile pack you know the forcep and cable will not fail.
- The feel and handling performance is the same as a reusable forcep.
- Cable fits Vallylab machines. Permanent reusable adapters available for Erbe/Conmed/Bovie/Martin/Other. These adapters do not have to be reprocessed, they live on the machine.

Description	Product Code
BIPOLAR CABLE—AARON FOR A942 & A1250s (EACH)	DAA827V
ELECTROSURGICAL PENCIL - AARON POWER CONTROL REUSEABLE #2 BUTTON (EACH)	DAESPR
PENCIL SHEATH AARON DISPOSABLE (B/100)	DAA910
PENCIL SHEATH AARON DISPOSABLE #X-LONG (B/25)	DAA910ST

Nylon Coated Vaginal Speculum

- Nylon coating resists abrasion and provides insulation
- Fully autoclaveable and easy to care for
- Non conductive for use during electrosurgery



Description	Product Code
VAGINAL SPECULUM—NYLON COATED SMALL GRAVES W/ DISP EVAC TUBE (EACH)	DA381606
VAGINAL SPECULUM—NYLON COATED MEDIUM GRAVES W/ DISP EVAC TUBE (EACH)	DA381608
VAGINAL SPECULUM—NYLON COATED LARGE GRAVES W/ DISP EVAC TUBE (EACH)	DA381610
STRAWS—ULTRASPEC #SUITABLE WITH VAGINAL SPECULUM FOR SMOKE EVAC (B/50)	DA920001

Bovie Smoke Shark Evacuator

Extremely effective at removing smoke, particles and odour at the surgical site. Protects clinicians and patients.

- A modern, light-weight, compact design allowing for simplified operation and storage
- Soft touch controls
- Packaged with both (DASF35) Extended-Life Filter and Pneumatic Footswitch.
- Extended-Life Filter offers up to 35 hours of performance
- Remaining filter life indicators
- Covered 3-port filter design accommodates a variety of tubing sizes [1/4"(6.4mm), 3/8"(9.5mm) and 7/8"(22mm)]
- Variable suction power and Quieter operation
- Operates in conjunction with optional Bovie Remote Activation Switch (SERS2)



Description	Product Code
SMOKE EVACUATOR—BOVIE SMOKE SHARK II (EACH)	DASE02220
SMOKE EVACUATOR FILTER—BOVIE SMOKE SHARK II #35 HOUR (EACH)	DASF35
TUBING—NON-STERILE 22 X 1830MM (B/24)	DA786T
REDUCER FITTING TUBE—22MM TO 6.4MM X 610MM L (B/10)	DASERF
ARMSTAND—BOVIE SMOKE SHARK II EVACUATOR (EACH)	DASEAS

Bovie Change-A-Tip Cautery

- A highly portable and easy to use cauterizing tool.
- To operate, depress the finger switch. Current flows through the burner tip.
- Reusable plastic handle with switch
- Batteries and disposable burner-tip can be replaced.
- Supplied with 1 non-sterile H101 fine tip and 2 x AA batteries
- Additional tips are sterile packed, 10 per box. Each tip has sterile drape to cover the power handle.



Description	Product Code
BATTERY CAUTERY UNIT - BOVIE CHANGE-A-TIP W/ TIP (EACH)	DAHIT1
CAUTERY DISPOSABLE TIP - HIGH TEMP FINE (B/10)	DAH101
CAUTERY DISPOSABLE TIP - HIGH TEMP LOOP (B/10)	DAH103
CAUTERY DISPOSABLE TIP - HIGH TEMP 50.8MM FLEX LOOP (B/10)	DAH105
CAUTERY DISPOSABLE TIP - HIGH TEMP 50.8MM FLEX FINE (B/10)	DAH106
CAUTERY DISPOSABLE TIP - HIGH TEMP CIRCULAR (B/10)	DAH111

Brymill Cry-Ac 0.5L Cryospray

Made in the USA and the number 1 choice of dermatologists worldwide, the Brymill cryosurgery device is unsurpassed for safety, versatility and precise control.

- 80 to 100 Freezes at an average of 10 seconds each
- 29cm tall, weighs 895g when full
- Ergonomic design. New contoured finger trigger for maximum control and comfort
- Suits left or right hand users
- Superior Valve system operated for longer without clogging up with contaminants
- Durable stainless steel and solid brass construction
- 2 Year warranty
- Anti-topple flange for extra stability

Kit Includes:

- Cry-Ac Spray Unit
- Five spray apertures with varying aperture sizes (1 x A, 1 x B, 2 x C & 1 x D)
- 28mm Straight spray needle, 20 gauge
- 78mm Curved spray extension, 20 gauge
- Tip holder
- Instruction and technique booklet



Open Spray Attachments

Aperture Spray Tips—A 1mm, B 0.8mm, C 0.6mm and D 0.4,,
Straight and Bent Spray Extension

Cryoplate—Transparent plastic plate with 4 conical openings—3, 5, 8 and 10mm. Limits the spread of the freeze zone and protects sensitive areas.

Closed Probes Attachments

Closed Probes provide precise control and excellent visibility Used to treat eyeline margins and difficult spots.

Round Mini Probes—1,2,3,4,5 & 6mm

Rectangular Mini Probes—1.5 x 3mm and 2.5 x 5mm.

Description	Product Code
CRYOSPRAY - BRYMILL CRY-AC 0.5L W/ APERTURES (SET)	DAB700
ACCESSORY—BRYMILL CRYOPLATE (EACH)	DAB303
ACCESSORY—BRYMILL CRYOTOTE #INSULATED BAG (EACH)	DAB601

Brymill Cry-Ac 150ML Cryospray

Made in the USA and the number 1 choice of dermatologists worldwide. The CRY-BABY® is a 150ml version of the larger Brymill units perfect for accurate, controlled treatment of smaller common lesions.

- 2-3 Hours static holding time
- 11.4cm tall, weighs 120g when full
- Same features as 0.5L version with same items in the kit.



Description	Product Code
CRYOSPRAY - BRYMILL CRY-AC 150ML W/ APERTURES (SET)	DAB400

Liquid Nitrogen Storage and Handling

- 10L and 20L Dewar Storage Containers for Liquid Nitrogen
- 10L and 20L Dewar Extractors
- 10L and 20L Dewar Ladles



These items are price on application—please enquire with customer service if interested.

CryOmega Disposable Cryopen

Clinically effective direct application of liquid nitrous oxide (-89°C) in a convenient and economical device. CryOmega® dispenses liquid nitrous oxide (-89°C) the coldest portable cryogen and it is cleared for 16 common lesions including verrucae, skin tags, and seborrheic keratoses. CryOmega® is an ideal choice in portable cryosurgery for physicians looking for convenience and great first-time results.

EFFECTIVE - Delivers a continuous, direct spray of nitrous oxide (-89°C) to treatment site. With a lower operating temperature compared to canister-based systems, provides clinically effective freezing with pinpoint precision.

CONVENIENT - Simple, hassle-free operation. Activate the internal cartridge once and the CryOmega® device is ready for use. No need for additional accessories or replacement cartridges. 30 to 45 five-second sprays per device.

ECONOMICAL - Low-cost, effective cryosurgery in any clinical setting.

DISPOSABLE - Simply dispose of the entire unit when the internal cartridge is empty.

Single Pack Contains:

1 Standard Pen with 16g of N2O cryogen

1 Instructions for Use

Twin Pack Contains:

2 Standard Pens with 16g of N2O cryogen

1 Instructions for use

Description	Product Code
CryOmega® Single Pack INTL	ECO160-2001
CryOmega® Twin Pack INTL	ECO160-2004



Sage Disposable Instruments

- Blue Handles for easy recognition
- Sharp each time—sharp, sterile and ready to go
- Sterile Guarantee
- Sold as single units

Description	Product Code
FORCEPS - S/S STERILE DRESSING BLOCK END BLUNT 12.5CM (EACH)	SA11011015
FORCEPS - S/S STERILE ADSON DRESSING FINE 12CM (EACH)	SA11011018
FORCEPS - S/S STERILE ADSON TISSUE 1:2 (EACH)	SA11011019
FORCEPS - S/S STERILE MOSQUITO STR 12.5CM (EACH)	SA11011020
FORCEPS - S/S STERILE MCINDOE DRESSING 15CM (EACH)	SA11011021
FORCEPS - S/S STERILE MOSQUITO CVD 12.5CM (EACH)	SA11011022
FORCEPS - S/S STERILE KELLY HEMOSTAT 14CM (EACH)	SA11011028
FORCEPS - S/S STERILE SPLINTER 12.5CM (EACH)	SA11011029
FORCEPS - S/S STERILE LABORDE DIALATING #TRIBLADE (EACH)	SA11011031
FORCEPS - S/S STERILE MAGILL 24CM (EACH)	SA11011084
FORCEPS - S/S STERILE CRILE ARTERY 1:2 14CM (EACH)	SA11011134
NEEDLEHOLDER - S/S MAYO STR 20CM (EACH)	SA11011156
NEEDLEHOLDERS - S/S STERILE MAYO-HEGAR 15CM (EACH)	SA11011023
SCISSORS - S/S STERILE IRIS 11.5CM (EACH)	SA11011024
SCISSORS - S/S STERILE SH/BL 12.5CM (EACH)	SA11011025
SCISSORS - S/S MAYO SH/SH STR 17CM (EACH)	SA11011151

Surgical Blades Swann-Morton

- Made from Carbon Steel
- Gamma Sterilised
- Box of 100
- Blades 06 to 16, 40, Std Stitch Cutter fit No. 3 Handles and 18 to 36 fit No. 4 Handles



Description	Product Code
SURGICAL BLADES STANDARD NO. 10 0201 (B/100)	BF6950
SURGICAL BLADES STANDARD NO. 11 0203 (B/100)	BF6951
SURGICAL BLADES STANDARD NO. 12 0204 (B/100)	BF6952
SURGICAL BLADES STANDARD NO. 15 0205 (B/100)	BF6955
SURGICAL BLADES STANDARD NO. 20 0206 (B/100)	BF6960
SURGICAL BLADES STANDARD NO. 21 0207 (B/100)	BF6961
SURGICAL BLADES STANDARD NO. 22 0208 (B/100)	BF6962
SURGICAL BLADES STANDARD NO. 23 0210 (B/100)	BF6963
SURGICAL BLADES STANDARD NO. 24 0211 (B/100)	BF6964
SURGICAL BLADES STANDARD NO. 10A 0202 (B/100)	BF6968
STITCH CUTTERS STANDARD 0420 (B/100)	BF6999

Scalpel Handles

- Stainless Steel Reusable Scalpel Handles
- L.R.I. Ward Quality Instruments

Description	Product Code
SCALPEL HANDLE - L.R.I. NO. 3 (EACH)	RH9250
SCALPEL HANDLE - L.R.I. NO. 4 (EACH)	RH9251

Disposable Scalpels Swann-Morton

- Made from Stainless Steel with Polystyrene Handle.
- Gamma Sterilised
- Box of 10—The sterile scalpels are individually packaged into a transparent thermos-formed tray



Description	Product Code
SCALPELS & HANDLES DISPOSABLE NO. 10 0501 (B/10)	BF6910
SCALPELS & HANDLES DISPOSABLE NO. 11 0503 (B/10)	BF6911
SCALPELS & HANDLES DISPOSABLE NO. 15 0505 (B/10)	BF6915
SCALPELS & HANDLES DISPOSABLE NO. 20 0506 (B/10)	BF6920
SCALPELS & HANDLES DISPOSABLE NO. 22 0508 (B/10)	BF6922
SCALPELS & HANDLES DISPOSABLE NO. 23 0510 (B/10)	BF6923

Disposable Scalpels Remover Swann-Morton



The **BladeFLASK** scalpel blade remover is simple to use with a quick, single-handed “push, click, and pull” action to remove and contain used scalpel blades. Used together with a hands-free passing technique, the Qlicksmart BladeFLASK can prevent up to 5 times more injuries than a safety scalpel.

Description	Product Code
QLICKSMART SCALPEL BLADE REMOVAL SYSTEM (EACH)	VH1390030
QLICKSMART UNIVERSAL BRACKET (EACH)	VH1390027

Reusable Surgical Instruments

Liberty Stainless Steel Instruments—Made in Pakistan to German specifications
 German Made and Finished AMS Stainless Steel Theatre Quality Instruments are also available to order.
 Please note this is a limited range—please enquire if instrument is not listed.

Artery Forceps	
Description	Liberty Product Code
Crile Straight 14cm	LIAFCS14
Crile Curved 14cm	LIAFCC14
Crile Renkin Str 16cm	LIAFCRS16
Crile Renkin Crv 16cm	LIAFCRC16
Hartmann Straight 10cm	LIAFHS10
Hartmann Curved 10cm	LIAFHC10
Kelly Straight 14cm	LIAFKS14
Kelly Curved 14cm	LIAFKC14
Mosquito Str 12.5cm	LIAFMS12
Mosquito Curved 12.5cm	LIAFMC12
Mosquito Straight 1/2 Teeth 12.5cm	LIAFMST12
Mosquito Curved 1/2 Teeth 12.5cm	LIAFMCT12
Mosquito Micro Straight 12.5cm	LIAFMMS12
Mosquito Micro Curved 12.5cm	LIAFMMC12
Ochsner-Kocher Str 20cm	LIAFOKS20
Ochsner-Kocher Curved 20cm	LIAFOKC20
Rochester Pean Straight 18cm	LIAFRPS18
Rochester Pean Curved 18cm	LIAFRPC18
Spencer Wells Straight 18cm	LIAFSWS18
Spencer Wells Curved 18cm	LIAFSWC18

Dilators	
Description	NOPA Theatre Product Code
St Marks Rectal Small 12-19mm	RH910
St Marks Rectal Large 15-25mm	RH9101

Aural Forceps	
Description	Liberty Product Code
Hartmann Croc Action 11cm	LIAFHC11
Hartmann Croc Action 14cm	LIAFHC14
Tilley 14cm	LIAFT14
Tilley 16cm	LIAFT16

Chalazion and Cilia Forceps	
Description	Liberty Product Code
Lambert Chalazion 9cm	LILCF09
Beer Cilia 9cm	LICFB09
Douglas Cilia 9cm	LILCFD09

Chiropody Pliers	
Description	Liberty Product Code
Plier Straight Blade 14cm	LICPSB14
Plier Curved Blade 14cm	LICPSS14
Cuticle Cutter 11cm	LICPCN

Curettes	
Description	Liberty Product Code
Buck Ear Blunt 0 15cm	LICEBB0
Buck Ear Blunt 00 15cm	LICEBB00
Buck Ear Blunt 1 15cm	LICEBB1
Buck Ear Sharp 0 15cm	LICEBS0
Buck Ear Sharp 00 15cm	LICEBS00
Buck Ear Sharp 1 15cm	LICEBS1
Williger Bone 00+0 14cm	LICBW00/0
Williger Bone 0+1 14cm	LICBW0/1
Williger Bone 1+2 14cm	LICBW1/2
Excavator Double Ended w/ Holes 14cm	LICEDE14

Reusable Surgical Instruments

Dressing Forceps	
Description	Liberty Product Code
Blunt 13cm	LIDFB13
Pointed 13cm	LIDFP13
Blunt 16cm	LIDFB16
Blunt 18cm	LIDFB18
Blunt 20cm	LIDFB20
Blunt 25cm	LIDFB25
Blunt 30cm	LIDFB30

Ear Syringe and Eye Magnet	
Description	Liberty Product Code
Ear Syringe 100ml O Ring Style w/ shield and 3 nozzles	RH9256
Eye Magnet and Loop	LIHEMAL

Files	
Description	Liberty Product Code
Blacks Files Angled Fine	LIBFAF
Blacks Files Angled Medium	LIBFAM
Blacks Files Straight Fine	LIBFSF
Blacks Files Straight Medium	LIBFSM

Introducing Forceps	
Description	Liberty Product Code
Magill Infant 15cm	LIFM15
Magill Child 20cm	LIFM20
Magill Adult 25cm	LIFM25

Graefe Iris Forceps	
Description	Liberty Product Code
Iris Graefe Half Curved 10cm	LIFGHC10
Iris Graefe Half Curved 10cm 1/2 Teeth	LIFGHCT10

Nasal Specula	
Description	Liberty Product Code
Thudicum #0	RH9344
Thudicum #1	RH9345
Thudicum #2	RH9346
Thudicum #3	RH9347
Thudicum #4	RH9348

Needle holders	
Description	Liberty Product Code
Crilewood 15cm	LINHCW15
Crilewood 15cm TC	LINHCW15TC
Derf 13cm	LINHD13
Derf 13cm TC	LINHD13TC
DeBakey 20cm TC	LINHD20TC
Kilner 13cm TC	LINHK13TC
Mayo Hegar 14cm	LINMH14
Mayo Hegar 14cm TC	LINMH14TC
Mayo Hegar 16cm	LINMH16
May Hegar 16cm TC	LINMH16TC
Mayo Hegar 20cm	LINMH20
Mayo Hegar 20cm TC	LINMH20TC

Pin Wheels and Reflex Hammer	
Description	Liberty Product Code
Wartenburg Pin Wheel	RHWPW
Nat Hosp Pat Aus Reflex Hammer	RH9312
Taylor's Neuro Reflex Hammer 18cm	RH9313

Reusable Surgical Instruments

Probes	
Description	Liberty Product Code
Farrell Cotton Applicator 13cm	LIPFCA13
Probe 2mm 14cm	LIP13
Probe 2mm 16cm	LIP16
Probe 2mm 18cm	LIP18
Jobson-Horne Probe & Applicator 18cm	LIPJH18

Ring Saw	
Description	L.R.I. Product Code
Finger Ring Saw w/ blade	LIRS
Blade for Finger Ring Saw	LIRSSB

Scissors	
Description	Liberty Product Code
Lister Bandage 14cm	LISBL140
Iris Straight 9cm	LISIS09
Iris Curved 9cm	LISIC09
Iris Straight 11cm	LISIS11
Iris Curved 11cm	LISIC11
Iris Straight 12cm	LISIS12
Iris Curved 12cm	LISIC12
Kilner Straight 12cm	LISKS12
Kilner Curved 12cm	LISKC12
Braun Stadler Episiotomy 22cm	LISBS22
Mayo Straight 15cm	LISMS15
Mayo Curved 15cm	LISMC15
Mayo Straight 17cm	LISMS17
Mayo Curved 17cm	LISMC17
Mayo Straight 20cm	LISMS20
Mayo Curved 20cm	LISMC20

Scissors	
Description	Liberty Product Code
Metzenbaum Straight 14cm	LISMES14
Metzenbaum Curved 14cm	LISMEC14
Metzenbaum Straight 16cm	LISMES16
Metzenbaum Curved 16cm	LISMEC16
Metzenbaum Straight 18cm	LISMES18
Metzenbaum Curved 18cm	LISMEC18
Metzenbaum Straight 20cm	LISMES20
Metzenbaum Curved 20cm	LISMEC20
Metzenbaum Straight 23cm	LISMES23
Metzenbaum Curved 23cm	LISMEC23
Spencer Suture 9cm	LISS09
Spencer Suture 11cm	LISS11
First Aid 13cm Straight SH/SH	LISSSHSF125
First Aid 13cm Straight SH/BL	LISSHBLF125
First Aid 13cm Straight BL/BL	LISBHLH125
Operating SH/SH Straight 13cm	LISSSHS13
Operating SH/BI Straight 13cm	LISSHBL13
Operating BI/BI Straight 13cm	LISBLBL13
Operating BI/BI Curved 13cm	LISBLBLC13
Operating SH/SH Straight 16cm	LISSSHS16
Operating SH/BI Straight 16cm	LISSHBL16
Operating BI/BI Straight 16cm	LISBLBL16
Operating SH/SH Straight 18cm	LISSSHS18

Reusable Surgical Instruments

Scissors	
Description	L.R.I Product Code
Universal 18cm w/ Auto-clavable Plastic Black	LISU180B

Sinus Forceps	
Description	Liberty Product Code
Lister 16cm	LISFL16

Splinter Forceps	
Description	Liberty Product Code
Feichenfeldt 11cm	LISFF11
Hunter Straight 11cm	LISFHS11
Hunter Curved 11cm	LISFHC11
Fine Splinter w/ Magnifier	LISFWM

Sponge Holding Forceps	
Description	Liberty Product Code
Foerester Straight 18cm	LISHFF18
Foerester Straight 24cm	LISHFF24

Sterilisation Forceps and Pins	
Description	Liberty Product Code
Cheatles Light 26cm	LISFCLP260
Cheatles Heavy 28cm	LISFCLP280

Tenaculum Forceps	
Description	Liberty Product Code
Schroeder 24cm	LIUFTS25
Skene 24cm	LIUFTS24

Tissue Forceps	
Description	Liberty Product Code
Micro Adson 12cm	LITFAM12
Micro Adson 1/2 T 12cm	LITFAMT12
Adson 12cm	LITFAP12
Adson 1/2T 12cm	LITFAT12
Adson Braun 7/7T 12cm	LITFAB12
Babcock 20cm	LITFB20
Gillies 1/2T 15cm	LITFGT15
Lanes 1/2T 15cm	LITFLT151/2
Lanes 2/3T 15cm	LITFLT152/3
Potts Smith 18cm	LITFPS18
Potts Smith 18cm 1/2T 18cm	LITFPST18

Towel Clips	
Description	Liberty Product Code
Backhaus 13cm	LITCB13
Backhaus 13cm w/ Ball & Socket	LITCBWS13
Backhaus 13cm w/ Plates	LITCBWP

Watchmakers Forceps	
Description	Liberty Product Code
Watchmakers #3	LIWMF3
Watchmakers #4	LIWMF4
Watchmakers #5	LIWMF5
Watchmakers #7	LIWMF7