

General Purpose Medical Carts and Trolleys

Product range for Medical Carts and Trolleys includes:

- Dressing, Instrument and Perspex Drawer Trolleys.
- Stainless Steel Transport Trolleys and Plastic Tub Trolley.
- Chart Trolleys
- Flow Chart Table and I.V. Poles
- Cylinder Trolleys
- Customisable Options and Powder Coated Colours Chart



Suppliers of quality products to the Healthcare Industry est 1889.
Australian Owned & Operated.

Web dbhealth.com.au

Email customerservice@denyer.com.au

Freecall: 1800 451 446

ABN: 37 000 459 174

Emery Dressing Trolleys

Available in Standard Stainless Steel or Stainless Steel with Powder coated drawer fronts.

- Precision ball bearing non breaking swivel castors enable the weight of the unit to be absorbed by the bearing. This design ensures that the unit will be able to be moved with minimal effort.
- HQSS11, HQSS12 and HQSS14 models of dressing trolleys have 4 x 100mm Zinc plated non breaking swivel casters. HQSS12, HQSS13 and HQSS15 models of dressing trolleys come with 4 x 125mm Zinc plated non-braking swivel casters
- Full Extension Ball Bearing drawer runners provide easy gliding of drawer movement, this allows full access to the contents of the drawer with ease. These drawers will not fly open when the trolley is moved.
- Made from 304 Welded Stainless Steel
- For safety reasons, our products have full safety folds that are standard on all our products.
- Large curved handles are entirely smooth to minimise infection control and make grasping easier
- Drawer sized 120mm deep
- Dressing Trolleys come in 3 standard size widths 50cm, 75cm and 100cm x 50cm depth and 90cm in height. Half Width Dressing Trolleys come in 2 width of 80cm and 96 cm x 50cm depth and 90cm height
- Durable and reliable medical equipment
- Customisable options available

Dressing Trolleys w/ 1 full width drawer

Trolley Width available x 50cm (D) x 90cm (H)	Product Code
DRESSING TROLLEY - EMERY 2 LVL S/S 50CM W/ 1 DRAWER (EACH)	HQSS10
DRESSING TROLLEY - EMERY 2 LVL S/S 75CM W/ 1 DRAWER (EACH)	HQSS10-1
DRESSING TROLLEY - EMERY 2 LVL S/S 100CM W/ 1 DRAWER (EACH)	HQSS10-2



Dressing Trolleys w/ 2 full width drawers

Trolley Width available x 50 (D) x 90cm (H)	Standard Stainless Steel
DRESSING TROLLEY - EMERY 2 LVL S/S 50CM W/ 2 DRAWERS (EACH)	HQSS11
DRESSING TROLLEY - EMERY 2 LVL S/S 75CM W/ 2 DRAWERS (EACH)	HQSS11-1
DRESSING TROLLEY - EMERY 2 LVL S/S 100CM W/ 2 DRAWERS (EACH)	HQSS11-2



Dressing Trolleys w/ 3 full width drawers

Trolley Width available x 50 (D) x 90cm (H)	Standard Stainless Steel
DRESSING TROLLEY - EMERY 2 LVL S/S 50CM W/ 3 DRAWERS (EACH)	HQSS12
DRESSING TROLLEY - EMERY 2 LVL S/S 75CM W/ 3 DRAWERS (EACH)	HQSS12-1
DRESSING TROLLEY - EMERY 2 LVL S/S 100CM W/ 3 DRAWERS (EACH)	HQSS12-2



Dressing Trolleys w/ 4 full width drawer

Trolley Width available ** (W) x 50 (D) x 90cm (H)	Standard Stainless Steel
DRESSING TROLLEY - EMERY 2 LVL S/S 50CM W/ 4 DRAWERS (EACH)	HQSS13
DRESSING TROLLEY - EMERY 2 LVL S/S 75CM W/ 4 DRAWERS (EACH)	HQSS13-1
DRESSING TROLLEY - EMERY 2 LVL S/S 100CM W/ 4 DRAWERS (EACH)	HQSS13-2



Dressing Trolleys w/ 2 1/2 width drawers

DRESSING TROLLEY - EMERY 2 LVL S/S 96CM W/ 2 1/2 WIDTH DRAWERS (EACH)	HQSS14
DRESSING TROLLEY - EMERY 2 LVL S/S 96CM W/ 2 1/2 WIDTH DRAWERS W/ LOCKING CASTERS (EACH)	HQSS14-1
DRESSING TROLLEY - EMERY 2 LVL S/S 80CM W/ 2 1/2 WIDTH DRAWERS (EACH)	HQSS14-3



Dressing Trolleys w/ 4 1/2 width drawers

DRESSING TROLLEY - EMERY 2 LVL S/S 96CM W/ 4 1/2 WIDTH DRAWERS (EACH)	HQSS15
DRESSING TROLLEY - EMERY 2 LVL S/S 80CM W/ 4 1/2 WIDTH DRAWERS (EACH)	HQSS15-3



Perspex Drawer Trolleys

In medical settings and hospitals, efficient organization and storage are paramount. Introducing our cutting-edge **Medical Perspex Drawer Trolley**, designed to streamline storage solutions while ensuring easy visual identification. These trolleys are meticulously crafted to meet the unique demands of medical environments.

- Available with 1 or 2 drawers
- 50cm (W) x 50cm (D) 90cm (H)
- Ergonomic handles on drawers
- Shelf with 3 Rails on sides
- Stainless steel welder frame
- 4 x 100mm Swivel casters
- Customisable

PERSPEX DRAWER TROLLEY - EMERY 1 LVL S/S 50CM W/ 1 DRAWER (EACH)	HQSS10VP
PERSPEX DRAWER TROLLEY - EMERY 1 LVL S/S 50CM W/ 2 DRAWERS (EACH)	HQSS11VPS



Pathology Trolley

- Welded 304 grade Stainless Steel frame
- Curved ergonomic handles on drawers
- 4 x 150mm Swivel Caster with 2 Locking
- 2 or 3 130mm deep Drawers with locking mechanism
- Comes with Perspex writing shelf
- Customisable options available

Model HQSP624

3 Drawers, 2 x Shelves - 1 x top shelf with 3 x Rails; 1 x Support Shelf under drawers;
1 x Perspex Writing Shelf

Model HQSP624-2

2 drawers, 2 shelves - Top Shelf with 3 rails; 1 x Perspex Writing Shelf



PATHOLOGY TROLLEY - EMERY 3 DRAWERS	HQSP624-1
PATHOLOGY TROLLEY - EMERY 2 DRAWERS	HQSP624-2

Emery Instrument Trolleys

- Made from 304 Welded Stainless Steel
- 4 x 100mm Zinc Plated Precision ball bearing castors enable the weight of the unit to be absorbed by the bearing. This design ensures that the unit will be able to be moved with minimal effort.
- Trolleys have full safety folds that are standard on all our products.
- Durable and reliable medical equipment
- Bottom Shelves are 30cm from the ground
- Available in 4 widths 50cm, 75cm, 100cm and 115cm X 50cm (D) x 90cm (H)
- Customisable options available

Description	Product code
F.T. INSTRUMENT TROLLEY - EMERY 1 SHELF S/S 50CM (EACH)	HQSS20
F.T. INSTRUMENT TROLLEY - EMERY 1 SHELF S/S 75CM (EACH)	HQSS21
F.T. INSTRUMENT TROLLEY - EMERY 1 SHELF S/S 100CM (EACH)	HQSS22
F.T. INSTRUMENT TROLLEY - EMERY 1 SHELF S/S 115CM (EACH)	HQSS23
F.T. INSTRUMENT TROLLEY - EMERY 2 SHELVES S/S 50CM (EACH)	HQSS24
F.T. INSTRUMENT TROLLEY - EMERY 2 SHELVES S/S 75CM (EACH)	HQSS25
F.T. INSTRUMENT TROLLEY - EMERY 2 SHELVES S/S 100CM (EACH)	HQSS26
F.T. INSTRUMENT TROLLEY - EMERY 2 SHELVES S/S 115CM (EACH)	HQSS27



HQSS20

Emery Instrument Trolleys with Side Rails

With the same features as the Flat Top Instrument Trolleys but with added side rails

- Available in 3 widths 50cm, 75cm, 100cm
- Customisable options available

INSTRUMENT TROLLEY - EMERY 2 SHELVES W/ RAILS S/S 50CM (EACH)	HQSS256
INSTRUMENT TROLLEY - EMERY 2 SHELVES W/ RAILS S/S 75CM (EACH)	HQSS256-1
INSTRUMENT TROLLEY - EMERY 2 SHELVES W/ RAILS S/S 100CM (EACH)	HQSS256-2



Stainless Steel Transport Trolleys

- Welded stainless steel frames
- Full safety folds
- 4 x Duffers as standard
- Push handles on both ends
- Precision ball bearing castors
- Available with either 3 or 4 shelves
- Recommended weight bearing capacity of 10kg per shelf



HQSS103



HQSS103.1 with Lipped Shelves

HQSS103

- 100cm (w) x 60cm (d) x 115cm (h)
- 3 x lipped shelves
- 4 x 100mm Swivel Castors

HQSS103-1

- 115cm (w) x 70cm (d) x 140cm (h)
- 4 flat shelves - Lipped shelves available on request
- 4 x 200mm Casters 2 x Brake, 2 x Directional



HQSP327

HQSP327 - Set Up Trolley

- 90cm (W) x 55cm (D) x 145cm (H)
- 3 shelves with rails, 3 x Top Rails, 1 x Extra Rail on Top Shelf
- 4 x 125mm Castors, 2 x Locking, 2 x Directional
- Customisable

HQSP288 - Drinks Trolley

- 117cm (W) x 70cm (D) x 112cm (H)
- 3 lipped shelves,
- push handles on both ends,
- 4 casters - 2 swivel, 2 locking
- Customisable



HQSP288

HQSP403 - Lipped Trolley

- 90cm (W) x 50cm (D) x 90cm (H)
- 2 lipped shelves with 3 rails on each shelf
- 1 push handle; 2 door handles
- 4 x 100mm Castors
- 4 buffers
- Customisable



HQSP403

CLEARING TROLLEY - EMERY S/S #3 SHELVES # 100 X 60 X 115CM (EACH)	HQSS103
CLEARING TROLLEY - EMERY S/S #4 SHELVES # 115 X 70 X 140CM (EACH)	HQSS103-1
SET UP TROLLEY - EMERY S/S #3 SHELVES W/RAIL 90 X 55 X 145CM (EACH)	HQSS327
DRINKS TROLLEY - EMERY S/S #3 LIPPED SHELVES 117 X 70 X 112CM (EACH)	HQSS288
LIPPED TROLLEY - EMERY S/S #2 LIPPED SHELVES W/ RAIL 90 X 50 X 90CM (EACH)	HQSS403

Plastic Bin Trolley

The Bin Trolley is designed to hold 2 x Bins.

- Push handle on end.
- Bin Dimensions: 64.5 (W) x 41 (L) x 21cm (H).
- Bins w/ lids are included with trolley.
- Wheels: 10cm standard w/ 2 braking
- Customisable



S/S BIN TRLY - EMERY S/S W/ 2 BINS #72.5 X 42.5 X 100CM (EACH)	HQSS443
--	---------

Ring Binder / Chart Station Trolley

Chart Station Trolley

- Welded stainless steel frame with high density compact laminate chart bays
- Dimensions: 125 (L) x 45 (W) x 110cm (H)
- 4 x 12.5cm Caster - 2 Locking and 2 swivel
- Individual bay sizes are 9cm (W) x 40cm (D) x 45cm (H)
- Customisable



Chart Trolley

- Powder coated Stainless Steel
- 85cm (W) x 40cm (D) x 112cm (H)
- 2 x Angled Shelves
- 4 x 100mm Non Breaking Casters
- Customisable.



CHART STATION TROLLEY - EMERY 30 S/S & HD LAMINATE (EACH)	HQSP650
CHART TROLLEY - EMERY S/S P/C EACH)	HQSP674

HQSP674 with
HQAO-CAST

DB Health has a variety of A4 Patient Polypropylene 2D-Ring Folders available in a range of colours and widths
Patient Polypropylene Folders can be found in DBH2401 - Patient PP Ring Binders and Accessories

Oxford Chart Trolley

Chart Trolley are designed to easily move a large number of charts as required. These trolleys when ordered come with the requisite number of chart holders. There are 3 models of Chart Trolleys which are 12, 24 and 36.

- Shelves: 1 x 25cm clearance from bottom shelf to bottom of chart holders.
- Wheels: 12.5cm standard w/ 2 braking and wall bumpers
- Optional locking lid
- Available in Stainless Steel or Powder Coated

CHART TROLLEY - OXFORD S/S W/ 12 CHART HOLDERS #41.2 X 50 X 90CM (EACH)	OX10-1100
CHART TROLLEY - OXFORD P/C W/ 12 CHART HOLDERS #41.2 X 50 X 90CM (EACH)	OX55-1100
CHART TROLLEY - OXFORD S/S W/ 24 CHART HOLDERS #75.2 X 50 X 90cm(EACH)	OX10-1200
CHART TROLLEY - OXFORD P/C W/ 24 CHART HOLDERS #75.2 X 50 X 90cm(EACH)	OX55-1200
CHART TROLLEY - OXFORD S/S W/ 36 CHART HOLDERS #112 X 50 X 90cm(EACH)	OX10-1300
CHART TROLLEY - OXFORD P/C W/ 36 CHART HOLDERS #112 X 50 X 90cm(EACH)	OX55-1300
ACCESSORY - OXFORD CHART HOLDER S/S (EACH)	OX10-1400
ACCESSORY - OXFORD CHART HOLDER P/C (EACH)	OX10-1500



OX55-1100

Flow Chart Table

- The working surface pivots to enable the user to select the most comfortable working angle
- Precision ball bearing castors enable the weight of the unit to be absorbed by the bearing.
- Compact Laminate Top – White 880 x 670 x (900-1200)mm.
- 30mm lip on lower side of the 880mm.
- 3 x Stainless Steel File Clamp across top of desk .
- 4 x 100mm Casters , 2 x Locking 2 x Swivel.
- 1 x file clamp on right hand side.
- Pivot Locking Knob (both sides).
- Gas Strut Height Adjustable.
- Wire basket mounted underneath desk.
- Customisable.



FLOW CHART TABLE - EMERY LARGE S/S W/ BASKET (EACH)

HQSP562-2

Interstackable I.V. Poles

Benefits

- These IV poles are space saving as they have the ability to nest together and can fit 10 of these poles in the space of 3 normal poles to maximise storage.
- IV Head is height adjustable.
- Extra low centre of gravity for increased stability.
- Available with weighted or non-weighted bases .
- Comes with 2 Prong Pigs Tail Head.
- 5 x 75mm Standard Casters
- Customisable options available.

IV POLE NESTABLE W/PIGS TAIL HEAD S/S # INTERNAL WEIGHTED BASE (EACH)	HQSS04IWBNEST
IV POLE NESTABLE W/PIGS TAIL HEAD S/S (EACH)	HQSS04NEST



IV Poles

These IV Poles are available with:

- IV heads - 2 Prong Pigtailed Head.
- Welded Stainless Steel or Welded Stainless Steel Weighted Base, —53cm diameter with 5 75mm casters.
- Height adjustable IV Pole from 133cm to 243cm.
- 6kg Internal Weights Base on weighted models.
- Customisable options available.

IV POLE W/ PIGS TAIL HEAD S/S (EACH)	HQSS04
IV POLE W/ PIGS TAIL HEAD S/S #INTERNAL WEIGHTED BASE (EACH)	HQSS04IWB



Cylinder Bottle Stands and Trolleys

Trolley for size C and D sized cylinder bottles

- Caster 4 x 125mm Swivel, 4 x Buffers,
- Push handles on both ends

Trolley caddy's for single C,D, E and G sized cylinder

- 1 Push handle
- 2 x 150mm fixed wheels
- 1 strap to secure cylinder

Portable stands

- Portable and for the use on the side of a patient bed or side table
- Hook to attach to trolley
- Holds 1 x C size cylinder or 1 x CD size cylinder
- Ring has internal diameter of 125mm with 9.5mm stainless steel



HQSP623-1P



HQSP623P



HQSP623C6



HQSP82D

Description	Code
CYLINDER HOLDER TROLLEY - EMERY 6 UNIT MOBILE SIZE C 49 x 33.5 x 100CM (H) (EACH)	HQSP623C6
CYLINDER HOLDER TROLLEY - EMERY 12 UNIT MOBILE SIZE C 640 x 49 x 100CM (H) (EACH)	HQSP623C12
CYLINDER HOLDER TROLLEY - EMERY 6 UNIT MOBILE SIZE D 61 x 42 x 115CM (H) (EACH)	HQSP623D6
SINGLE CYLINDER CADDY TROLLEYS	
CYLINDER HOLDER TROLLEY - EMERY SIZE C SIZE S/S 22 x 22 x 90CM (H) (EACH)	HQSP82C
CYLINDER HOLDER TROLLEY - EMERY SIZE D SIZE S/S 29.5 x 30 x 91.5CM (H) (EACH)	HQSP82D
CYLINDER HOLDER TROLLEY - EMERY SIZE E SIZE S/S 29.5 x 30 x 110CM (H) (EACH)	HQSP82E
CYLINDER HOLDER TROLLEY - EMERY G SIZE S/S 29.5 x 30 x 91.5CM (H) (EACH)	HQSP82G
PORTABLE / BED RAIL MOUNTABLE	
CYLINDER HOLDER PORTABLE C SIZE CYLINDER S/S 20 x 20 x 55CM (H) (EACH)	HQSP623P
CYLINDER HOLDER PORTABLE CD SIZE CYLINDER S/S 14.5 x 17 x 58CM (H) (EACH)	HQSP623-1P
WALL MOUNTING BRACKET FOR HQSP623P AND HQSP623-1P	HQSP623-WB

If a stationary unit is required please inform customer service when requesting your quote

Customisable Options for Emery Trolleys

Upgrade Options	Product Code
Shelf Rails Bottom	HQAO-SRB
Braking Caster	HQAO-CAST
Bucket holder LHS Mounted	HQAO-BUCKET HOLDER
Chart Holder***	HQAO-CH
Drawer Dividers / Trays	HQAO-PDD001
Drawer Dividers / Trays	HQAO-PDD001-1
Drawer Dividers / Trays	HQAO-PDD002
Drawer Dividers / Trays	HQAO-PDD005
Glove Box Holder*** (3 x Ansell Boxes)	HQAO-GLOVE BH
IV Pole Mount on corner	HQAO-IV-TROLLEY
IV Pole Flat Top Trolley Mount on corner	HQAO-IV TROLLEY FT
Keyed Lock - RHS Mounted**	HQAO-LC****
Magnifying Lamp	HQAO-MAG LIGHT TROLLEY
Mesh Drawers Bottom	HQAO-MESH DRAWER
Mesh Base SS Catheter Holder***	HQAO-CH-150-800
Push Handle	HQAO-PH
Drop down Shelf RHS Mount	HQAO-DDS
Size C Gas Bottle Holder	HQAO-GAS BOTTLE HOLDER
Size C Gas Bottle Holder Portable***	HQAO-GBHP

*HQAO-PDD001

Drawer Dividers to fit Standard 500 x 500 Unit
3mm Clear Polycarbonate to fit One drawer
2 of 400 x 105 high dividers with 5 slots
2 of 365 x 105 high dividers with 5 slots
Creates 9 Compartments

*HQAO-PDD001-1

Drawer Dividers to fit Standard 500 x 500 Unit
3mm Clear Polycarbonate to fit One drawer
1 of 400 x 105 high dividers with 5 slots
1 of 365 x 105 high dividers with 5 slots
Creates 4 Compartments

*HQAO-PDD002

3mm Clear Drawer Dividers – T Divider
Fit drawers on Standard 500 x 500 unit
1 of 390 x 115 high T dividers
Run front to back only
Creates 2 Compartments

*HQAO-PDD005

Drawer Dividers to fit Standard 750 x 500 Unit
3mm Clear Polycarbonate to fit One drawer
2 of 665 x 105 high dividers with 5 slots
4 of 365 x 105 high dividers with 5 slots
Creates 15 Compartments



HQAO-PDD001

All accessories are added during manufacture

* Please select required Drawer Dividers

** Fitted to Right hand side unless another accessory is on the right hand side

*** For Trolleys with Side Rails only

**** Locking System can be for all drawers or a combination of drawers. Please state how many drawers require to have the locking system when placing your order.

Emery I.V. Pole Customisable Options

Upgrade Options	Product Code
Powder Coated IWB I.V. Pole Base	HQAO-PCIVIWB
2 Hook Rectangular Rams Head	HQAO-U2 HOOKS
4 Hook Rams Head	HQAO-RH4
5 Hook Rams Head	HQAO-RH5
IV Push Handle	HQAO-IV PUSH HANDLE 3
IV Push Handle Quarter Ring 380mm Diameter	HQAO-IV PUSH HANDLE
IV Size C gas bottle hold	HQSS04IWB-GBH
Mid Bag Hooks (x2)	HQAO-MBH-MID BAG HOOK
Height Adjustable Basket to Fit on IV Poles	HQAO-IV-BASKET



HQAO-PCIVIWB
For internally weight base



HQAO-IV PUSH
HANDLE 3



HQAO-IV PUSH HANDLE



HQAO-U2 HOOKS



HQAO-RH4



HQSS04IWB-GBH
"C" size gas bottle holder. 520 w x 530 d x (1350-
2300mm) h. For internal weighted base



HQAO-MBH-MID BAG
HOOK x2



HQAO-IV-BASKET
Height Adjustable Basket to Fit on IV
Poles
Dimensions: 330 x 200 x 40mm H
Construction: Stainless Steel / Welded

Emery Powdercoat Colours

PEARL WHITE



BIRCH WHITE



LEMON YELLOW



PINK BLUSH



FLAME RED



VIPER GREEN



POMMEL BLUE



FRENCH BLUE



DEEP VIOLET

