

Maternity Department Products

Product range includes:

- Soft Foam Wristbands,
- Print on Demand wristbands
- Mother/Father and Infant Wristbands
- Barcode Ready Wristbands and Tags
- Insert Wristband
- Wristband Extender Tether
- Foot Printer
- Family Wristbands and Stickers
- Breast Milk Labels
- Patient Chart, Drug line and Nursing Labels
- Baby Scales and Measuring Mat
- Bassinets, baby Baths and Mobile Cots
- Products for Consideration



Web dbhealth.com.au

Email customerservice@denyer.com.au

Freecall: 1800 451 446

ABN: 37 000 459 174

Suppliers of quality products to the Healthcare Industry Est 1889.
Australian Owned & Operated.

Precision Foam Wristbands—Neonatal and Newborn

The Precision Wristband helps meet the special needs of patients with sensitive skin for Neonatal and Premie Infants as well as elderly patients. Available in 4 sizes, this moveable ID band offers the highest quality, convenience, and comfort.

- Constructed of soft, comfortable foam with hook & loop Velcro® style closure - No cutting necessary to remove.
- Blank ID panel accepts handwritten information or label and Adhesive shield protects critical patient ID information and barcodes.
- Easily repositioned on wrist or ankle to accommodate IV site changes
- Ensures ID band is present on limb instead of taped to incubator, guaranteeing positive patient identification.

Description	Product Code
ID BAND NEONATAL NEWBORN SOFT FOAM BLUE-WHITE W/SHIELD (B/12)	PD523-11
ID BAND NEONATAL PREMIE SOFT FOAM BLUE-WHITE W/ SHIELD (B/12)	PD533-11

	Length	Max Label Size
Neonatal Newborn	6.5cm	5.7 x 2.5cm
Neonatal Premie	13.5cm	5.7 x 2.5cm



PD533-11

Scanband DR Direct Thermal Wristbands

Best-selling Direct Thermal wristband on the market to withstand the harshest conditions in the hospital environment.

- Provides a high-quality bar code for scanning, accurate patient ID, Supports text, linear and 2D bar codes, photos and graphics.
- Top quality tamper-evident adhesive closure
- Unbeatable print quality and Extremely easy to print and apply to patients
- Features a soft material, ideal for patients with sensitive skin.

Product	Wristband Dimensions	Print Area	Core Size	Black Bar	Wound	Outer Diameter	Closure
PD8180-**	17.8cm x 1.9cm	7.3 x 1.9cm	2.5cm	At perforation	Out	9cm	Adhesive

Description	Product Code
DT WRISTBAND - SCANBAND SOFT INFANT ADH CLOSURE WHITE (B/2 ROLLS OF 200)	PD8180-11
DT WRISTBAND - SCANBAND SOFT INFANT ADH CLOSURE RED (B/2 ROLLS OF 200)	PD8180-16



PD8180-11 Series—Infant

Precision® Infant Direct Thermal Wristbands

Superior imprint quality, maximum resistance to solvents, and seamless printing with the patented PermaPrint® surface. Also provides a flat surface for best bar code scanning results, eliminating curvature issues.

Easy to apply adhesive closure ensures safety and security

Antimicrobial additive protects wristband surface against tested non-pathogenic bacteria*

Compatible with most direct thermal printers

PD7710 has Black Bar at perforation



Description	Product Code
ID BAND - PRECISION INFANT DIRECT THERMAL W/ ADH CLOSURE WHITE #2 ROLLS (B/400)	PD7701-11
ID BAND - PRECISION INFANT DIRECT THERMAL W/ ADH CLOSURE WHITE #2 ROLLS (B/400)	PD7710-11

Speedi-Print® Infant Wristbands

Soft vinyl material and Latex Free
 SecurSnap closure to prevent transfer
 ClearImage® surface protects handwritten and embossed information
 Information Area: 57 x 19mm and Wristband Size: 169 x 19mm



Description	Product Code
ID BAND SPEEDI-PRINT INFANT IMPRINT WHITE (B/250)	PD510-11

Speedi-Print® Mother-Infant Wristbands

Soft vinyl material and Latex Free
 SecurSnap closure to prevent transfer
 ClearImage® surface protects handwriting and embossed information
 Serialized with alpha-numeric characters
 DuoBand - Set of 2 Wristbands - 1 Adult/Paed and 1 Infant with the same alpha-numerial serial number.



Information Area: 60 x 19mm and Adult/Paed Wristband Size: 237 x 25mm
 Infant Information Area: 47 x 17mm and Infant Wristband Size: 169 x 19mm

Description	Product Code
ID BAND SPEEDI-PRINT DUO-BAND WRITE-ON WHITE (B/150)	PD220S-11

Veri-Colour® Infant Wristbands

Soft vinyl material - Latex Free
 SecurSnap closure
 Insert cards are included (Booklet with 288 cards)
 Information Area: 52 x 14mm and Wristband Size: 170 x 18mm




NSW Contract

Description	Product Code
ID BAND VERI-COLOR INSERT INFANT CLEAR (B/250)	PD110-10
ID BAND VERI-COLOR INSERT INFANT WHITE (B/250)	PD110-11

Veri-Colour® Mother-Infant Wristbands

Soft vinyl material - Latex Free
 SecurSnap closure
 Available in 2, 3 or 4 part - each set has the same serial number - Serialized with alpha-numeric characters
 Insert cards are included
 PD220, PD320 and PD324 Series Adult/Paed Information Area: 60 x 19mm
 Adult/Paed Wristband Size: 237 x 25mm Infant Information Area: 47 x 16mm
 Infant Wristband Size: 169 x 19mm



Description	Product Code
ID BAND VERI-COLOR DUO-BAND INSERT WHITE (B/150)	PD220-11
ID BAND VERI-COLOR TRI-BAND INSERT WHITE (B/100)	PD320-11
ID BAND VERI-COLOR QUAD-BAND INSERT WHITE (B/100)	PD324-11

Conf-ID-Ent™ Newborn Bar Code Tags

For on-demand printing with a direct thermal printer. Unique, double-sided design and adhesive-free strap provides a flat bar code surface for easy access and scanning without disturbing the patient. Tags are made from soft material with rounded corners so they are gentle on your most sensitive patient's skin.

- Double-sided tag with strap allows scanning in any position, without disturbing the patient
- Easy-to-use and attaches to any wristband
- Durable, lightweight and comfortable design
- Supports text, bar codes, and photos
- Compatible with most direct thermal printers
- Patient ID and bar code on both sides for easy access and scanning

Product Number	Core Size	Wound	Qty./Box	Colour
WBTAG1XN-9	1"	In	500	White
WBTAG1XT-9	1"	Out	500	White
WBTAG15N-9	1-1/2"	In	500	White
WBTAG15T-9	1-1/2"	Out	500	White
WBTAG3XN-9	3"	In	1,000	White
WBTAG3XT-9	3"	Out	1,000	White



TenderCare® Thermal Mother/Father/Infant Wristbands

TenderCare® Thermal Newborn ID System

Positively identify newborn and parent(s) with ultra-soft thermal wristbands.

TenderCare® provides single-step, on-demand printing of bar-coded, ultra-soft, thermal wristbands to ensure that every newborn is properly identified and matched to their parent(s).

- Ultra-soft material for delicate newborn skin
- Tamper-resistant adhesive closure
- Fast, accurate, on-demand printing with a thermal printer
- Superior quality, seamless printing, and maximum resistance to water, soap, and hand sanitizer with the patented PermaPrint® surface
- Antimicrobial additive protects wristband surface against tested non-pathogenic bacteria*
- TenderCare® Thermal Mother/Father/Infant

Exclusive, patented 4-part wristband system positively identifies baby with mother and father or designated family

Wristband Length: Adult: 27.94cm / Infant: 20.3cm

Information Area Dimensions: Adult: 8.25cm x 2.5cm / Infant: 5cm x 2.2cm




Scanband® Z Infant Direct Thermal Wristbands

Soft and smooth infant wristbands are compatible with Zebra and other direct thermal printers.

Enjoy seamless printing with the patented PermaPrint® surface designed for durability and resistance to moisture, alcohol, and hand sanitizer to preserve patient data and bar codes.

Antimicrobial additive protects wristband surface against tested non-pathogenic bacteria*

Indexing bar is located at the perforation to ensure proper alignment and calibration



Precision® Thermal Bar Code ID Tags

Simple, print-on-demand tags turn any Mother/Father/Infant wristband* into a bar code wristband to improve patient safety and efficiency. Available in several formats to fit your application needs.

- Compatible with TenderCare® Infant ID Tether found below on this page.

Product Number	Description	Size	Core Size	Wound	Sets/Roll	Rolls/Box
PD7600-11	Mother/Father/Infant	Adult: 2.86cm L x 3.17cm H Infant: 2.86cm L x 1.6cm H	1.5"	In	200	2
PD7601-11	Mother/Father/Infant	Adult: 2.86cm L x 3.17cm H Infant: 2.86cm L x 1.6cm H	3"	In	200	2
PD7610-11	Infant Bar Code ID Tag	2.86cm L x 1.6cm H	1.5"	In	500	2
PD7612-11	Infant Bar Code ID Tag	2.86cm L x 1.6cm H	3"	In	500	2
PD7630-11	Paediatric Bar Code Tag	2.86cm L x 1.2cm H	1.5"	In	500	2



PD7600-11



PD7610-11



PD7630-11

Safeguard® Sealident® Insert Infant Wristbands

Tri-laminate material and latex free

Adhesive closure

Water-resistant Sealident insert pocket (insert pocket has flap with adhesive making insert pocket water resistant)

Insert cards are included



PD811-11—Information Area: 55 x 9mm and Wristband Size: 217 x 13mm

ID BAND SAFEGUARD INFANT INSERT WHITE (B/250)	PD811-11
---	----------

TenderCare® Infant ID Tether

Simple method to access identification tag for swaddled newborns

Ultra-soft and durable material for patient with sensitive skin.

Easy application plus tamper evident closure for patient safety.

Compatible with any wristband and ID tag combination.

Wristband Dimensions: 317mm x 11mm

Description	Product Code
INFANT TETHER - TENDERCARE WHITE (B/100)	PD7780-11
INFANT TETHER - TENDERCARE YELLOW DUCKS (B/100)	PD7780-82



Kleen-Print Disposable Foot Printer

Kleen-Print Foot Printer ensures that ink does not come into contact with baby, mother or nurse.

Designed with micro-thin film which prints on any type paper stock.

The foot printer includes 3 panels for Mother's Right Index, Infant's Left Footprint (or palm print) and Infant's Right Footprint (or palm print).



Description	Product Code
FOOT PRINTER - KLEEN-PRINT DISPOSABLE (B/25)	PD3100

SnugFit® Wristbands

Newborns typically lose water weight within the first 72 hours after birth.

SnugFit® features a unique adjustable closure to keep a "snug fit" at all times for safe and secure identification.

- Soft vinyl material with no rough edges
- Closure is re-adjustable to size needed
- Available in Imprinter and Insert styles (Insert cards included)
- Serialized with alphanumeric characters



Bundle of Joy® Wristbands and Stickers

Bundle of Joy wristbands

Proud family members deserve fun and colourful wristbands to celebrate baby's arrival and their new role. Wristbands are durable enough to last for days to let everyone know the exciting news. One size fits all. 250 wristbands per box.

Bundle of Joy Stickers

Cute and playful stickers announce baby's arrival! Made with glossy, high quality stock. 5cm diameter / 250 per roll



Breast Milk Labels

These secondary breast milk labels help communicate additional information such as whether the milk is fresh or thawed, date and time milk was thawed, and if any additives have been fortified in the milk.

BREAST MILK LABEL

INFANT NAME: _____

MRN: _____

COLLECTED DATE: _____ TIME: _____

THAWED: DATE: _____ TIME: _____

ADDITIVES: _____

BREAST MILK:

Fresh _____ Frozen _____

THAWING: Date: _____ Time: _____

ADDITIVES: Date: _____ Time: _____

Breast Milk

Freezer to Refrigerator Date _____ Time _____

Fortified No Yes Date _____ Time _____

By _____

INFANT'S LAST NAME:

Collected Date _____ Time _____

Freezer To Refrigerator Date _____ Time _____

Fortified No Yes Date _____ Time _____

Cal _____

BREAST MILK

BREAST MILK

BREASTMILK FOR 24 HOUR SUPPLY FEEDING

Patient _____

Date/Time Thawed _____

Signature _____

BREASTMILK FOR FREEZER STORAGE

Patient _____

Date _____

Signature _____

BREAST MILK

Time _____

Date _____

Expires _____

BREAST MILK

Fresh _____ Frozen _____

THAWING: Date: _____ Time: _____

ADDITIVES: Date: _____ Time: _____

Patient Chart, Drug Line and Nursing Labels

BED NO.

Smudge Proof Flex with Removable Adhesive
Size: 3.8 x 3.49cm Roll of 200

TMNPFBED-1 LABELS - "BED NO." WHITE #3.8X3.5CM #REMOV (R/200)
Only available in White

ROOM NO. PATIENT DOCTOR

Smudge Proof Flex with Removable Adhesive
Size: 13.65 x 3.49cm Roll of 200

TMNPF2286-** LABELS - PATIENT CHART PRE-PRINTED (COLOUR) #REMOV # 13.6 X 3.5CM (R/200)

TMNPF2286-** LABELS - PATIENT CHART BLANK (COLOUR) #REMOV #13.6X3.5CM (R/200)

** Colour code chart for TMNPF and TMNPFB

1 White	2 Yellow	3 Green	4 Pink	5 Red	6 Orange	7 Blue	8 Tan	9 Salmon
10 Copper	11 Rose	12 Pink	13 Lavender	14 Grey	15 Lime	16 Aqua	17 Chartreuse	20 Silver

SYRINGE MEDICATION

PATIENT _____

DRUG _____

AMT _____ /ml

PREPARED DATE _____ TIME _____

PREPARED BY _____

EXPIRES DATE _____ TIME _____

Code: TMSY7983

LABEL - "SYRINGE MEDICATION" YELLOW (R/333)

NAME ALERT

Code: TMFR1207

LKABEL - "NAME ALERT" FL RED (R/760)

ALLERGIES

Code: TMFR1435

LABEL - "ALLERGIES" FL RED (R/400)

MEDICATION ADDED TO I.V.

Time _____ Date _____ By _____

Drug _____ Quantity _____

Expires: Date _____ Time _____

Code: TMFR7329

LABEL - "MEDICATION ADDED TO I.V" FL RED (R/400)

MEDICATION ADDED

Solution _____ DRUG ADDITIVE _____ AMOUNT _____

ADDED BY: _____ DATE _____ TIME _____

PREPARED ROOM NO. _____ RATE _____ SPECIAL INFORMATION _____ ML/HR _____

DISCARD TIME _____ DATE _____

Code: TMFR17401

LABEL - "MEDICATION ADDED" FL RED (R/333)

ALLERGIC TO:

Code: TMFR1431

LABEL - "ALLERGIC TO" FL RED (R/420)

CAUTION EPIDURAL NARCOTIC

Code: TMFP9665

LABEL - "CAUTION EPIDURAL NARCOTIC" FL PINK (R/420)

Epidural BC2744

Code: TMFY59745608

LABEL - "EPIDURAL HIGH ALERT" FL YELLOW (R/1000)

CESAREAN

Code: TM59705025

LABEL - "CESAREAN" FL YELLOW (R/1000)

More Medication, Drug Line Labels and Nursing labels can be found in our Nursing Labels and Time Tape brochure reference number DBH2403

Seca Baby Electronic Scales

Seca 334 Baby Scale is designed as a reliable portable baby scale with a integrating weighing tray larger and is more comfortable than other models.

Seca 354 Baby Scale is designed as a dual function portable baby / child scale. At a touch of a button, the tray can be detached quickly and easily from the base leaving you with robust flat scales for children weighing up to 20 kg. As the two elements are securely connected with one another, it is absolutely safe to lift and transport the scales.

Seca 374 Baby Scale is designed as a dual function portable baby / child scale with extra large weighing tray. Seca 374 is also compatible with 360 wireless system - allowing measurements to be sent wireless to the wireless printer or to analytics 115 on your PC.



PN334

Model	Seca 334	Seca 354	Seca 374
Capacity & Graduation	20kg - 5gm up to 10kg and 10gm for 10kg plus	20kg - 10gm up to 10kg and 20gm for 10kg plus	20kg - 5gm up to 10kg and 10gm for 10kg plus
Functions	Automatic switch-off, TARE, Auto-HOLD, BMI Function	Breast Milk Intake, TARE, Auto-HOLD and BMI Function	TARE, Breast Milk Intake Function, auto-HOLD, auto-CLEAR, damping, automatic switch-off, SEND/PRINT
Weight & Dimensions	2.8kg - 63.8 x 30 x 10.5cm	2.3kg - 55.2 x 33.2 x 15.6cm	3.6kg - 62 x 35.8 x 19cm
Power	Batteries	Batteries	Batteries
Options	400 Power Adaptor, 232 Measuring Rod and 428 Carrying Case	413 Carrying Case	400 Power Adaptor, 233 Measuring Rod, 418 Head and Foot Positioners and 360 compatible

Description	Product Code
SCALE - DIGITAL BABY SCALE 20KG (EACH)	PN334
SCALE - DIGITAL BABY SCALE 20KG (EACH)	PN354
SCALE - DIGITAL BABY SCALE WITH X-LARGE TRAY 20KG (EACH)	PN374
SWITCHMODE POWER SUPPLY UNIT—SECA ADAPTOR FOR 703, 634, 869, 334 AND 374 (EACH)	PN400
OPTION - SECA MEASURING ROD FOR BABY SCALES (EACH)	PN232
CARRY BAG FOR BABY SCALES #334 (EACH)	PN428
TRANSPORT CASE—SECA FOR BABY SCALE 354 (EACH)	PN413
MEASURING ROD—SECA FOR BABY SCALE 374 #35-80CM (EACH)	PN233
HEAD AND FOOT POSITIONER—SECA FOR BABY SCALE 374 (EACH)	PN418



PN354



PN374

Seca Baby Measuring Mat

- Seca 210 Measuring Mat is kind to the skin, washable and foldable.
- The light weight, space-saving solution for easy, precise measurement of the length of babies and toddlers while lying down.
- The fixed head piece and the easy to slide foot positioner make it simple to use.
- Measuring range: 10 to 99cm
- Weight: 580g

Description	Product Code
MEASURING MAT - SECA BABY 10 - 99CM (EACH)	PN210



PN210

Baby Bassinet Trolley

The Baby Bassinet Trolley is available in 2 sizes Single or Twins.

- Mattress base height: 87.5 cm.
- Wheels: 7.5cm standard with 2 locking.
- Clear acrylic bassinet tub and 2 storage trays.
- Trendelenberg and Reverse Trendelenberg positioning.
- 5cm Mediflex covered mattress.
- Available in Powder Coated and Stainless Steel

BASSINET TRLY - OXFORD SINGLE S/S #80 X 42 X 106.5CM (EACH)	OX20-3815
BASSINET TRLY - OXFORD SINGLE P/C #80 X 42 X 106.5CM (EACH)	OX20-3812
BASSINET TRLY - OXFORD TWINS S/S #80 X 64 X 106.5CM (EACH)	OX20-3842
BASSINET TRLY - OXFORD TWINS P/C #80 X 42 X 106.5CM (EACH)	OX20-3840
BASSINETT MATTRESS - OXFORD SINGLE MEDIFLEX #66 X 32 X 5CM (EACH)	OX20-4100
BASSINET ACRYLIC TUB - OXFORD SINGLE #71 X 36 X 19-29CM (EACH)	OX20-4000
BASSINET TRAY - OXFORD WHITE PLASTIC #36 X 32 X 8CM (EACH)	OX20-4200
BASSINET ACRYLIC TUB - OXFORD TWINS (EACH)	OX20-4001
BASSINETT MATTRESS - OXFORD TWINS MEDIFLEX #66 X 50 X 5CM (EACH)	OX20-4101



Single Bassinette Bath Trolley

The Single Bassinette Bath Trolleys provides a cost effective alternative to Baby Baths without height adjustment.

- Bassinette features: clear acrylic twin tub with plug/waste drain hose and valve
- Tray feature: removable 2 x 390 x 350 x 80mm deep, white acrylic
- Tub height: 820mm high
- Manufactured from 304 stainless steel tube
- Castors standard: 75mm single rubber with 2x braking

BASSINET BATH TRLY - OXFORD SINGLE S/S W/ TUBING FOR BATH (EACH)	OX20-3850
--	-----------



Hi-Lo Baby Bath

The Hi-Lo Baby Bath Trolley is designed specifically for the purpose of giving baby baths without the need to empty heavy bathtubs and to be able to adjust the bath to a comfortable height. Electronically adjustable height from 63cm to 93cm

- Wheels: 10 cm standard with 2 locking.
- Clear acrylic bassinet single tub with drainage and plug.
- Controlled via Foot pedal control.
- Available in Powder Coated Only.

BABY BATH - OXFORD HI-LO (EACH)	OX20-3860
---------------------------------	-----------



Mobile Cots

The range of Mobile Cot are designed and manufactured to meet AS2130.

- Light weight side rails with squeeze grip handle for easy single handed operation.
- On each corner Socket with adjustable knob for either IV Poles / Detachable C Size Oxygen Bottle Holder can placed as required (IV Poles and Oxygen Bottle Holder is optional).
- Bottom shelf for additional storage.
- 10cm wall bumpers on each corners.
- Trendelenberg Angle Gauge
- Wheels: 12.5cm standard with 2 locking (diagonally positioned - expect for King Cot).
- Mediflex covered mattress—7.5cm high



OX20-3802

	Trendelenberg Adjust	Mattress Height	End Panels	Extras
Neo-Natal Cot	Both Ends Manual	86cm Height	Fixed end panels	
Premier Cot	Both Ends Manual	86cm Height	Fixed end panels	
King Cot	Both Ends Electrical	95.5cm Height	Fixed end panels	Central Locking Castors
Quattro Cot	Both Ends Manual	86cm Height	Removable end panels	
King Quattro Cot	Both Ends Electrical	95.5cm Height	Removable end panels	Central Locking Castors



OX20-3811



OX20-3800



OX20-3804

Cot Accessories

Include Telescopic IV Poles 76cm long with Ramshead and Detachable C Size Oxygen Bottle Holders.



OX13-3012

Description	Code
MOBILE COT - OXFORD NEO-NATAL P/C #102 X 68.5 X 145CM (EACH)	OX20-3802
MOBILE COT - OXFORD PREMIER P/C #143 X 79 X 145CM (EACH)	OX20-3800
MOBILE COT - OXFORD KING ELECTRIC P/C #143 X 79 X 152.5CM (EACH)	OX20-3804
MOBILE COT - OXFORD QUATTRO P/C #143 X 79 X 145CM (EACH)	OX20-3811
MOBILE COT - OXFORD KING ELECTRIC QUATTRO P/C #143 X 79 X 152.5CM (EACH)	OX20-3818

Products for Consideration

ECO Gel is a very good echo-transmission gel for ultrasound. Crystal-looking, it is salt-free and without string effect. It does not contain formaldehyde and, thanks to its particular structure, it helps the moving probe, thus making the operator's work easier.

Available in: 260ml bottle

<p>ECO Gel is: Push-vacuum-produced and it is therefore free from air microbubbles Non irritant Water-soluble Drying resistant Odourless Easily removable Alcohol Free</p>	<p>ECO Gel does not: Stain Grease Ionize Oxidize over time</p>
--	---



NSW Contract

ECO Gel does not contain any active ingredients that might damage the probe transducer. The used dyes and preservatives are non-toxic and strictly in accordance with Australian regulations.

ULTRASOUND GEL - GEL ECO CLEAR 260ML (EACH)	G008A
---	-------

Fresh+Clean 3-in-1 Air Freshener is a versatile surface cleaner, disinfectant and air freshener. It is designed for easy use. Initially designed and approved for use in Aircraft Fresh+Clean 3-in-1 Air Freshener can be used as both a room spray and surface spray. It also counteracts and nullifies unpleasant odours from the toilet and kitchen / laundry areas.

FEATURES & BENEFITS

- Non-flammable
- Non-toxic and non-corrosive formulation – will not cause corrosion
- Suitable for use in around food preparation areas
- Available in a range of regionally inspired fragrances, including natural oil and hypoallergenic formulations

APPLICATION AREAS

- Commercial grade disinfectant for use in galleys, kitchen facilities, toilets and air-conditioning vents to kill germs
- Effective surface cleaner for use in hand basins, toilets, cupboards, mirrors and hard surfaces.
- Air freshener and deodoriser eliminating odours and bad smells
- Can be used on surface areas that are in contact with food



Ki-San is a ready to use 2-in-1 high performance disinfectant cleaner for all hard surfaces. Kills viruses and bacteria, it provides broad spectrum disinfection in industrial, healthcare environments, shopping centres, schools, universities, and household settings.

ARTG listed **Ki-San** kills SARS-CoV-2 (COVID-19), Norovirus, Adenovirus, Influenza A (H1N1), Murine hepatitis, *Escherichia coli* (*E coli*), *Pseudomonas aeruginosa*, *Staphylococcus aureus* (*Golden Staph*), *Proteus vulgaris*, *K.pneumonia*, *E.faecalis* and *Acinetobacter*.

Ki-San will not harm or damage fabrics, plastic, metals, rubber or other surfaces. It is safe and suitable for use on high touch point areas such tabletops, door handles, light switches, desks, toilets and toilet doors, taps, TV remotes, kitchen surfaces and cupboard handles. Hard surface disinfectant only. Use undiluted. Do not mix with detergents or other chemicals. Not to be used on skin. Not to be used on medical devices or other therapeutic goods.



NSW Contract

Products for Consideration

Spirit Deluxe Dual Head Stethoscopes

Available in Adult 47mm, Paediatric 38mm and Infant 32mm (Head ø)
 Stainless Steel Dual Head Design with Floating Diaphragm and deep cone shaped bell
 Ultraflex™ binaural PVC tubing available in 20 colours.
 Includes spare diaphragm, retaining ring and soft ear tips.



Description	Product Code
STETHOSCOPE - SPIRIT DELUXE INFANT D/H BLACK (EACH)	STS607BLACK

Long Life Re-usable Polypropylene Bed Pan.

Ideal for patients confined to bed for:

- Medical
- Post surgery
- Geriatric
- Fall risk
- Maternity



NSW Contract

Description	Carton Qty	Product code
BEDPAN - SELINA SUPER-SLIPPER 335 X 245MM (EACH)	25	WASPN

Warwickware Holloware - Reuseable Polypropylene Products

Features.

- Strengthened sides with strong but lightweight design
- Manufactured from high quality polypropylene
- Autoclavable up to 135 ° C

Product Range includes:-

Gallipots,
 Jugs,
 Kidney Dishes,
 Lotion Bowls,
 Instrument Trays,
 Wash Bowls,
 Vomit Bowls.



WAKD200



WAIT1813



WAJG1000



WAWB3112



WAVB300

Aspire Happy Sleeper Chair

Seamlessly transforming from a snug lounge chair into a supportive sleeping surface, this chair caters to the needs of loved ones who wish to stay close to their hospitalised family members. The chair is upholstered in medical-grade silicone vinyl, meeting rigorous infection control standards while delivering unparalleled comfort and functionality.



SLEEPER CHAIR - ASPIRE HAPPY BEIGE VINYL (EACH)	AD198810
---	----------