



est. 1889

Aged Care Catalogue



Table of Contents

- 1. Featured Products..... 3
- 2. Mobility Equipment 5
- 3. Bedroom Products 7
- 4. Bathroom Products 10
- 5. Healthcare Carts & Trolleys 11
- 6. Medical Devices 18
- 7. Consumables 19
- 8. Infection Control 24
- 7. Patient Identification 26

About Us

Denyer Bros. Pty. Ltd was originally founded in 1889 by two brothers to manufacture surgical instruments. The business expanded from its original base in Melbourne with the Sydney operation becoming independently owned.

In 2000, Denyer Brothers re-named itself with the new trading name of DB Health. This change was to market & characterise the now broader scope of business activities and reflect the more current evolution in today's environment, compared to its original beginnings in 1889.

Today, DB Health remains a pivotal member in the health care community, with its experience and supply solutions strongly focused around patient safety. Quality products including Monitoring Electrodes, BetterShield Faceshield, and Electric Beds for Aged Care, GelFast Hand Sanitiser, Patient Identification Bands, Silent Knight Pill Crushers and Skin Markers reflect just a sample of the varied product range available.

DB Health is still one of the few Australian medical wholesalers to be 100% family owned and operated.

1. FEATURED PRODUCTS

Silent Knight Pill Crusher

The Silent Knight Pill crusher proven performer for over 10 years

- Eliminates the risk of cross contamination
- The pouches are designed to leave no residue behind.
- Quickly and evenly crushes tablets
- No mess and no clean up after



Code	Description	Qty
NS9999	Silent Knight Pill Crusher	Each
NS9997	Deluxe pouches	B/1000



The Grip-n-Assist

The Grip-n-Assist transfer belt makes every day living easier for patients and caregivers alike.

- One size fits all
- Sturdy, so slip grips
- Can be worn by patient or caregiver
- Can be worn backwards or forwards
- Four ergonomic grips provide versatility
- Patented support bones behind each handle
- Medicare and private health insurance reimbursable



Product Code	Description
LCPGNAW001	Grip N Assist Professional Belt #White (Each)

CoaguChek XS Plus PT/INR Meter

CoaguChek XS is designed to test Warfarin patients for PT/INR with results in less than 1 minute.

- Accurate, precise results in about one minute
- Small sample size (8µl) from a finger stick

Product Code	Description
BO06644511001	CoaguChek XS Plus Aged Care Kit (Each)
BO04625315017	INR Test Strips - PT CoaguChek XS (B/48)



BO04696522190	INR Control – PT CoaguChek XS (Each)
BO03603539195	Lancet - Accu-Chek Safe-T Pro Plus (B/200)

GelFAST anywhere

- sanitise hands instantly in between and during patient visits
- wearable & super convenient
- swappable refills
- kills 99.99% of germs
- made from natural plant based materials, orange scented and moisturising
- improves compliance with the 5 moments of hand hygiene



Product code	Description
MX11011	GelfastAnywear Starter Kit
MX2100	GelfastAnywear Refills
MX5101	GelfastAnywear Holster



Hydrasyn25

Hydrasyn25 is a specially formulated moisturiser that helps treat and prevent cracked heels, fissures, and severe dry skin conditions

- Pure silk amino acid softens and nourishes hard, dry, calloused heels and feet.
- Alpha Hydroxy Acid (6%) – helps exfoliate old, dry skin and allows for the growth of healthy, new, more pliable skin.
- UREA (25%) – penetrates deep below the surface, causing skin cells to attract, absorb, and hold moisture better.
- Non-greasy, fragrance-free formulation.

Product Code	Description
LMPHS25	Hydrasyn25 heel cream 150GM (Tube)



The Spirit Deluxe Stethoscope

Features stainless steel dual head design & dual ear soft tips, Superflex™ binaural tubing, Colour co-ordinated retaining & non-chill rings, Includes a spare Diaphragm & ear tips, Write-on name tag, Unique to the Spirit Deluxe range,

Code	Description	Qty



STS601-**	Deluxe Stethoscope Adult	Each
-----------	--------------------------	------

2. MOBILITY EQUIPMENT

Patient Lifting Equipment:

Designed and manufactured in Australia, each of these products states full compliance with Australian standards: ASA ISO 10535-2002 & ISO 10535:1988.

SONATA 150 is a lightweight compact mobile patient lifting hoist with innovative & ergonomic design and an easy to use every day lifter.

ALTO Lift 200 is a medium sized heavy duty mobile patient lifting hoist with a compact ergonomic design.

Code	Description
AESON-E	Patient Lifter- Sonata150kg Hoist w/ electric leg spread (Each)
AEALTO	Patient Lifter - Alto 200kg Hoist w/ electric leg spread (Each)



SALSA 200 is designed as an active standing aid allowing a single carer to raise a patient into a comfortable secure standing position. The SALSA is ideal for transferring from bed to chair, or chair to commode. Allied to suit the SALSA is a padded sling which allows superior access to remove clothing for toileting purposes.

- *Height & angle adjustable*
- *contoured kneo pad offers maximum lower limb support*
- *Removable footplate*
- *Patient friendly lift arms enable multiple grip positions during lift*



Walking frame folding

Features adjustable height, light aluminium frame, folds flat, foam handles, and 2.6kg unit weight

- Handle height: 780-930mm
- Base size: 650 x 500mm



Coopers walking frame folding

Features adjustable height, folds flat, contoured plastic hand grips and 2kg unit weight

Adult

- Handle height: 780-860mm
- Base size: 640 x 580mm

Tall:

- Handle height: 830-910mm
- Base size: 640 x 560mm



Self-propelled Combi Wheelchair

Features: Swing away removable leg rests, calf strap, removable armrests, quality steel folding frame, breathable padded nylon upholstery and rear backrest upholstery storage pocket

- Net weight: 14.5kg
- Max user weight: 100kg
- Dimensions: seat height 520mm; back height: 400mm;
- Push handle height: 910mm; overall width: 690mm



Vito Plus Transit Wheelchair

Features hand breaks, swingaway removable leg rests, heel loops, folding frame, quick release wheels, rear storage pocket, and breathable nylon upholstery

- Gross weight: 11.7kg
- Max user weight: 110kg
- Seat width x depth x height: 46 x 40 x 48cm



Product Code	Description
ALB202	Wheelchair - Care-Quip Vito Plus Transit (Each)

3. BEDROOM PRODUCTS

Motiva Electric Bed

The Motiva has been designed to provide a high quality and cost effective bed with an attractive resident room solution.

- Highly mobile with 125mm Castors and central locking
- 4 Section Patient Surface with Aero Surface
- Motor controlled function for backrest, thigh rest and height adjustment
- Height adjustment range of 36 – 78cm
- Hand control with lockable functions
- Available with 2 style of side rail: Twin Full Length (09) or Split 4-section (10)
- Max User Weight: 190kg
- Patient Surface area is 200 x 87cm with potential additional 20cm extension for tall patients
- 3 year Warranty



G10 Motiva



Bed Options: IV Poles, Urine Bag Holder, Bed Control Holder, lifting Pole and Hand Grip, Mattress Bolster for 20cm extension, Patient Surface Extension + Side rail Upholestry

D09 Motiva

Product Code	Description
WB10-0550	Bed – Motiva Electric (Each)

Carisma

The Carisma nursing home bed reflects the combination of cosiness and individual nursing comfort.

- Divided side protection
- Hybrid lifting column
- Motor controlled function for backrest, thigh rest and height adjustment
- Height adjustment from 36-80cm
- 2x2 castor breaks
- 225 kg safe working load
- 3 year warranty



Product Code	Description
WB10-0513	Bed – CarismaAge Care (Each)

Viskomatt Mattress


The ViskoMatt Memory Foam Mattress is a 2-layered mattress, upper layer made of 4 cm viscoelastic thermo-foam, lower layer made of 8 cm cold foam RG 40.

- Designed for patients under 150kg with medium risk.
- Cambio (Aqua Plus on one side and Terry Cloth on other w/ zipper) or Aqua Plus (waterproof w/ zipper) Covers
- Mattress Dimensions: 200 x 86 x 12cm

Product Code	Description
WB50-0085	Mattress - ViskoMatt Aqua Plus (Each)
WB50-0085CAM	Mattress - ViskoMattCambio (Each)



ViskoMatt Memory Foam



Cambio Cover

Bedside Cabinet - Attend

The Bedside Cabinet Attend is an Australian Made Cabinet that can be stained to match your décor.

- Available with 2 small drawers + 1 large drawer or in the configuration pictured.
- Available with or without drawer locks.
- Mounted on 5cm castors

Product Code	Description
CARVBS	Bedside Cabinet - Attend w/ 2 drawers and empty top shelf (Each)
CARVBS2	Bedside Cabinet - Attend w/ 3 drawers (Each)



Over bed Tables - Macs

Mac's Metal craft over bed table is manufactured from heavy tubular mild steel and has laminated table top with raised edges (only for tilt top)

- Easily height adjustable and simple to manoeuvre, the Mac's server table is not merely stylish.
- Tilt or standard top
- 4 castors of Ø50mm
- Overall dimensions (mm): 1160H (max) /790H (min) x 400W x 790L
- Available in various finishes, including laminated melamine, timber and vinyl wrap* with colour to match any existing furniture.

Product Code	Description
MACOB	Overbed Table W/ Tilt Top & C Base (Each)
MACOBN	Overbed Table W/O Tilt Top & C Base (EACH)



Overbed Table with Tilt Top & C Base

Deluxe Bed with Air Comfort Seating System

The Deluxe Bed is designed for use as a pressure reducing day bed.

- Total surface air pressure system.
- No gaps in seating surface.
- Adjustable seat and back (Infinite locking).
- Adjustable wings and drop-down arms
- (Multiple positions).
- Adjustable leg rest (infinite locking).
- Detachable/Adjustable head rest.
- Locking castors (Standard and Directional).
- Full metal frame.
- Replaceable cover.
- Options: Flip-open footrest (Adjustable), Meal tray, Incontinence cover, Seat belt, Harness, Buggy wheels or large castors.
- Outside Dimensions- Outside: Back width – 71cm, Seat width – 69cm, Sitting Length – 104cm.
- Inside Dimensions: Back width – 46cm, Seat width – 51cm, Back height – 71cm. Suit patients from 162cm – 188cm.



This product is available in various sizes along with Tilting, Electric and Bariatric versions.

Product Code	Description
AD198200	Day Bed – Air Comfort Deluxe (Each)

Nordic Mobicline Reclining Mobile Chair

The Nordic Mobicline Reclining Mobile Chair is versatile chair designed for use in Aged Care Facilities includes large activity tray and slide out footrest.

- Available in either fabric or vinyl with a variety of colour options.
- Fully welded tubular frame
- Dimensions: Overall width 70cm, Seat depth 49cm, Seat height (from footrest) 44.5cm, Seat to top of back 65cm, Seat to top of arm 19cm and Length (extended) 141cm



Product Code	Description
AD194700	Mobile Reclining Chair - Nordic Mobilcline #Std Colours (Each)

4. BATHROOM PRODUCTS

Kerry Mobile Shower Trolley

This Shower Trolley allows a user with limited mobility to shower and bathe with the assistance of a carer. Hydraulically operated fully height adjustable with dual pedals on both sides to suit a range of washing heights

- Stainless steel frame and push handle with aluminium bed frame
- 4 lockable castors for manoeuvrability and stability
- Dimensions: 185 x 80 x 50-95.5cm
- Weight: 63kg
- Max User Weight: 150kg



Product Code	Description
KH600	Shower Patient Trolley - Kerry Mobile P/C #Standard (Each)

Care Quip Mobile Shower Commode

The Mobile Shower Commode is designed to assist transportation of the patient / resident for showering and toileting.

- Available in durable stainless steel or powdercoated steel frame
- Padded commode seat (choice of open or closed)
- Sliding Footrest + swingback arms for side transfers
- PVC slatted backrest + 125mm castors w/ 2 locking
- Seat Dimensions (W x D x H): 44.5 x 43 x 59cm
- Overall width: 51cm
- Max User Weight: 120kg
- Options: Wide range of options including Folding Safety Arms, Lock Down Arms, IV Pole, Oxygen Bottle Holder, Pan Carrier, Bowl Carrier, etc.



Product Code	Description
ALB1020	Shower Commode - Mobile P/C (Each)

5. HEALTHCARE CARTS & TROLLEYS

Oxford Instrument & Dressing Trolleys

- Drawer Options: full width or half width drawers with roller bearer slides (half width only available on trolleys 80 & 90cm wide).
- Trolley Widths: From 50cm to 120cm (Standard Depth – 50cm and Height – 90cm).
- Castors: 7.5cm standard.
- Options: Drawer Locks, IV Poles, Sharp Container Brackets, etc

Product Code	Description
OXS18-0250	Instrument Trolley –50cm Stainless Steel (Each)
OXS18-0450	Dressing Trolley – 50cm Stainless Steel 1 drawer (Each)
OXS18-0650	Dressing Trolley – 50cm Stainless Steel 2 drawer (Each)
OXS18-0950	Dressing Trolley –50cm Stainless Steel 3 drawer (Each)



OXS18-0450



OXS18-0250

Oxford 30 Dosage Medication Trolley

Oxford 30 Dosage Medication Trolley is available in either powder coated or stainless versions.

- Central Locking + Wheels: 12.5cm standard with 2 locking
- Drawers: 3 x 10cm and 1 x 25cm deep
- Top 3 drawers have removable aluminium compartments with 3 dividers.

Product Code	Description
OX18-4560	Medication Trolley – Oxford 30 Dosage Stainless Steel (Each)
OX55-4560	Medication Trolley – Oxford 30 Dosage Powder Coat (Each)



OX55-4560

Kerry Webster Medication Trolleys

The Kerry Flexicart Webster Medication Trolleys are designed to be used with the Webster Medication System.

- ABS Plastic Removable Top
- Central Locking Drawers
- Wide Wheel Base with 125mm Castors
- Push Handle + Light weight Durable Construction
- Dimensions: KH220WEBSTER – 63.5 x 48.5 x 104.5cm
- Dimensions: KH222WEBSTER – 63.5 x 48.5 x 101.5cm
- Drawer Dimensions: 54 x 42 x 26cm
- Options: Folding Shelves, Waste Bins, Sharp Brackets, etc.



KH220WEBSTER + Folding Shelf + Waste Bin



KH222WEBSTER + Folding Shelf + Waste Bin

Product Code	Description
KH220WEBSTER	Medication Trolley - Kerry Webster Flexicart 3 Drawer P/C (Each)
KH222WEBSTER	Medication Trolley - Kerry Webster Flexicart 2 Drawer P/C (Each)

Oxford Linen Skips

Oxford Linen Skips are constructed from stainless steel and are available in 2 models with or without foot operated lid.

- Lid Colours: White, Red, Blue, Yellow and Green
- Linen Bags: Available in same colour as lids

Product Code	Description
OX26-1605	Linen Skip – Slimline Single w/o lid (Each)
OX26-1615	Linen Skip – Slimline Double w/o lid (Each)
OX26-1670	Linen Skip – Slimline Triple w/o lid (Each)
OX26-1645	Linen Skip – Slimline Single w/ foot operated lid (Each)
OX26-1655	Linen Skip – Slimline Double w/ foot operated lid (Each)
OX26-1680	Linen Skip – Slimline Triple w/ foot operated lid (Each)



OX26-1605 OX26-1670



OX26-1645 OX26-1680

Oxford Mobile Waste Bag Holder

Oxford Waste Holders are constructed of stainless steel with a powder coated lid

- Solid platform to support bags
- Incorporated clip to hold bags

Product Code	Description
OX10-1000	Waste Bag Holder – Oxford w/ foot operated lid (Each)



NS9000Y	Waste Bags - Yellow 70 X 54cm 40um LDPE (Each)
---------	--

Kerry Wet & Dry Laundry Trolleys

The Kerry Wet & Dry Laundry Trolleys are a mobile moist linen trolley made from a galvanized plus zinc rich epoxy coated frame.

- 125mm Heavy Duty Nylon Swivel Castors + Spring loaded bases
- Options: Push Handle for Medium & Large Sizes, and Spring Covers for Medium and Large Sizes.
- Dimensions: KH550B – 80 x 60 x 80cm
- Dimensions: KH551B – 105 x 70 x 85cm
- Dimensions: KH552B – 125 x 70 x 85cm



KH551B – Medium with Spring Loaded Base

Product Code	Description
KH550B	Wet & Dry Trolley - Kerry Small P/C w/ base (Each)
KH551B	Wet & Dry Trolley - Kerry Med P/C w/ base (Each)
KH552B	Wet & Dry Trolley - Kerry Large P/C w/ base (Each)

Rotoplas Linen Exchange Trolleys

The Rotoplas Linen Exchange Trolley™ has been innovatively designed for both the delivery of clean linen and collection of soiled linen.

- Two fixed, two swivel 150mm quality blue rubber castors.
- Fitted 3mm Galvanised base plate.
- Ergonomic cut out handle slots.
- Options: Vertical bolt on handles, Weather proof rip stop PVC cover, Tow Hitch for use with Electric Tug units, Hygienic Poly Shelves, Company Logo, Serial Numbering
- Dimensions: RTLET1 – 114 x 64 x 147cm
- Dimensions: RTLET2 – 105 x 71 x 145cm
- Dimensions: RTLET3 – 105 x 71 x 130cm



RTLET1 + Optional Shelf



RTLET2



RTLET3 + Optional

Shelves can be used to slot in front to form large bin or used as shelves.

Product Code	Description
RTLET1	Linen Exchange Trolley – Rotoplas #1 Standard (Each)
RTLET2	Linen Exchange Trolley – Rotoplas #1 Standard (Each)
RTLET3	Linen Exchange Trolley – Rotoplas #1 Standard (Each)
RTSLET1	Option – Poly Shelve (Each)

Rotoplas Rectangular Tubs Trolleys

The Rotoplas Rectangular Tubs are available for multiple uses from Laundry to Waste mobile trolleys

- Available in 3 separate parts: Rectangular Tub, Trolley, Rising Bases for 350L & 450L only
- Options: Lids and Tub Handle (bolted on)
- Dimensions: RTTTRT150 – 150L – 81 x 57 x 50cm
- Dimensions: RTTTRT250 – 250L – 115 x 63.5 x 67cm
- Dimensions: RTTTRT350 – 350L – 102 x 72 x 62.5cm
- Dimensions: RTTTRT350 – 450L – 132 x 72 x 62.5cm
- Add 13cm to the height if adding a trolley.

250L Tub has straight sides while the rest have tapered. The Rotoplas Rectangular Tubs are available for multiple uses from Laundry to Waste mobile trolleys



150L Tub



250L Tub



350L Tub



450L Tub

Product Code	Description
RTTTRT150	Rectangular Tub – Rotoplas 150L (Each)
RTTTRT250	Rectangular Tub – Rotoplas 250L (Each)
RTTTRT350	Rectangular Tub – Rotoplas 350L (Each)
RTTTRT450	Rectangular Tub – Rotoplas 450L (Each)
RTTTRT150	Trolley for Rectangular Tub – Rotoplas 150L (Each)
RTTTRT250	Trolley for Rectangular Tub – Rotoplas 250L (Each)
RTTTRT350	Trolley for Rectangular Tub – Rotoplas 350L (Each)
RTTTRT450	Trolley for Rectangular Tub – Rotoplas 450L (Each)
RTRBTRT350	Rising Base for Rectangular Tub – Rotoplas 350L (Each)
RTRBTRT450	Rising Base for Rectangular Tub – Rotoplas 450L (Each)

Bedpan and Bottle Holder

Stainless steel

Holds 4 or 6 bottles

Dimensions: 55-80x56cm high

Product Code	Description
OX30-1020	Bedpan and bottle rack- Oxford s/s.o 4ea



Oxygen Bottle Trolley

'C' size, 'D' size, 'G' size

Wheels: 15cm heavy duty

Product Code	Description
HQSP82C	Oxygen bottle holder size c mobile s/s (each)



HQSP82D	Oxygen bottle holder size d mobile s/s (each)
HQSP82G	Oxygen bottle holder size g mobile s/s (each)

IV Stands

On heavy cast alloy or nylon 65cm base

Telescopic height adjustment from 133-243cm

5 different style heads: 4 prong head, inside hook,

Outside hook, Ramshead, carol pole

Available in alloy or nylon base

Wheels: 5cm standard



Product Code	Description
OX10-3699	IV Stand- oxford 4 prong head w/ nylon base (ea)
OX10-3700	IV Stand- oxford o/s hook head w/ nylon base (ea)
OX10-3701	IV Stand- oxford ramshead w/ nylon base (ea)
OX10 3772	IV Stand- oxford carol head w/ nylon base (ea)

Chair Scale

Seca 952 chair scales for weighing while seated

Capacity of 200kg and graduation of 100 grams. Fold up arm rests, foot rests, wheel brakes, smooth running rollers, and a handle. TARE, auto-hold, and automatic switch off function.



Product Code	Description
PN954	Seca Scale digital chair 200kg (ea)

Wheel Chair Scale

Seca 676 electronic wheelchair scales with hand rail and transport castors

Stable, functional and mobile, patients can be weighed by sitting in a chair or wheelchair, and also ideal for weighing heavy persons. The scale can be folded to save space.

Capacity 360kg, graduation 100g, TARE, prew-TARE, auto hold, BMI, damping functions.



Scale with Hand Rails

Seca 644 multifunctional handrail scales



Low leveled platform and 300kg weight bearing capacity.
hand rail and a chair can be placed on platform if patient prefers
function memories detect storage weight such as a chair.

Integrated
to sit. 3



Seca 634 platform and bariatric scales

Capacity 360 kg. Graduation 100g, hold and auto-hold function.

Seca 750 robust mechanical floor scales

Stable and comfortable powdercoated scales.

Capacity, 150kg, graduation 1kg.



Office furniture

Wallaroo 2000

Wallaroo wall mounted workstations holds patient charts, electronic devices.

Has an electronic keyless lock and indicator lights signal when workstation is locked or open.



Deluxe Gas Lift Chair

Control seat height with the gas lift lever, raises from 53cm to 72cm

Adjustable height foot ring

Back rest fully adjustable from 10 to 16 degrees

Cast alloy base 65cm diameter

Available with non slip glides as an alternative to castors

Base available in alloy or nylon



For a full, free furniture catalogue: Free Call 1800 451 446

6. MEDICAL DEVICES

Spirit Sphygmomanometers

The Spirit Aneroid Sphygmomanometers are available in a Mobile Unit, Wall Unit and Portable Unit. These all come with a standard Adult size cuff but other sizes are available.

Product Code	Description
SP110	Sphyg - Spirit Aneroid Portable (Each)
SPCK141A	Sphyg - Spirit Aneroid Wall Mounted w/ Basket (Set)
SPCK142	Sphyg - Spirit Aneroid on Mobile Stand (Set)



CAS740 Vital Signs Monitor

The Casmed Monitor offers leading edge technology in a rugged, portable, lightweight unit.

Models available (with or without mobile stand and basket)

740-1 – MAXNIBP

740-2 – MAXNIBP + spO2

740-3 – MAXNIBP + spO2 + Temperature (SureTemp)

Product Code	Description
CAS740-2NL-M	Vital Signs Monitor - CAS 740 w/ Maxnibp&spO2 w/ Stand #240V (Each)
CAS740-3NL-M	Vital Signs Monitor - CAS 740 w/ Maxnibp, spO2 and temperature w/ Stand #240V (Each)



Accu-Chek Go Glucometer is designed for use on the go.

5 second test time, small 1.5µL blood sample, expired test strip warning, alarm clock function, target range indicator - warning indicator, hygienic test strip ejector, 300 test memory

Code	Description	Qty
BO03603806186	Accu-Chek Go includes meter, test strips & lancet	Kit
BO0353651386	Accu-Chek Go Test strips	B/50
BO03603539195	Accu-Chek Safe-T Pro lancets	B/200



Accu-Chek Safe-T Pro Lancet
- 3 adjustment depth settings

AccutrendPlus Meter is an all-in-one hand held device for measuring major CVD risk parameters which are blood glucose, cholesterol, troglycerides and lactate.

Hand-held point of care device, easy to use, fast results, holds up to 100 results can be stored with date, time, and additional information for each parameter.



Code	Description	Qty
BO05196418001	Accutrend Plus	Kit
BO11418262171	Accutrend Cholesterol test strips	B/25
BO11443054187	AccutrendGlucosel test strips	B/50
BO11538144016	Accutrend Triglyceride test strips	B/25

7. CONSUMABLES

Urinalysis Tests

Glucose Testing

Code	Description	Features
BO1443054	Glucometer test strips- Accutrend glucose	pack of 50

Kidney/Cardiovascular testing

Code	Description	Features
BO1544039	Micral Test	Results in 1 minute

Urine Testing

Code	Description	Qty
AM2300	Multistix 10 SG	Strips 10 SG
AM2301	Multistix 9	Strips 100 per vial
AM2304	Multistix 8 SG	Strips 100 per vial
AM2308	Multistix 5	Strips 100 per vial
AM2855	Uristix	Strips 100 per box



Multistix 10 SG

Faecal Screening

Faecal occult blood tests enhance detection of hidden blood in stool specimens

Code	Description	Features
CM800	Hemoccult SENSE	For colorectal cancer tests



Hemoccult SENSE

Catheters & Urinary Bags

Axis Urine Bags

Axis Urological Drainage Bags include the standard urine drainage bags A1 & A2, Urimaax Leg Bags and Urimaax Night Drainage Bags.

A1 & A2 Urological Drainage Bags:

- Latex Free + 120cm Regular Bore kink resistant tubing
- Stepped Catheter Connector + Anti-reflux valve
- Push-Pull bottom drain valve + Graduated front



A1 Urological Drainage Bags

Urimaax Leg Urine Bags

- Latex Free + Kink resistant tubing + Anti-Reflux valve
- Velcro Straps + Flip tap bottom drain
- 500-750ml Bag Size and tubing from 60-500mm



Urine Leg Bags

Urimaax Night Drainage

Product Code	Description
SA10013101	Urine Drainage Bag – Axis A1 2000ml #Non-Sterile (B/200)
SA10013102	Urine Drainage Bag – Axis A2 2000ml # Sterile (B/75)
SA10013200	Plastic Hanger – Axis #Fits A1 / A2 and Urimaax Night Drainage Bag (B/100)
SA10013110	Urine Drainage Bag – Axis Night 2000ml #Non-Sterile (B/250)
SA10013305	Urine Leg Bag – Urimaax 500ml w/ 50cm inlet tube #Sterile (B/80)
SA10013405	Urine Leg Bag – Urimaax 750ml w/ 50cm inlet tube #Sterile (B/80)



Urimaax Night Urine Bags

Warwickware – Urinals

- Comfortable
- Easy to clean
- Polypropylene material



WAUML1000



WAUF1000



WASPN

Code	Description	Size	Qty
WAUM1000	Male Urinal	1000ml	Each
WAUF1000	Female Urinal	1000ml	Each
WASPN	Bedpan	335 x 245mm	Each

Greatcare Vomit Bags

Disposable 1.5L Vomit / Emesis Bag with moulded plastic opening and twist lock closure.

Product Code	Description
GCG305001	Vomit Bag - Disposable Plastic Bag w/ Moulded Plastic Opening (B/50)



Zeomed Spill Kits

Zeomed spill kits have been designed to meet infection control and workcover guidelines. These Spill Kits are designed for Body Fluids, Biohazards and Cytotoxic Body Fluids.

Product Code	Description
ZF305	Spill Kit - Zeomed Biohazard (Kit)
ZF306	Spill Kit - Zeomed Body Fluid (Kit)
ZF307	Spill Kit - Zeomed Cytotoxic Body Fluid (Kit)
ZF304	Body Fluid Absorbent Powder 200gm Sachets (Each)



Contents of ZF306

Sharps Container

Product Code	Description
BD5487	Sharps container- bd guardian 1.4L
BD4588	Sharps container- bd guardian 3.1L



Safeskin* Synthetic

Powder Free Vinyl examination Glove. Length 240mm. Latex Free

Code	Size	Gloves/Box	Boxes/Case
KC55031	Small	100	B/10
KC55032	Medium	100	B/10
KC55033	Large	100	B/10



Proshield Loop & Soft Mask

The Proshield soft mask is ideal for clinicians working in general care or allied health environments where exposure to blood or body fluids is not a concern.

- » Soft hypoallergenic liner
- » Submicron filtration.
- » Triple pleat design for easy breathing.
- » Soft and durable earloops for Proshield Loop.

Code	Description	Quantity
SN2111444A	Proshield Soft Mask	B/50
SN211444A	Proshield Loop Mask	B/50



Hand Towels & Toilet Tissue

Medical Towel

Versa Towel is made from high quality, white 2 ply crepe wadding used across a spectrum of medical areas.

Code	Description	Size	Sheets/Roll	Qty/Case
KC4210	Versatowel-small roll	24.5 x 41.5cm	100	16 rolls
KC4220	Versatowel-large roll	49 x 41.5cm	100	8 rolls
KC4921	Dispenser	Suits KC4210	-	1 dispenser
KC4922	Dispenser	Suits KC4220	-	1 dispenser



Kleenex Facial Tissues & Deluxe Toilet Tissue

Code	Description	Qty/Pack	Qty/Case
KC4715	Facial Tissues	200	24
KC4735	Deluxe Toilet Paper 2ply	400 sheets	48 rolls



Dermatological solutions

No Rinse hair and body

Code	Description	Size	Qty
NR00200	No rinse shampoo	473ml	Each
NR00540	No rinse hair conditioner	250ml	Each
NR00910	No rinse body bath	473ml	Each



No rinse shampoo



Ego body wash



Selsun

Ego Body wash

Code	Description	Size	Qty
EG74500	Ego QV wash w/pump	500ml	Each

Medicated shampoo

Code	Description	Size	Qty
CH2525	Shampoo medicated Selsun	125ml	Each

Wound Care

Medium Crepe Bandages – Q Hospital

Q Hospital Bandages are a cotton crepe bandage designed for general purpose use. The bandage is made from a lightweight open fabric composition with non-fraying edges.

Code	Description
QB1605	Bandage- crepe hospital 5cm x 1.6m (Pkt 12)
QB1675	Bandage- crepe hospital 7.5cm X 1.6m(Pkt 12)
QB1610	Bandage- crepe hospital 10cm X 1.6m(Pkt 12)
QB1615	Bandage- crepe hospital 15cm X 1.6m(Pkt 12)



Elastic Adhesive Bandage – Tensoplast

Tensoplast is an elastic adhesive bandage that provides strong compression and long-lasting support. Manufactured completely from cotton, the bandage is ideal for treatment of all conditions that require limitation of movement.

Code	Description
SN1002	Bandage elastoplasts adhesive 5cm x 2.5m (roll)
SN1003	Bandage elastoplasts adhesive 7.5cm X 2.5m (roll)
SN1004	Bandage elastoplasts adhesive 10cmX 2.5m (roll)



8. INFECTION CONTROL

Disinfectants

Viraclean is a neutral detergent as well as a hospital grade disinfectant. Suitable for cleaning & disinfection of most surfaces

Code	Description	Size	Qty
WH008	Viraclean disinfectant	5Lt	Each
WH009	Viraclean disinfectant squeeze	500ml	Each



Triclosan- Microshield

Rapid disinfectant and hard surface cleaner

Code	Description	Qty
JJ60342	500ml spray pump	Each
JJ61363	5Lt	Each



Teriwipers have a number of applications eg. Surgeons scrub towel, tray liner & a general wiper. They are individually folded, are of scrim construction, maintain absorbency, have a strong wet strength, minimal lint generation & are easy dispensing.

Code	Description	Case	Total
KC4465	Teriwipers SML32 x 32cm	6pktsx75wipers	450 wipers



Isowipes are isopropyl alcohol impregnated bactericidal wipes that offer cost effective bactericidal disinfection of clean medical surfaces.

Code	Description	Size	Case	Total
KC6835	Isowipe Canisters	42 x 14.3cm	12x75wipers	900 wipes



Tuffie Wipes- detergent wipes are perfect for using anywhere where general cleanliness is required quickly and conveniently. They form an integral part of Infection Control practices.
Simple, effective way of removing dirt and grime

Code	Description	Per Canister
ZF225	Tuffie Wipes	225

MicroshieldRange

Pale pink liquid with cologne. Use as a rapid hand antiseptic, instant skin and hand disinfectant.

Code	Description	Size	Qty
JJ61356	MicroshieldHandrub	125ml	Each
JJ61357	MicroshieldHandrub	500ml	Each
JJ60238	MicroshieldHandwash	5Lt	Each
JJ60244	Microshield 2	5Lt	Each
JJ60239	Microshield 5 Handwash (5%chl.com)	5Lt	Each



Microshield handrub



Microshield handwash



Microshield 2 Chlorhexidine skin cleanser



Microshield 5 Chlorhexidine concentrate

7. PATIENT IDENTIFICATION

SoftGuard®

Soft tube design will not scratch delicate skin. Includes a clear sleeve design to protect patient information. Insert cards included.

Product Code	Material	Closure Type	Band Size	Label Size	Bands per Box
PD190-10	Soft tube design	Adhesive	Adult 25.4cm (I)	16.5cm x 0.7cm	50

NOTE: Available in clear (10) **ONLY**

Veri-Colour®



This band is Australia's most popular. It is made from soft vinyl for comfort and utilises a plastic sleeve design for insert cards. Cards included.

Product Code	Material	Closure Type	Band Size	Insert Size	Bands per Box
PD140-**	Soft vinyl	Adhesive	Adult 33cm (I)	8.8cm x 2.2cm	250
PD120-**	Soft vinyl	Adhesive	Adult/Pediatric 33cm (I)	5.2cm x 2cm	250
PD110-**	Soft vinyl	Adhesive	Infant 16.5cm (I)	5.2cm x 1.4cm	250

11

16

NOTE: Also available in clear (10) and alternative colours

PD140-16



Contact us for a full range of adult, short stay, RFID, bar code, and alert identification wristbands

Free Call 1800 451 446



est. 1889

Contact Us:

Phone: (02) 9638 2100

Free Call: 1800 451 446

Fax: (02) 9638 2700

Email: customerservice@denyer.com.au

Website: www.dbhealth.com.au

Address:

PO Box 7182 Silverwater NSW 1811

Unit 17, 38-46 South Street, Rydalmere
NSW 2116

DB1712