


# Holloware

Product range for Medical Holloware includes:

Warwick Gallipots, Lotion Bowls, Wash Bowls, Kidney Dishes, Instrument Trays, Compartment Instrument Trays, Jugs, Bedpans, Slipper Bedpans, Urinals, Vomit Bowls and Denture Cups, Drinking / Feeding Cups, Mesh / Perforated CSSD Instrument Trays, Silicone Instrument Protection Pads, Instrument Tray Tags, Scope Containers, Disinfection Sets (Soaking), Ident-Alert CSSD Tags and SteriLite Holloware.



Suppliers of quality products to the Healthcare Industry est 1889.  
Australian Owned & Operated.

Web [dbhealth.com.au](http://dbhealth.com.au)

Email [customerservice@denyer.com.au](mailto:customerservice@denyer.com.au)

Freecall: 1800 451 446

ABN: 37 000 459 174

## Warwick Polypropylene Reusable Holloware

- Long Life Re-usable Polypropylene Holloware is Designed for up to 1000 washes and sterilisations -
- Lightweight, strong and durable.
- Improved removal of blood residue compared to Stainless Steel.
- Grip feature moulded below the rim to improve holding.
- Smooth edges and surfaces eliminate paper rips during packing and sterilization.
- Note: scouring and abrasion will reduce the performance of products.
- Can be sterilised by Autoclave, Ethylene Oxide or Hydrogen Plasma and Autoclavable up to 135C
- CE Mark – Class 1 Medical Device

### Warwick Gallipots and Lotion Bowls

- Non drip with strengthened design for secure grip and anti vacuum pips on the bases.

Description	Carton Qty	Product Code
GALLIPOT 35ML #40 X 30MM (EACH)	100	WAGP40
GALLIPOT 95ML #60 X 40MM (EACH)	300	WAGP60
GALLIPOT 280ML #80 X 60MM (EACH)	80	WAGP80
LOTION BOWL 280ML #100 X 45MM (EACH)	200	WALB100
LOTION BOWL 900ML #150 X 70MM (EACH)	60	WALB150
LOTION BOWL 2L #200 X 90MM (EACH)	30	WALB200
LOTION BOWL 3L #250 X 105MM (EACH)	20	WALB250



WALB100



WAGP60

### Warwick Kidney Dishes

- Non drip with strengthened design for secure grip and graduations

Description	Carton Qty	Product Code
KIDNEY DISH 300ML 150 X 75 X 35MM (EACH)	210	WAKD150
KIDNEY DISH 500ML 200 X 98 X 45MM (EACH)	150	WAKD200
KIDNEY DISH 750ML 250 X 130 X 55MM (EACH)	60	WAKD250
KIDNEY DISH 1.5L 300 X 144 X 65MM (EACH)	25	WAKD300



WAKD150

### Warwick Kidney Dishes

- Non drip with strengthened design for secure grip and graduations

Description	Carton Qty	Product Code
WASH BOWL 3L 310 X 118MM (EACH)	25	WAWB3112
WASH BOWL 5L 350 X 120MM (EACH)	25	WAWB3513



WAWB3513

## Warwick Instrument Trays

The Warwickware Autoclavable Instrument Trays are available in a range of sizes with the following features.

- Strengthened sides with strong but lightweight design with features such as — Flat and Ribbed bases, Smooth edges and corners.
- Small trays suitable for ophthalmic, chiropody and dental sets
- Large trays suitable for general use in surgical sets including heavy orthopaedic.
- WAIT1913 and WAIT2715 come with hole for Tray Tracability Tag.
- Can be sterilised by Autoclave, Ethylene Oxide or Hydrogen Plasma and Autoclavable up to 135C.
- CE Mark – Class 1 Medical Device
- Designed for up to 1000 washes and sterilisations - scouring and abrasive cleaning will reduce the performance of products.



WAIT3025



WAIT1813



WAIT1915

Description	Carton Qty	Order Code
INSTRUMENT TRAY 181 X 134 X 54MM #PLAIN BASE (EACH)	40	WAIT1813
INSTRUMENT TRAY 195 X 135 X 23MM #PLAIN BASE (EACH)	115	WAIT1913
INSTRUMENT TRAY 200 X 150 X 51MM #RIBBED BASE (EACH)	60	WAIT2015
INSTRUMENT TRAY 200 X 150 X 51MM #PLAIN BASE (EACH)	60	WAIT2015PB
INSTRUMENT TRAY 270 X 150 X 30MM #PLAIN BASE (EACH)	60	WAIT2715
INSTRUMENT TRAY 300 X 250 X 52MM #RIBBED BASE (EACH)	25	WAIT3025
INSTRUMENT TRAY 300 X 250 X 52MM #PLAIN BASE (EACH)	25	WAIT3025PB
INSTRUMENT TRAY 420 X 305 X 75MM #RIBBED BASE (EACH)	12	WAIT4030
LID FOR INSTRUMENT TRAY #IT2715 (EACH)	10	WAITL271

## Warwick Compartment Instrument Trays

The Warwickware Autoclavable Compartment Instrument Trays are available in a range of sizes with the following features.

- Strengthened sides with strong but lightweight design with features such as— Smooth edges and corners, Integrated gallipots to maximise space for instruments and accessories.
- Can be sterilised by Autoclave, Ethylene Oxide or Hydrogen Plasma and Autoclavable up to 135C.
- CE Mark – Class 1 Medical Device
- Designed for up to 1000 washes and sterilisations - scouring and abrasive cleaning will reduce the performance of products.



WACIT1814



WACIT2718

Description	Carton Qty	Code
COMPARTMENT INSTRUMENT TRAY 181 X 134 X 54MM (EACH)	38	WACIT1814
COMPARTMENT INSTRUMENT TRAY 205 X 130 X 30MM (EACH)	120	WACIT2013
COMPARTMENT INSTRUMENT TRAY 270 X 180 X 41MM (EACH)	25	WACIT2718
COMPARTMENT INSTRUMENT TRAY - LIGHT 220 X 150 X 24MM (EACH)	160	WALCT2215

## Warwick Jugs

Long Life Re-usable Polypropylene Jugs are designed for up to 1000 washes and sterilisations.

- Clear, easy to read graduated measurements.
- Lightweight, strong, durable and easy to handle.
- Can also be washed in automatic washer disinfectant or by hand using detergents and disinfectants.
- Can be sterilised by Autoclave, Ethylene Oxide or Hydrogen Plasma and Autoclavable up to 135C.
- CE Mark – Class 1 Medical Device



WAJG1000



WAJG2000

Description	Carton Qty	Product code
JUG NATURAL 500ML (EACH)	50	WAJG500
JUG NATURAL 1000ML (EACH)	45	WAJG1000
JUG NATURAL 2000ML (EACH)	40	WAJG2000

## Warwick Holloware Spa

- Spa Products are designed for a lengthy service - guaranteed for 1000 washes and sterilizations.
- Scouring and abrasive cleaners will reduce the performance of products.
- Can be cleaned in automatic washer disinfectors or suitable for hand washing using detergents and disinfectants.

This range encompasses Urinals, Bedpans and Drinkware.

### Warwick Spa Bedpans

Long Life Re-usable Polypropylene Bed Pan, Ideal for patients confined to bed due to illness or after surgery

- Lightweight, strong and durable - with a weight tolerance up to 250kg (WASPN has a max user limit of 160kg).
- Can be cleaned in automatic washer disinfectors or hand washed. Both WAVPN and WASPN can be autoclaved.
- Smooth material surfaces ensure ease of use and cleaning
- CE Mark – Class 1 Medical Device and Autoclaveable.

Description	Carton Qty	Product code
BEDPAN - WARWICK BUDGET 370 X 310 X 90MM (EACH)	30	WABPN
BEDPAN - VECTOR W/ LID 335 X 308 X 108MM (EACH)	12	WAVPN
BEDPAN - SELINA SUPER-SLIPPER 335 X 245MM (EACH)	25	WASPN
BEDPAN - SLIPPER LARGE 449 X 309MM (EACH)	12	WAFPN



WABPN



WASPN



WAFPN



WAVPN



## Warwick Spa Male and Female Urinals

- Re-usable Polyethylene Female and Male Urinals
- Comfortable and easy to clean.
- Manual washing (maximum 90C) and processing in Washer Disinfector Machines.

WAUF1000



WAUMC1000

Description	Carton Qty	Product code
FEMALE URINAL - CYGNET BOTTLE W/ HANDLE 1000ML (EACH)	30	WAUF1000
URINAL MALE - TOPPER BOTTLE W/ HANDLE 1200ML (EACH)	45	WAUMC1000

## Warwick Spa Vomit Bowls and Denture Cups

- Re-usable Polypropylene Vomit Bowl and Denture Cups
- Manual washing (maximum 90C) and processing in Washer Disinfector Machines.

Description	Carton Qty	Product code
VOMIT BOWL W/ HANDLE 300 X 65MM 1500ML (EACH)	25	WAVB300
DENTURE CUP W/ HINGED LID & LIFT OUT STAINER (EACH)	75	WADNT



WAVB300



WADNT

## Warwick Drinking and Feeding Cups

- Easy to use, clean and disinfect, meaning both patient and carer can be confident that products are safe to re-use.
- Beakers and Drinking Cups are wide, easy to grip, non-scratch, drop-proof and have a non-drip rim.
- Handled beakers are ideal for patients with unsteady or infirm hands.
- All Spa products can be cleaned in automatic washer disinfectors compliant with ISO 15883. Also suitable for hand washing using detergents and disinfectants.



WABFC250



WABFH250

Description	Carton Qty	Product code
BEAKER FEEDER / DRINKING CUP 250ML (EACH)	300	WABFC250
BEAKER FEEDER / DRINKING CUP 250ML W/ HANDLES (EACH)	140	WABFH250
LID W/ NARROW SPOUT #USE WITH BFH250/BFC250 #4X4MM SPOUT (EACH)	300	WABNS
LID W/ MEDIUM SPOUT #USE WITH BFH250/BFC250 #7X4MM SPOUT (EACH)	300	WABMS
LID W/ WIDE SPOUT #USE WITH BFH250/BFC250 #13X8MM SPOUT (EACH)	300	WABWS
LID RIM FEEDER #USE WITH BFH250/BFC250 (EACH)	300	WABRF



WABNS



WABMS



WABWS





WABRF



# Warwick Mesh and Perforated Instrument Trays

Manufactured from a high quality polypropylene material that is approved for use in medical devices.

- CE Mark – Class 1 Medical Device
- Low Cost, Robust and hard wearing—Guaranteed for 1000 cycles
- 50% lighter than Stainless Steel & easier to clean off blood protein
- Can be sterilised using Autoclave, ETO or Hydrogen Plasma
- Smooth corners & surfaces reduce the risk of piercing wraps.
- Perforated or mesh base for improved water penetration and steam circulation,
- Select trays also have perforations in the side walls
- Strengthened sides
- Tag holes for tray tags.



Description	Carton Qty	Product Code
INSTRUMENT TRAY - MESH 200 X 180 X 30MM (EACH)	50	WAMIT1515
INSTRUMENT TRAY - MESH 200 X 150 X 51MM (EACH)	60	WAMIT2015
INSTRUMENT TRAY - MESH 275 X 110 X 30MM (EACH)	60	WAMIT2709
INSTRUMENT TRAY - MESH 270 X 150 X 30MM (EACH)	60	WAMIT2715
INSTRUMENT TRAY - MESH 280 X 170 X 55MM (EACH)	25	WAMIT2817
INSTRUMENT TRAY - MESH 300 X 250 X 52MM (EACH)	24	WAMIT3025
INSTRUMENT TRAY - MESH 424 X 305 X 75MM (EACH)	12	WAMIT4030
INSTRUMENT TRAY - PERFORATED 300 X 250 X 52MM (EACH)	25	WAPIT3025
INSTRUMENT TRAY - PERFORATED 424 X 305 X 75MM (EACH)	12	WAPIT4030
DENTAL INSTRUMENT TRAY - 255 X 255X 50MM (EACH)	20	WADINT2525
DENTAL INSTRUMENT TRAY - 480 X 255 X 50MM (EACH)	25	WADINT4825



## Warwick Silicone Protection Pads and Rolls

- Designed to protect instruments during washing, sterilisation and transportation.
- Pliable pad with perforations for maximum water penetration and steam circulation with long fingers for instrument positioning.

Description	Product Code
SILICONE PAD - PINK #230 X 130MM (B/6)	WASP2313
SILICONE PAD - PINK #280 X 230MM (B/6)	WASP2823
SILICONE PAD - PINK #270 X190MM (B/6)	WASP2719
SILICONE ROLL - PINK 1M #230 X1000MM (EACH)	WASR10



## Warwick Tray Tags with S/S Rings

- Designed to enable tracking of trays .
- Tray Tags are autoclavable with area for barcode or label.




Description	Product Code
ATTACHMENT RING FOR TRAY TAG S/S (EACH)	WATR30
TRAY TAG - WARWICK 105 X 40MM WHITE (EACH)	WATAG010
TRAY TAG - WARWICK 105 X 40MM YELLOW (EACH)	WATAG015
TRAY TAG - WARWICK 105 X 40MM BLUE(EACH)	WATAG020
TRAY TAG - WARWICK 105 X 40MM LIGHT GREEN(EACH)	WATAG022
TRAY TAG - WARWICK 105 X 40MM ORANGE (EACH)	WATAG025
TRAY TAG - WARWICK 105 X 40MM RED(EACH)	WATAG030
TRAY TAG - WARWICK 105 X 40MM GREEN (EACH)	WATAG050
TRAY TAG - WARWICK 105 X 40MM GREY (EACH)	WATAG055
TRAY TAG - WARWICK 105 X 40MM LAVENDAR(EACH)	WATAG060
TRAY TAG - WARWICK 105 X 40MM TAN(EACH)	WATAG065
TRAY TAG - WARWICK 105 X 40MM BLACK (EACH)	WATAG070

## Warwick Container for Long Instruments and Scopes

- Designed to process up to 3 long instruments at a time.
- Each Container can also be stacked on each other
- Mesh base, sides & lid allow water penetration and steam circulation.
- Moveable silicone grips to hold instruments securely in place.
- Rapid water drainage with raised floor and internal curves.
- Positioning lock rib to centralise lid.
- Handle area designed to help lid removal.

Description	Product Code
DISINFECTION SET - CONTAINER + LID FOR RIGID SCOPES / LONG INSTRUMENTS (EACH)	WASMT5510





## Warwick Disinfection Sets

- Suitable for decontamination / manual sterilization of instruments and equipment by liquid soaking & manual cleaning.
- Can be used to disinfect both sterile and non-sterile equipment.
- Strong and durable - container resistant to most hospital cleaning/disinfectant chemicals.
- Recommended for low temperature sterilization processes.
- Smooth edges and surfaces eliminate paper rips during packing and sterilization.
- 4L / 12L Sizes can be used for cleaning Transoesophageal Echocardiography, Transvaginal Ultrasound and Transrectal Ultrasound Probes.



WADSSW1000



WADSS2711



WADSS4215

Description	Product Code
DISINFECTION SET - WARWICK BOX W/ STRAINER AND LID WHITE 1L #26 X 13 X 6.5CM (EACH)	WADSSW1000
DISINFECTION SET - WARWICK BOX W/ STRAINER AND LID BLUE 3L # 25.7 X 18 X 12.8CM (EACH)	WADSS2711
DISINFECTION SET - WARWICK BOX W/ STRAINER AND LID BLUE 4L # 42.6 X 15.6 X 13.8CM (EACH)	WADSS4215
DISINFECTION SET - WARWICK TANK W/ STRAINER AND LID BLUE 12L # 89 X 20 X 14.5CM (EACH)	WADSS8515



WADSS8515

## Ident-Alert CSSD Tags

- Easy to apply to tray handles, scopes, or surgical instruments
- Durable and lab-tested to withstand autoclave and common disinfectant solutions
- Cost-effective and versatile for a variety of applications
- Multipurpose Tags—Write critical ID information or alerts directly on the tag with a permanent marker or apply autoclavable labels or autoclave tape
- Message Tags—Special alerts communicate tray type and status
- Scope Tags—Printed with helpful fields to streamline and document the entire scope sterilization process

Description	Product Code
CSSD TAGS - IDENT-ALERT MULTIPURPOSE WHITE (CTN/1000)	PDTSTAG-01
CSSD TAGS - IDENT-ALERT MULTIPURPOSE BLUE (CTN/1000)	PDTSTAG-02
CSSD TAGS - IDENT-ALERT MULTIPURPOSE GREEN (CTN/1000)	PDTSTAG-03
CSSD TAGS - IDENT-ALERT MESSAGE PINK BROKEN (CTN/1000)	PDTSTAG-04
CSSD TAGS - IDENT-ALERT MESSAGE BLUE MISSING (CTN/1000)	PDTSTAG-05
CSSD TAGS - IDENT-ALERT MESSAGE GREEN VENDOR (CTN/1000)	PDTSTAG-06
CSSD TAGS - IDENT-ALERT MESSAGE YELLOW LOANER (CTN/1000)	PDTSTAG-07
CSSD TAGS - IDENT-ALERT MESSAGE RED PRIORITY (CTN/1000)	PDTSTAG-08
CSSD TAGS - IDENT-ALERT SCOPE WHITE DATE & SERIAL (CTN/1000)	PDTSTAG-09
CSSD TAGS - IDENT-ALERT SCOPE WHITE CLEANED BY (CTN/1000)	PDTSTAG-10
CSSD TAGS - IDENT-ALERT SCOPE WHITE SCOPE REPAIR (CTN/1000)	PDTSTAG-11





## SteriLite Holloware

The SteriLite Holloware Range is a selection of single use polypropylene holloware. Design to speed up packing and sterilising processes. SteriLite is ideal for use in Sterile Pack Manufacturing, Decontamination Services, Operating Rooms and Clinics due to its low cost, product design and performance.

- Low cost and disposable
- No need to invest time and money into retrieval and decontamination after use.
- Clean and ready for sterilisation
- Can be sterilised by steam autoclaving at 135C, Ethylene Oxide (EO) and Hydrogen Peroxide (H2O2).
- Can be used as standalone items or as a supplementary component.
- Widely used in sterile services for packing instrument sets and procedure kits
- Used non-sterile in hospital wards and Clinics.

### SteriLite Compartment and Packing Trays



Description	Product Code
PACKING TRAY—STERILITE 43X33X6.4CM (CTN/63)	WAST4330
COMPARTMENT TRAY—STERILITE 4 SECTIONS 22X15X2.4CM (CTN/160)	WALCT2215S

### SteriLite Gallipots



WAGPS60



WAGPS120



WAGP100

Description	Product Code
GALLIPOT—STERILITE 60ML 6.6 X 3CM (CTN/2000)	WAGPS60
GALLIPOT—STERILITE 120ML 8.5 X 3CM (CTN/1000)	WAGPS120
GALLIPOT—STERILITE 250ML 10 X 5.3CM (CTN/750)	WAGP100
GALLIPOT—STERILITE 500ML 13 X 6.6CM (CTN/735)	WAGP130



WAGP130

## SteriLite Sponge and Theatre Bowl



WASB170



WASB200

Description	Product Code
SPONGE BOWL—STERILITE 1000ML 17X8.4CM (CTN/340)	WASB170
SPONGE BOWL—STERILITE 2000ML 21X9CM (CTN/220)	WASB200
THEATRE BOWL—STERILITE 5000ML 35.2X11.6CM (CTN/50)	WATB350



WATB350

## SteriLite Kidney Dishes



WAKD250



WAKD300

Description	Product Code
KIDNEY DISH—STERILITE 750ML 25.2X10X4.8CM (CTN/450)	WAKD250
KIDNEY DISH—STERILITE 1500ML 30X14X5.7CM (CTN/135)	WAKD300

---

## SteriLite Jug and Quiver

---



WAJ1000



WAQV30

Description	Product Code
JUG—STERILITE 1L (CTN/100)	WAJ1000
QUIVER—STERILITE 19 X 6.5CM (CTN/75)	WAQV30