

Gynaecological Products

Product range for Gynaecological Products includes:

Pregnancy Tests, Uriscreen UTI Test.

Intra-Uterine Devices, Lubrication Gel, Disposable Vaginal Specula, Reusable Vaginal Specula, Vaginal Glass Dilators, Pessary Rings, Ultrasound Gel and Sterile Ultrasound Gel



Suppliers of quality products to the Healthcare Industry est 1889.
Australian Owned & Operated.

Web dbhealth.com.au

Email customerservice@denyer.com.au

Freecall: 1800 451 446

ABN: 37 000 459 174

Prgenancy Test - Bio Cel Xcel kit

BIO-CEL XCEL PREGNANCY TEST

Rapid immunoassay for the qualitative detection of Human Chronic Gonadotropin (hCG) in urine for the early detection of pregnancy.

- 2nd generation, one step pregnancy test
 - Rapid results in 3 minutes
 - Easy to read colour
 - Built-in test control
 - Room temperature storage and long shelf life
 - Self contained reagents - no mixing
- Single or batch testing - easy and rapid.



Description	Product Code
PREGNANCY TEST - BIO-CEL EXCEL KIT (B/35)	MEBIO

UriScreen UTI Rapid Urine Test

Rapid Screening Test for detection of Urinary Tract Infection.

- Result: 1 - 2 minutes
- High sensitivity and high predictive rate
- Detects both bacteriuria and pyuria
- No instrumentation needed
- Room temperature storage and long shelf life
- User friendly and inexpensive



Description	Product Code
URINE SCREENING TEST - URISCREEN (B/20)	MEURI

Intra-Uterine Devices

Load 375 IUD—This device is surrounded by a copper wire coil of 375mm and fitted with a monofilament Nylon string at the base of the device. The load 375 IUD is available in standard size only.

IUD Thread Retriever—The Emmett IUD Thread Retriever is designed to retrieve the lost threads if they are within the endocervical canal or uterine cavity. The procedure is simple, painless, atraumatic – generally no analgesia nor local anaesthesia is required – and can be carried out in a Doctor's surgery or Clinic setting. This saves the patient time and inconvenience and may avoid the need for hysteroscopy.



MEPEL400565



MELOAD375

Description	Order Code
IUD - LOAD 375 (EACH)	MELOAD375
IUD THREAD RETRIEVER - PELICAN STERILE (B/25)	MEPEL400565

L-Gel Surgical Lubricant

L-Gel Surgical Lubricant is suitable for personal and medical / surgical use. This lubricant is water soluble, non staining, greaseless and non-irritating. Available in 500gm pump pack, 100gm tube or 3 gm sterile sachet.



ME1002



ME1004

Product	Product Code
L-GEL SURGICAL LUBRICANT STERILE 3 GRAM SACHET (B/150)	MG1002
L-GEL SURGICAL LUBRICANT 100GM TUBE (EACH)	MG1003
L-GEL SURGICAL LUBRICANT 500GM PUMP PACK (EACH)	MG1004

Med Gyn Vaginal Specula with LED Light Disposable Sterile

The MedGyn disposable Vaginal Speculum is made with clear polystyrene for easy, one-hand operation and insertion. Features cordless LED illumination for improved visualization of exam area and greater convenience.

- Bright LED light with built-in battery – no cord to tangle.
- Rounded, smooth edges ensure patient comfort.
- Durable plastic construction.
- Easily locks in several positions to facilitate examination and save time.



Product	Product Code
VAGINAL SPECULA—MEDGYN DISP W/ LED LIGHT SML (B/50)	MEMG020225
VAGINAL SPECULA—MEDGYN DISP W/ LED LIGHT MED(B/50)	MEMG020226
VAGINAL SPECULA—MEDGYN DISP W/ LED LIGHT LGE (B/50)	MEMG020227

KleenSpec® Disposable Vaginal Specula

- Clear acrylic allows complete visualization of the vaginal canal
- Wide, ergonomic handle improves balance for easier manipulation during exams
- Available in a broad range of sizes, for all types of patients
- Look for the color-coded Welch Allyn handles
- Smooth acrylic will not pinch or bind tissue for improved patient satisfaction
- Strong grip and yoke designed to be more flexible than low-cost substitutes




4 Models pictured from X-Small, Small, Medium and Large.

Description	Product code
VAGINAL SPECULUM - PREMIUM 590 DISPOS SML (B/24)	WE59000
VAGINAL SPECULUM - PREMIUM 590 DISPOS MED (B/24)	WE59001
VAGINAL SPECULUM - PREMIUM 590 DISPOS LRG (B/18)	WE59004

Livingstone EOS Disposable Vaginal Specula's

- Made from clear plastic
- With easy to use and simple lever operation
- Disposable and for single use to avoid cross infections
- Provides maximum strength and tested up to 16lbs pressure
- Individually packed in hygienically sealed plastic
- Non toxic, pyrogen free and FDA(USA) approved



Descriptions	Product code
SPECULUM VAGINAL - E.O.S. DISPOSABLE PLASTIC SMALL (PKT/10)	LVSPECS
SPECULUM VAGINAL - E.O.S. DISPOSABLE PLASTIC MEDIUM (PKT/10)	LVSPECM
SPECULUM VAGINAL - E.O.S. DISPOSABLE PLASTIC LARGE (PKT/10)	LVSPECL

Reusable Vaginal Specula's

Liberty Pederson Reusable Vaginal Specula's.

Liberty Surgical Instruments are made from the finest French and Japanese stainless steel using German machinery and technology. This gives Liberty® the finish that professional people are looking for with the added benefit of low cost.

Liberty® instruments are made in Pakistan to German specifications, giving you a quality instrument at a lower price.




Descriptions	Product code
VAGINAL SPECULA - LIBERTY PEDERSON SML (EACH)	RH894
VAGINAL SPECULA - LIBERTY PEDERSON MED (EACH)	RH893
VAGINAL SPECULA - LIBERTY PEDERSON LGE (EACH)	RH892

Vaginal Glass Dilators

Made from custom blown glass and comes in 6 sizes.
Can be purchased individually or in a set of all 6 sizes.

Size	A	B	C	R
1	22mm	80mm	25mm	6.5
2	24mm	85mm	30mm	7mm
3	26mm	95mm	32mm	8mm
4	30mm	105mm	34mm	8mm
5	34mm	115mm	36mm	9mm
6	38mm	125mm	38mm	10mm

Description	Product Code
VAGINAL DILATOR GLASS - SET OF #1 - 6 (SET)	GG113
VAGINAL DILATOR GLASS #1 (EACH)	GG107
VAGINAL DILATOR GLASS #2 (EACH)	GG108
VAGINAL DILATOR GLASS #3 (EACH)	GG109
VAGINAL DILATOR GLASS #4 (EACH)	GG110
VAGINAL DILATOR GLASS #5 (EACH)	GG111
VAGINAL DILATOR GLASS #6 (EACH)	GG112



Pessary Rings

High quality flexible PVC pessary ring

- Latex free
- 12.5mm thick PVC ring
- Size ranges from 50mm to 110mm

Description	Product Code
PESSARY PORTIA VINYL RING 50MM (EACH)	PA366-50
PESSARY PORTIA VINYL RING 53MM (EACH)	PA366-53
PESSARY PORTIA VINYL RING 56MM (EACH)	PA366-56
PESSARY PORTIA VINYL RING 59MM (EACH)	PA366-59
PESSARY PORTIA VINYL RING 62MM (EACH)	PA366-62
PESSARY PORTIA VINYL RING 65MM (EACH)	PA366-65
PESSARY PORTIA VINYL RING 68MM (EACH)	PA366-68
PESSARY PORTIA VINYL RING 71MM (EACH)	PA366-71
PESSARY PORTIA VINYL RING 74MM (EACH)	PA366-74
PESSARY PORTIA VINYL RING 77MM (EACH)	PA366-77
PESSARY PORTIA VINYL RING 80MM (EACH)	PA366-80
PESSARY PORTIA VINYL RING 85MM (EACH)	PA366-85
PESSARY PORTIA VINYL RING 90MM (EACH)	PA366-90
PESSARY PORTIA VINYL RING 95MM (EACH)	PA366-95
PESSARY PORTIA VINYL RING 110MM (EACH)	PA366-110



Ultrasound Gel

ECO Gel is a very good echo-transmission gel for ultrasound. Crystal-looking, it is salt-free and without string effect. It does not contain formaldehyde and, thanks to its particular structure, it helps the moving probe, thus making the operator's work easier.

ECO Gel is:

- Push-vacuum-produced and it is therefore free from air microbubbles
- Non irritant
- Water-soluble
- Drying resistant
- Odourless
- Easily removable
- Alcohol Free

ECO Gel does not:

- Stain
- Grease
- Ionize
- Oxidize over time

ECO Gel does not contain any active ingredients that might damage the probe transducer. The used dyes and preservatives are non-toxic and strictly in accordance with Australian regulations. Available in: 260ml bottle



Product Code	Description
G008A	ULTRASOUND GEL - GEL ECO CLEAR 260ML (EACH)

Sterile Ultrasound Gel

Sterile Ultrasound Gel satchels are used in vaginal and intrauterine examinations, surgeries, ultrasonography applications for in vitro fertilization and biopsy taking. It is in a sterile disposable package of 20 ml. Sterilized by gamma ray.

ARTG Number: 346886

Product Code	Description
G020M	ULTRASOUND GEL – KONIX STERILE 20ML SATCHET (CTN/240)

