


Laundry and Storage

Product range for Medical Carts and Trolleys includes:

- Laundry and Waste carts and Trolleys,
- Rectangular Tubs
- Linen Exchange Trolleys
- Wire Baskets
- Shelving



Suppliers of quality products to the Healthcare Industry est 1889.
Australian Owned & Operated.

Web dbhealth.com.au

Email customerservice@denyer.com.au

Freecall: 1800 451 446

ABN: 37 000 459 174

Clean Linen Carts

The Clean Linen Trolley has been designed to move large amount of clean linen as required by healthcare staff.

- Light weight construction.
- Three shelves with 30cm clearance
- Push handle at both ends.
- Castors: 12.5cm standard swivel.
- Available in Stainless Steel or Powder Coated.

CLEAN LINEN TROLLEY - OXFORD S/S 3 LVL #113 X 60 X 110CM (EACH)	OX26-1300
CLEAN LINEN TROLLEY - OXFORD P/C 3 LVL #113 X 60 X 110CM (EACH)	OX55-1300



OX26-1300

Mobile Linen Storage Carts

- Strong, lightweight fully welded construction
- Shelf clearance: 350mm (top shelf 1440mm high)
- Shelf feature: 4x sloping shelves at 2° and 3x lower shelves with centre partitions
- Perforated sheet metal on 3 sides
- Optional vinyl or clear PVC cover with zippers at the front (sold separate)
- Manufactured from 304 stainless steel or from mild steel with zinc coating
- Push handle
- Castors standard: 200mm heavy duty with 2x fixed swivel/2x braking and bumpers

MOBILE LINEN STORAGE S/S 4 LVL #145 X 61 X 140CM (EACH)	OX23-3114
MOBILE LINEN STORAGE P/C 4 LVL #145 X 61 X 140CM (EACH)	OX55-3114



OX55-3114

Tallboy Linen Exchange Trolley

Linen Exchange Trolley has been innovatively designed for both the delivery of clean linen and collection of soiled linen in healthcare facilities. The design of the Linen Exchange Trolley allows use for bulk collection and linen storage.

Model TLET1 Specifications:

- Tallboy Linen Exchange Trolley 1 with 3mm Galvanised Base Plate and 150mm Castors (2 fixed and 2 swivel) — 114 x 64 x 170cm
- Poly Shelf / Bin Front 1125 x 600 x 15mm
- Available in Blue, Yellow, Green and Red (Standard Colours—a wide range of colours are available for larger purchases).
- Options Available: Rip Stop PVC Clip-on Cover, Vertical Stainless Tubular Handle, Polyethylene Shelf/Bin Front, Powder Coated Wire Shelf, Stainless Steel Wire Shelf, Document Holder, Tow Hitch Kit and Galvanised Hanging Rail, custom marking options such as logo prints / serial number, etc.

These LET are made to order a number of models are available with different configurations to suit different uses, or can be customised if you have a specific request please contact customer service.

LINEN EXCHANGE TROLLEY— TALLBOY #1 W/ 2 POLYSHELF & GAL BASE PLATE AND 15CM WHEELS #115 X 64 X 170CM (EACH)	RTP-LET1-SET1
---	---------------



TLET1 with PVC Cover



TLET1 with Shelves to form bulk collection



TLET1 with Poly shelves for

Oxford Linen Skips

The Linen Skips are available in Without Lid models - Single, Double or Triple and Foot Operated Lid - Single Double or Triple.

- Rectangular open ring accommodates a standard drawer string linen bag
- Base of skip designed to support the linen bag on removal.
- Wheels: 7cm standard single
- Available in Stainless Steel Only.

Linen Skip Bags Specifications:

- Impermeable water-resistant sl300 machine washable polyester
- Pull drawcord with cord lock
- 5 colours— White, Red, Blue, Yellow and Green



OX26-1681



OX26-1645



OX26-1605



OX26-1655



OX26-1615



OX26-1680



OX26-1670

LINEN SKIP - OXFORD SINGLE W/O LID #35 X 45 x 87CM(EACH)	OX26-1605
LINEN SKIP - OXFORD DOUBLE W/O LID #35 X 45 x 87CM(EACH)	OX26-1615
LINEN SKIP - OXFORD TRIPLE W/O LID #35 X 45 x 87CM(EACH)	OX26-1670
LINEN SKIP - OXFORD SINGLE W/ FOOT OP. LID #38 X 59.5 X 87CM (EACH)	OX26-1645
LINEN SKIP - OXFORD DOUBLE W/ FOOT OP. LID #38 X 59.5 X 87CM (EACH)	OX26-1655
LINEN SKIP - OXFORD TRIPLE W/ FOOT OP. LID #38 X 59.5 X 87CM (EACH)	OX26-1680
LINEN SKIP BAG - OXFORD WHITE (EACH)	OX26-1660
LINEN SKIP BAG - OXFORD RED (EACH)	OX26-1661
LINEN SKIP BAG - OXFORD BLUE(EACH)	OX26-1662
LINEN SKIP BAG - OXFORD YELLOW (EACH)	OX26-1663
LINEN SKIP BAG - OXFORD GREEN (EACH)	OX26-1664
LINEN SKIP REPLACEMENT LID - OXFORD SLIMLINE (EACH)	OX26-1681

Kerry Wet and Dry Laundry Trolley

The Kerry Wet and Dry Laundry Trolley is designed for use in institutional laundries.

- Galvanised plus rich zinc epoxy powder coated.
- Spring loaded base designed to minimise bending. Base hangs from 4 springs which are attached to the top of the trolley.
- Light weight.
- Full wraparound PVC bumper.
- Moist linen is able to breathe.
- Wheel: 12.5 cm Single Swivel on Small and 15cm Single Swivel on Medium and Large Trolleys.



KH550B—Wet & Dry Trolley with Spring Base

WET & DRY LINEN TROLLEY - KERRY SMALL P/C W/ BASE #80 X 60 X 80CM (EACH)	KH550B
WET & DRY LINEN TROLLEY - KERRY MEDIUM P/C W/ BASE #105 X 70 X 85CM (EACH)	KH551B
WET & DRY LINEN TROLLEY - KERRY LARGE P/C W/ SPRING BASE #125 X 70 X 85CM (EACH)	KH552B

Oxford Waste Bag Holders

The Waste Bag Holders have been designed as a mobile waste receptacle. Easy to operate.

- Designed for use with wide plastic bags. Bag holder dimensions: 35 x 23.5cm
- Solid platform to support bags.
- Incorporated clips holds any bag securely.

Waste Bag Size - NS9000Y is a short bag - the base does not touch the base plate.
Recommended Bag Size - 70cm wide x 90cm deep.

WASTE BAG HOLDER - OXFORD W/FOOT OP. LID S/S #28 x 41 x 81CM (EACH)	OX10-1000
WASTE BAGS - YELLOW 70 X 54CM 40UM LDPE (B/500)	NS9000Y



OX11-1000

Rectangular Tub Moist Linen Trolleys

Rotoplas Rectangular Tub Moist Linen Trolley is designed for easy of use when a significant volume of Moist Linen is removed for cleaning. The moist linen trolley is made up of two main parts a 3mm galvanised dolly base with fitted wheels (2 fixed and 2 swivel 125mm rubber) and a polyurethane bin.

Options available include rising base:

- Coil Spring Rising Base with Poly Top - Spring Capacity of 25, 50, 75 or 100kg (also available with galvanised top).
- Hanging Spring Rising Base
- Bolt on Handles and Lids.




Powdercoated Handle



RTP-SRT250D

RECTULANGULAR TUB—150L TAPERED SIDES W/ MOBILE BASE #81 X 57 X 65CM (EACH)	RTP-TRT150D
RECTULANGULAR TUB—250L STRAIGHT SIDES W/ MOBILE BASE #115 X 63.5 X 67.5CM (EACH)	RTP-SRT250D
RECTULANGULAR TUB—340L STRAIGHT SIDES W/ MOBILE BASE #99 X 63.5 X 82CM (EACH)	RTP-SRT340D
RECTULANGULAR TUB—450L STRAIGHT SIDES W/ MOBILE BASE #123 X 73 X 82CM (EACH)	RTP-SRT450D
OPTION—COIL SPRING BASE W/ POLYTOP FOR 150L RECTANGULAR TUB (EACH)	RTP-CSBTRT150P
OPTION—COIL SPRING BASE W/ POLYTOP FOR 250L RECTANGULAR TUB (EACH)	RTP-CSBSRT250P
OPTION—COIL SPRING BASE W/ POLYTOP FOR 340L RECTANGULAR TUB (EACH)	RTP-CSBSRT340P
OPTION—COIL SPRING BASE W/ POLYTOP FOR 450L RECTANGULAR TUB (EACH)	RTP-CSBSRT450P



Hanging Spring Rising Base



Coil Spring Base w/ galvanised top

Wire Baskets

MEDIMESH provides innovative storage solutions for your health facility. The wire shelving can be quickly and easily changed to meet your needs. The range of functional accessories offers a vast variety of configuration possibilities.

MEDIMESH requires little or no tools for assembly.

MEDIMESH design minimises dust collection on the surface and easy to wipe down whilst maintaining excellent airflow. Static shelving, mobile shelving and wall framed storage for hanging wire mesh baskets are available.

These storage solutions 304 Stainless Steel.



HQSSBSK413021



HQSSBSK413010



Basket wall frame HQSSBRT



Product code	Description
HQSSBRT	BASKET WALL FRAME - S/S #485X485MM #7 RAILS (EACH)
HQSSBSK413010	BASKET - S/S MEDIMESH #410X300X105
HQSSBSK413021	BASKET - S/S MEDIMESH #410X300X215
HQSSBSK414510	BASKET - S/S MEDIMESH #410X450X105
HQSSBSK414521	BASKET - S/S MEDIMESH #410X450X215
HQSSDIV3010	DIVIDER - S/S MEDIMESH #300X100MM (EACH)
HQSSDIV3020	DIVIDER - S/S MEDIMESH #300X215MM (EACH)
HQSSDIV4510	DIVIDER - S/S MEDIMESH #450X100MM (EACH)
HQSSDIV4520	DIVIDER - S/S MEDIMESH #450X215MM (EACH)

Chrome Wire Baskets

FIT ON LOUVRE PANELS OR SHELVES

Versatile baskets to replace the classic plastic bin, can hook on existing Louvre panels, sit on a bench or back-to-back on shelves

HIGH QUALITY CHROME WIRE

Manufactured with high quality chrome wire ensures a high quality and durable basket

UNIQUE SOFT EDGED DESIGN

Unique soft edged design - removes harsh edges - won't damage sterile stock and reduces potential hazards to employees

EASY CLEANING

Wire basket ensures dust falls through and doesn't collect in the base, allowing for easier cleaning

VERSATILE FOR STORING GOODS

Can be used in any industry to stock goods or products - hospitals, food & beverage, retail, commercial, warehouses


DIVIDERS & STACKING FRAMES AVAILABLE

Dividers available to adapt to the storage requirements and stacking frames allow baskets to be placed on top of each other

Chrome Wire Baskets

Ideal for storing sterile and non-sterile stock, protective equipment, consumables and other items


- Allows items to be visible and organised
- Customise with labels to identify goods stored
- Use stacking frames to stack baskets on top of each other for space saving
- Available in 4 sizes: small, medium, large and extra large (see back for dimensions)



Chrome Wire Bed baskets

Ideal for storing patient files, protection equipment, wound drains, monitoring equipment, small consumables and other items


- Attach securely to patient's bed or mount to wall
- Dividers available to customise and better organise items stored in basket
- Available in 5 sizes: wide in small, medium and large and tall in small and large (see back for dimensions)



Chrome Catheter Wire Basket

Protect and store catheters or other long equipment (endoscopes etc) from being damaged

- Can be placed on shelving unit or attach to louvre panel
- Ensures greater organisation of catheters
- Dividers available to customise items stored in basket
- Available in 3 sizes: horizontal small and large and vertical (see back for dimensions)



Small Maxi Basket				Medium Maxi Basket		Large Maxi Basket	Extra Large Maxi Basket
457mm	1	2	3	1	2	1	1
	4	5	6	3	4	2	2
	7	8	9				
914mm	10	11	12	5	6	3	3
	13	14	15	7	8	4	4
	16	17	18	9	10	5	
1372mm	19	20	21	11	12	6	5
	22	23	24				6
	25	26	27	13	14	7	
1829mm	28	29	30	15	16	8	7
	31	32	33				8
	34	35	36	17	18	9	

Product Codes	Description
NACWB DIVS	Divider for Large Chrome Wire Maxi Basket 245mm x 155mm
NACWB DIVSXL	Divider for Extra Large Chrome Wire Maxi Basket - 425mm x 210mm
NACWB PANS	Plastic Louvre Panel - 460mm x 460mm
NACWB DIVBBW	Divider for Wide Bed Baskets - 134mm x 143mm x 138mm
NACWB DIVBBH	Divider for High Bed Baskets - 134mm x 143mm x 238mm
NACWBSSBBS	Stainless Steel (202) Bracket for Wall Mounting (2 required for each basket)
Stacking frames for Chrome Wire Maxi Baskets	
NACWB20FR	Small stacking frame to suit Small Chrome Wire Maxi Basket
NACWB30FR	Medium stacking frame to suit Medium Chrome Wire Maxi Basket
NACWB40FR	Large stacking frame to suit Large Chrome Wire Maxi Basket
NACWB60FR	Extra Large stacking frame to suit Extra Large Chrome Wire Maxi Basket
NACWB20	Small Chrome Wire Maxi Basket - 225mm x 135mm x 125mm
NACWB30	Medium Chrome Wire Maxi Basket - 235mm x 205mm x 160mm
NACWB40	Large Chrome Wire Maxi Basket - 250mm x 420mm x 160mm
NACWB60	Extra Large Chrome Wire Maxi Basket - 400mm x 420mm x 210mm
NACWB15	Small Wide Bed Basket - 150mm x 300mm x 130mm
NACWB35	Medium Wide Bed Basket - 150mm x 500mm x 130mm
NACWB45	Large Wide Bed Basket - 150mm x 600mm x 130mm
NACWB55	Small Tall Bed Basket - 150mm x 350mm x 250mm
NACWB65	Large Tall Bed Basket - 150mm x 500mm x 250mm
NACWB25	Catheter Chrome Wire Basket Horizontal Small - 400mm x 130mm x 135mm
NACWB70	Catheter Chrome Wire Basket Horizontal Large - 600mm x 130mm x 135mm
NACWB80	Catheter Chrome Wire Basket Vertical - 150mm x 130mm x 600mm

Mobile Storage Trolleys and Stationary Shelving Units

Louvre Panel Trolley

Designed to hold wire or plastic baskets and includes bottom shelf with Stainless Steel construction.

Dimensions: 950w x 600d x 1200h

Shelves: 1 x Shelf

Castors: 4 x 100mm Non-Braked

Panels: 4 x 45.7 x 91.4cm Louvered Panels (2 mounted on each side).

Product code	Description
EMSP619-1	LOURVE PANEL TROLLEY S/S (EACH)
NACWBPAN-S	LOURVE PANEL - METAL SMALL #45.7X45.7CM (EACH)
NACWBPAN-M	LOURVE PANEL - METAL MEDIUM #45.7X91.4CM (EACH)
NACWBPAN-L	LOURVE PANEL - METAL LARGE #45.7X182.8CM (EACH)



Suture Mobile Trolley or Shelving—Available in a number of sizes and shelves with or without castors

122 x 46 x 139 or 159 or 174cm (W x D x H)

Each of the angled shelves has a height of 11cm. Usually these trolley have 4-5 levels.

For Mobile Trolley—2 swivel/directional castors, 2 brake castors are included.

Please contact DB Health regarding the size of Suture Shelving or Mobile Suture Trolley you are interested in and how many shelves you may wish to have.



* **Sliding baskets shelving unit and trolleys** are available in a number of configurations. The standard kit on feet (no castors) - 46cm x 122cm x 159cm (D x W x H) with 8 sliding baskets (2 sizes available—54 x 41 x 25cm or 54 x 41 x 18cm)

Castors can be added as an option.

Please contact DB Health regarding the type of sliding basket shelving you are interested in and any options you may wish to have added to it and which size of baskets you wish have included.



A variety of shelving units are available with or without casters!

Depending on your available space, height and other requirements, these units can be configured to your needs.

Standard stainless steel shelving units available also with either 3, 4 or 5 tier / shelves, these shelves come in lengths of 1000mm, 1250mm 1500mm or 1750mm and in 2 sizes either 450mm or 600mm shelf depth.

Height of the shelving units will depend on the space you have available. Additional tiers available on request depending on height of unit


If you require mobile shelving units, keep in mind an additional height of 170mm for the caster is required.

Please contact us with your requirements.

Oxford Powdercoat Colours

Below are the colours for Oxford Products along with white gloss. Please colours are approximate representations of the final powdercoating colour.

When selecting colour please advise drawer colour and body colour.



Standard Powdercoat Colours



Non-Standard Powdercoat Colours (additional cost may apply)

Flexicart Powdercoat Colours

Below are the standard colours for Flexicart Medication Trolleys along with White and Signal Red.

