

Sterilisation Equipment and consumables

Product range for Sterilisation Equipment and Consumables includes: Sterilisers, Sterile Pouches, Sterilisation Reels, Sterilisation Indicator Tapes, Chemical Indicator, Cleaning Brushes.



Suppliers of quality products to the Healthcare Industry est 1889.
Australian Owned & Operated.

Web dbhealth.com.au

Email customerservice@denyer.com.au

Freecall: 1800 451 446

ABN: 37 000 459 174

Sterilisers

Siltex Pratika B20 (20L) Printer

Able to convert from Type S to Type B sterilizing at the touch of a button, the Siltex Pratika B has been designed to combine all the benefits of a Class B sterilizer with the convenience of choice.

The Pratika B has pre-vacuum for more efficient sterilizing, fully programmable sterilizing and drying cycles, two chamber sizes to match your needs, a Bowie-Dick test cycle and clean and waste water tanks.

The features of this steriliser are:

- Type B or Type S sterilizing.
- 20 litre chamber.
- 121°C or 134°C sterilizing cycles.
- Closed door drying cycle.
- Digital temperature readout.
- User selectable cycle lengths.
- Printer option available.
- Thermal cut-out safety switches.
- External waste water dump tank 7 litre chamber.
- Inbuilt SD card for electronic backup and SD card reader.

Dimensions: 45 (W) x 39 (H) x 58cm (D).

Weight: 45kg.

Chamber capacity: 20 Litre - 24cm Diameter x 43cm (L) with 5 litre external dump tank.

240V, 10A power outlet required.

Sterilizer comes with circular rack system and 4 trays, tray handle, 2 wire racks, dump tank and technical manual.

Printers require one additional power point, and connect directly to the sterilizer.



Siltex Baby Sterilizer

The Siltex Baby sterilizer has been designed for a medical standard sterilizer that is simple to use, compact and affordable. It is a downward displacement steam sterilizer, with a closed door drying cycle and an optional printer. The features of this steriliser are:

- 7 litre chamber.
- 121°C or 134°C sterilizing cycles.
- Closed door drying cycle.
- Digital temperature readout.
- Analogue pressure gauge.
- User selectable cycle lengths.
- Printer option available.
- Thermal cut-out safety switches.
- External waste water dump tank.

Dimensions: 37 (W) x 39 (D) x 30cm (H).

Weight: 20kg.

Chamber capacity: 17cm Diameter x 27cm (L) with 5 litre external dump tank.

Standard 240V, 10A power outlet required.

Accessories includes circular rack system and 2 trays, tray handle, dump tank and technical manual.



Product Code	Description
SI910	STERILISER - SILTEX PRATIKA B20L W/ PRINTER (EACH)
SI710	STERILISER - SILTEX BABY 7L W/ PRINTER (EACH)
SI750650	OPTION - RACK FOR SILTEX STERILISER AC3D (EACH)

Medipack Sterilisation Pouches

The Medipack Self-Sealing Pouches are suitable for use in Steam and EO Gas Sterilisation.

- Sealed by hand.
- One-side clear film for easy identification of contents and the other medical grade paper.
- Integrated sterilisation indicator on medical grade paper side.



Product Code	Description
MEAP499	STERILISATION BAG - MEDIPACK 7 X 23CM (B/300)
MEAP500	STERILISATION BAG - MEDIPACK 9 X 13.5CM (B/300)
MEAP501	STERILISATION BAG - MEDIPACK 9 X 23CM (B/300)
MEAP502	STERILISATION BAG - MEDIPACK 13.5 X 26CM (B/200)
MEAP503	STERILISATION BAG - MEDIPACK 19 X 33CM (B/200)
MEAP504	STERILISATION BAG - MEDIPACK 23 X 36.5CM (B/200)

Medipack Sterilisation Reels

The Medipack Reels are suitable for use in Steam and EO Gas Sterilisation.

- One-side clear film for easy identification of contents and the other medical grade paper.
- Integrated sterilisation indicator on medical grade paper side.

Product Code	Description
MEPR002	STERILISATION REEL - MEDIPACK 5CM X 200M (ROLL)
MEPR003	STERILISATION REEL - MEDIPACK 7CM X 200M (ROLL)
MEPR004	STERILISATION REEL - MEDIPACK 10CM X 200M (ROLL)
MEPR005	STERILISATION REEL - MEDIPACK 15CM X 200M (ROLL)

Medipack Sterilisation Indicator Tapes

The Medipack Sterilisation Indicator Tapes are designed for use in Steam Sterilisers.

Product Code	Description
ME320	TAPE - STERILISATION INDICATOR 12MM X 55M (ROLL)
ME321	TAPE - STERILISATION INDICATOR 19MM X 55M (ROLL)
ME322	TAPE - STERILISATION INDICATOR 25MM X 55M (ROLL)



Titems Monitor HP Chemical Indicator - Steam

The Titems monitor HP are designed for use in steam sterilisers providing an integrated response to defined combinations of time, temperature and presence of moist saturated steam.

- These attributes are monitored to a point that is well in excess of biological spore kill where upon the indicator spot will turn black.
- Class 4 performance with easy interpretation.

Product Code	Description
AUS415	STERILISATION INDICATORS - TITEMS MONITOR HP SHORT STRIPS (PKT/125)



Titems Diagnostic Emulator Chemical Indicator - Steam

The Titems diagnostic emulators are designed for use in steam sterilisers providing an integrated response to various defined combinations of time, temperature, pressure and presence of steam.

- These indicators provide an easily identifiable visual confirmation of correct sterilisation conditions within the steriliser chamber, in addition to aiding in the detection of a variety of steriliser malfunctions and human factor occurrences.
- Class 6 performance with easy interpretation.
- Permanent colour change.

Product Code	Description
AUS508	STERILISATION INDICATORS - TITEMS DIAG. EMULATOR (B/250)



Instrument Cleaning Brushes for CSSD

Product Code	Description
EDCB20043	BRUSH - INSTRUMENT NYLON (PKT/3)
EDCB20060	BRUSH - BOTTLE/TUBE 5CM X 30CM (PKT/3)
LR2020	CLEANING BRUSH 21.5CM - BL 6 X BD 1.3CM (EACH)
LR2022	CLEANING BRUSH 21.5CM - BL 6.5 X BD 1.5CM (EACH)
LR2024	CLEANING BRUSH 21.5CM - BL 9.5 X BD 2CM (EACH)
LR2030	CLEANING BRUSH 26.5CM - BL 6.5 X BD 0.3CM (EACH)
LR2032	CLEANING BRUSH 43CM - BL 6.5 X BD 0.3CM (EACH)
LR2034	CLEANING BRUSH 28.5CM - BL 7.5 X BD 0.7CM (EACH)
LR2036	CLEANING BRUSH 26.5CM - BL 6.5 X BD 0.8CM (EACH)



EDCB20060



EDCB20043