

CSSD Equipment

Includes:-

Static benches, Electric Height Adjustable Tables, Ultrasonic Cleaning Table, Heat Sealing Station, Electric Height Adjustable Wrap Stations, Endoscopy Case Cart, CSSD Case Cart, CSSD Transport Trolley, Customisable Options for CSSD Equipment



DB Health also has a wide range of products of interest such as Reusable Warwick Instrument Trays , Heine Diagnostic equipment including Endoscopy kits and Disinfectants products from Callington Haven. Contact DB Health to find out more or to request our latest catalogues



Web dbhealth.com.au

Email customerservice@denyer.com.au

Freecall: 1800 451 446

ABN: 37 000 459 174

Suppliers of quality products to the Healthcare Industry Est 1889.
Australian Owned & Operated.

Static Benches

Static benches are made of 304 Grade Welded Stainless Steel
Benches come in a variety of lengths from 75cm to 2.25metres in increments of 25cm.
***When ordering please state the length required.
Standard size 75cm (d) x 90 (h) excluding splash back
15cm Splashback
2 shelves. Lower shelf is 30cm from the ground
Dry edge
Height adjustable feet
Customisable



Illustrated is HQSP547-1250

| | |
|--|-------------|
| STATIC DRY BENCH W/ SPLASHBACK - EMERY 70X90X*** | HQSP547-*** |
|--|-------------|

Height Adjustable table

240V Height Adjustable Stainless Steel Table
150cm (W) x 75cm (D) x 80-145cm (H)
2 x lower Stainless Steel Wire Mesh Shelf
2 x 30cm deep Top Solid Stainless Steel shelves 30cm and 60cm High
4 x 12.5cm casters, 2 locking and 2 swivel
Customisable



| | |
|---|---------|
| HEIGHT ADJUSTABLE TABLE 240V - EMERY S/S 150X75X80-145CM (EACH) | HQSP661 |
|---|---------|

Ultrasonic Cleaning Table

Ultrasonic Cleaning Bench
100cm (W) x 60cm (D) x 50cm (H)
Height adjustable feet
1 shelf
Shelf braced for 40kg capacity
Customisable



| | |
|--|-----------|
| ULTRASONIC CLENAING TABLE - EMERY S/S 100X60X50CM (EACH) | HQSS22USC |
|--|-----------|

Heat Sealing Station

Heat Sealing Station has a Compact Laminate Working Surface
2 x Compact Laminate Storage Shelf (1200 x 300).
140cm (W) x 80cm (D) x 80-145cm (H).
Full Extension Ball Bearing drawer runners.
240 volt height adjustable.
140cm (W) x 80cm (D) x 80-145cm (H).
4 x casters, 2 locking, 2 swivel.
Over Bench Light, located under top storage shelf
Standard colour is white.
Customisable



| | |
|--|-----------|
| HEAT SEALING STATION 240V HEIGHT ADJUSTABLE - EMERY S/S W/ COM. LAM. | HQSP534.1 |
|--|-----------|

Electric Wrap Station

Electric Wrap Station Height Adjustable
201.5cm (L) x 75cm (D) x 90-130cm (H)
2 lower Stainless Steel Mesh Shelves (46cm x 122cm)
1 x 12.5cm casters, 2 locking and 2 swivel
Maxi Bin Rail and 2 x Steri Peel Holders
Push handles on both ends.
Handles rounded and entirely smooth to minimise infection control and make grasping easier
Lower Shelf 30cm from the floor
Customisable



| |
|--|
| WRAP STATION HEIGHT ADJUSTABLE 240V ELECTRIC - EMERY S/S 201.5X75X90-130CM (EACH) |
|--|

| |
|-----------|
| HQSP536-1 |
|-----------|

Endoscopy Case cart

Endoscopy Case Cart 73.5cm (L) x 48.8cm (D) x 109.1cm (H)
1 full size, magnetic door retainer with Stainless Steel Handle
8 x Pull Out Solid Stainless Trays (665x410x30), Shelves removable and height adjustable
1x Rounded Push Handle (LHS of case cart)
12.5cm Synthetic castors, 2 x Directional and 2 x Locking
4 x Buffers
Rated to 120 degree Celsius
Customisable



| |
|--|
| ENDOSCOPY CASE CART - EMERY 73.5X48.8X109.1CM (EACH) |
|--|

| |
|------------|
| HQSP398-16 |
|------------|

CSSD Case Cart

CSSD Case Cart 92cm (W) x 65cm (D) x 97cm (H)
1 x Fixed Wire Shelf, 1 top and bottom Shelf, Bottom Shelf 30cm from ground.
2 x Rounded Push Handles
2 x Door Handles
4 x 125mm Castors, 2 x Locking, 2 x Directional
4 x Corner buffers



| |
|-------------------------------------|
| CASE CART - EMERY 92X65X97CM (EACH) |
|-------------------------------------|

| |
|---------|
| HQSP398 |
|---------|

CSSD Transport Trolley

Transport Trolley – For Trolley Wash System 80cm (W) x 65cm (D) x 125cm (H) (+ push handles)
2 x Shelves + Top Angled
2 x side by side container holders
2 x Vertical Curved Push Handles on RHS
4 x 12.5cm Castors, 2 x Directional, 2 x Locking
Left and Right Side Panels
Front and Back Open Doors
Rated to 120 degree Celsius



| |
|---|
| CSSD TRANSPORT TROLLEY - EMERY 80X65X125CM (EACH) |
|---|


| |
|---------|
| HQSP397 |
|---------|

CSSD Customisable Options



HQRACK010

Sterilisation Rack - 590 x 180 x 110
Approximately 45mm gaps between uprights
Minimum order quantity 10 units



HQRACK020

Wash Rack – Container Lid
640 x 270 x 200
Holds 4 Lids



HQSP1002

Dolly Trolley Stacker Floor Mounted
Holds 8 units
48cm (W) x 44cm (D) x 160cm (H)



HQAO-BACK RACK

added to 500 x 500 Trolley
Mounts up to 3 tilt bin units
Legs to Extend Vertically
Tilt bins are 600 wide so will overhang outside edge of trolley by 50mm on each side
Tilt Bins and Baskets sold separately



HQSP440-3

Stainless Steel Wrap Rack - 160cm (W) x 60cm (D) x 120cm (H)
Open bottom, 4 hanging rails.
4 x 12.5cm Castors, 2 x Locking, 2 x swivel,



HQAO-S-MAG-BRKT

Add on Magnifying Lamp Bracket for Splashback (Hooks On)
Hooks simply over the top of CSSD Sink or Work Bench Splashback
Easy to fit, Move or Remove



HQAO-SHELF-UNDERSINKRACK

Stainless Steel Under Sink Tray Rack
Suits Autoplas Trays
1 Column of 5 trays
Under Sink Tray Rack
Mounted Underneath Sink Top or Bench Top




HQAO-WRAP RACK

Three Bar Wrap Rack Extension
Fitted to Height Adjustable Tables During Manufacture
Dimension: 430 x 700 x 430H



HQAO-STERIPEEL HOOK

100 mm Steripeel Hook
Fitted to under side of shelves
Stainless Steel/Welded



Transport trolley shelf divider

HQSP397-3 DIVIDER-S

200 x 65

HQSP397-3 DIVIDER-L

650 x 65

Helps prevent tray movement on mesh shelves



HQAO-Shelf 4

Wrap Top Mesh SS Shelf
1500 x 300 Above Bench Mesh Shelf
400mm from working surface
If Unit Modified Fit Shelf to Working Surface Length



HQAO-STERIREEL RAIL

Steripeel Rail (Removable)
Fitted to Back of Wrap station

HQAO-STERIPEEL HOOK02

50 mm Steripeel Hook
Fitted to under side of shelves
Stainless Steel/Welded



HQAO-Roller 1

Product: Transfer Roller – SS (400mm)
Dimensions: 50mm diameter, 400mm long
Other: Fitted to Table During Manufacture
Construction: Stainless Steel/Welded



HQMTCOV

Thermal Trolley Covers
Made to measure



HQAO-DDS-500-400x1

Item: One Drop Down Shelf – added to trolley in construction
Right Hand Side of Trolley
Dimensions: 400w x 500d x 30h
Shelves: 1 x Shelves
Construction: Stainless Steel/Welded

HQAO-DDS-500-400x2

Item: Two Drop Down Shelves – added to trolley in construction
Dimensions: 400w x 500d x 30h
Shelves: 2 x Shelves (one each side of the trolley)
Construction: Stainless Steel/Welded

HQAO-Roller 2-600

Product: Transfer Roller – PP
Dimensions: 50mm diameter, 600mm long
Other: Fitted to Table During Manufacture
Construction: Polypropylene