



est. 1889

Hospital Catalogue





About Us

Denyer Bros. Pty. Ltd was originally founded in 1889 by two brothers to manufacture surgical instruments. The business expanded from its original base in Melbourne with the Sydney operation becoming independently owned.

In 2000, Denyer Brothers re-named itself with the new trading name of DB Health. This change was to market & characterise the now broader scope of business activities and reflect the more current evolution in today's environment, compared to its original beginnings in 1889.

Today, DB Health remains a pivotal member in the health care community, with its experience and supply solutions strongly focused around patient safety. Quality products including Monitoring Electrodes, Bettershield Faceshield, Electric Beds for Aged Care, gelFast Hand Sanitiser, Patient Identification Bands, Silent Knight Pill Crushers and Skin Markers reflect just a sample of the varied product range available.

DB Health is still one of the few Australian medical wholesalers to be 100% family owned and operated.

*Suppliers of quality products to the
healthcare industry.*

Table of Content

1	■	DB Health	2
2	■	Safety	4
	■	Patient Identification	5
	■	Labels	13
	■	Markers	14
	■	Pill Crusher	15
3	■	Infection Control	16
	■	GelFast Anywear	17
	■	BetterShield	17
4	■	Equipment and Instruments	18
	■	Diagnostic	19
	■	Trolleys	22
	■	Electrodes	25
	■	Stethoscopes	31
5	■	Other	32
	■	Cadaver Bags	33
	■	Warwick Ware	33

Safety



Patient Identification

Barcode Identification

Sentry® LabelBand™

Designed to provide a waterproof seal to prevent damage to identification labels. Non-Transferable, durable, strong and comfortable to wear.

Product Code	Material	Closure Type	Band Size	Label Size	Bands per Box
PD5080-11	Smooth, soft & non-abrasive	SecurSnap	Adult 27.9 (l) x 3.2cm (w)	2.5cm x 7.6cm	500
PD5090-11	Smooth, soft & non-abrasive	SecurSnap	Adult/Pediatric 27.3cm x 3.1cm	2.5cm x 5cm	500



NOTE: Also available in alternative colours



PD5080-11

DuraPrint®

Product Code	Material	Closure Type	Label Size	Labels per Sheet	Bands per Box
PD9425-00	Smooth, soft with contoured edges	Adhesive tamper-evident	5.7cm x 2.5cm	20	1000
PD9425-11	Smooth, soft with contoured edges	Adhesive tamper-evident	5.7cm x 2.5cm	20	1000

Utilises a patent design that protects patient data and barcodes. Simple and quick to use, without the worry of damaging or transferring information.

NOTE: Available with multi colour tags or all red as pictured right.



PD9425-11

Short Stay LabelBand®

Constructed with water-resistant Tyvek® this is a durable, strong, yet lightweight and comfortable band. A clear adhesive shield also protects the barcode providing accurate scanning.

Product Code	Material	Closure Type	Band Size	Label Size	Bands per Box
PD3080-**	Smooth, soft & non-abrasive	Adhesive tamper-evident	Adult 25.4 (l) x 2.5cm (w)	7.3cm x 2.2cm	500

11

16

NOTE: Also available in alternative colours.
MUST BE ORDERED IN.



PD3080-16

Insert Bands

Veri-Colour®

This band is Australia's most popular. It is made from soft vinyl for comfort and utilises a plastic sleeve design for insert cards. Cards included.

Product Code	Material	Closure Type	Band Size	Insert Size	Bands per Box
PD140-**	Soft vinyl	Adhesive	Adult 33cm (l)	8.8cm x 2.2cm	250
PD120-**	Soft vinyl	Adhesive	Adult/Pediatric 33cm (l)	5.2cm x 2cm	250
PD110-**	Soft vinyl	Adhesive	Infant 16.5cm (l)	5.2cm x 1.4cm	250

11

16

NOTE: Also available in clear (10) and alternative colours
Infant colours available: 10, 11, 12 & 18.



PD140-16

Safeguard® Sealident®

Made for water resistance to protect patient information. Insert cards get placed into a plastic sleeve casing. Cards included.

Product Code	Material	Closure Type	Band Size	Label Size	Bands per Box
PD841-**	Tri-laminate	Adhesive	Adult 33cm (l)	8.2cm x 1.5cm	250
PD821-**	Tri-laminate	Adhesive	Adult/Pediatric 33cm (l)	5cm x 1.5cm	250
PD811-11	Tri-laminate	Adhesive	Infant 21cm (l)	5cm x 0.7cm	250

11

16

NOTE: All but PD811-11 are available in alternative colours



PD841-16

SoftGuard®

Soft tube design will not scratch delicate skin. Includes a clear sleeve design to protect patient information. Insert cards included.

Product Code	Material	Closure Type	Band Size	Label Size	Bands per Box
PD190-10	Soft tube design	Adhesive	Adult 25.4cm (l)	16.5cm x 0.7cm	50

NOTE: Available in clear (10) **ONLY**



Mother-Infant


Veri-Colour® Mother/Child

Ideal for maternity use, coming in 4, 3 or 2 parts. Made for comfort.

Product Code	Material	Closure Type	Band Size	Insert Size	Bands per Box
4 Part PD324-10	Soft vinyl	SecurSnap	Adult 23.5cm (l) Infant 16.5cm (l)	6cm x 1.9cm(Adult) 4.7cm x 1.7cm(Infant)	100 sets of 4
3 Part PD320-**	Soft vinyl	SecurSnap	Adult 23.5cm (l) Infant 16.5cm (l)	6cm x 1.9cm(Adult) 4.7cm x 1.7cm(Infant)	100 sets of 3
2 Part PD220-**	Soft vinyl	SecurSnap	Adult 23.5cm (l) Infant 16.5cm (l)	6cm x 1.9cm(Adult) 4.7cm x 1.7cm(Infant)	150 sets of 2

11 12 18

NOTE: Also available in clear




PD220-12

Speedi-Print®

Ideal for maternity use. Comes in 2 parts. Made for comfort and larger printing. Can be used with embossing machines or handwriting.

Product Code	Material	Closure Type	Band Size	Insert Size	Bands per Box
2 Part PD220S-**	Soft vinyl	SecurSnap	Adult 23.5cm (l) Infant 16.5cm (l)	6cm x 1.9cm Adult 4.7cm x 1.7cm Infant	150 sets of 2

11 12 18



Write-On & Imprint

Sentry® SuperBand®

The 5040 as shown is water-resistant and imprinted on, whilst the 5070 is written on.

Product Code	Material	Closure Type	Band Size	Write-on/Imprint Area	Bands per Box
PD5040-**	Lightweight non-stretch polyester	SecurSnap	Adult 29.2cm (l)	9cm x 2.5cm	500
PD5070-**	Lightweight non-stretch polyester	SecurSnap	Adult 25.4cm (l)	6.6cm x 2.5cm	500

11 16

NOTE: PD5070 is not available in 11, however both bands are available in alternative colours



PD5040-16

Speedi-Print®

Special Clearimage surface protects handwritten and embossed information.

Product Code	Material	Closure Type	Band Size	Imprint Size	Bands per Box
PD510-**	Soft vinyl	SecurSnap	Infant 16.5cm (l)	5.7cm x 1.9cm	250

11 12 18



Short Stay

Short Stay® Band

Ideal for day-only admissions and procedures.

Product Code	Material	Closure Type	Band Size	Label Size	Bands per Box
PD3000-**	Tyvek	Adhesive tamper-evident	Adult 25.4cm (l)	10.1cm x 2.5cm	1000



NOTE: Also available in alternative colours



Short Stay® Visitor Band

Coming in a variety of bright bold colours, this band is ideal for knowing your visitors. Comes pre imprinted with ability to write on.

Product Code	Material	Closure Type	Band Size	Label Size	Bands per Box
PD3015-**	Tyvek	Adhesive tamper-evident	Adult 25.4cm (l)	25.4cm x 2.5cm	1000



Alert Bands

SuperBand® Wide

Extra large text with the use of primary and secondary colours provides clarification to clinicians.

Product Code	Material	Closure Type	Band Size	Label Print	Bands per Box	Colour
PDP470-37R	Lightweight non-stretch polyester	SecurSnap	Adult 26.3cm (l)	Restricted Limb	500	Day Glow Pink
PDP5040-14CR	Lightweight non-stretch polyester	SecurSnap	Adult 29.2cm (l)	Caution Radiation	500	Yellow
PDP5055-12	Lightweight non-stretch polyester	SecurSnap	Adult 29.2cm (l)	Limb Alert	500	Pink
PDP5070-16A	Lightweight non-stretch polyester	SecurSnap	Adult 25.4cm (l)	Allergy	500	Red
PDP470-35	Lightweight non-stretch polyester	SecurSnap	Adult 26.3cm (l)	Oxygen Alert	500	Aqua

NOTE: Colours vary according to band label



PDP470-37R

Neonatal

Precision® Neonatal Wristbands

Available in two sizes, this moveable band offers newborns and preemies the highest quality, convenience and comfort.

Product Code	Material	Closure Type	Band Size	Label Size	Bands per Box
PD534-11	Soft comfortable foam	Hook and loop Velcro	15.8cm (l)	6.6cm x 2.5cm	12
PD524-11	Soft comfortable foam	Hook and loop Velcro	12.7cm (l)	6.6cm x 2.5cm	12



Other Bands

Xtender Band™

Available in two sizes, this adjustable band offers newborns and preemies the highest quality, convenience and comfort. Can be used with most ID bands for additional length.

Product Code	Material	Closure Type	Band Size	Label Size	Bands per Box
PD5062-11	Durable water resistant	Secursnap	Adult 25.4cm (I)	N/A	250

NOTE: Available in White **ONLY**



Smart Band® RFID System

Radio Frequency Identification- transfer digital information via radio waves between a tag and a reader.

Product Code	Closure Type	Band Width	Core Size	Outer Ø	Bands per Box
PD7070E-11	Adhesive tamper-evident	3cm	3.8cm	120mm	125

NOTE: Available in White **ONLY**



Foot Printer

Kleen-print® Disposable Foot Printer

With a clean and easy to print system, ink does not come into contact with baby, mother or nurse.

Product Code	Sheets per Box
PD3100	25



Newborn Identification Forms

Product Code	Forms per Pad	Sheets per Box
PD2010	50	25

If you would like to see accessories to go with our wristbands please contact your local sales person or refer to our website: www.dbhealth.com.au.

Labels

Hospital Labels

Product Code	Label	Labels per Roll
TMSY7983	Syringe Medication	333
TMFR7329	Medication added to I.V	400
TMFR1431	Allergic to	420
TMF04885	This Tray is Incomplete	333
TMFR7401	Medication added	333



Pediatric Labels

Product Code	Label	Labels per Roll
SLS004	I'm A Great Kid	80
SLS060	Star Patient	100
SLS024	I Was Very Brave Today	100
SLS008	Beware I'm Under Repair	80
SLS007	Star Patient	50
SLS003	I've Been Brave Today	100
SLS005	Purr-fect Patient	100
SLS061	Perfect Patient	100



If you would like to see our full range of both labels and stickers please contact your local sales person or refer to our website: www.dbhealth.com.au.

Markers

Securline®

Product Code	Product type	Markers per Box
PD1000	Surgical Skin Marker	10
PD1400-**	Lab Marker	10
PD1451	Superfrost-/ Marker II	10



Surgical Skin Marker provides fine and thick tip for marking patients skin.



Lab Marker is permanent and writes even on damp labware. It is ideal for writing on plastic, metal, glass and porcelain labware.



Superfrost®/ Marker II is ultra permanent and is ideal for histology and cytology. May be erased if Alcohol and Xylene are combined.

Pill Crusher

Silent Knight Pill Crusher

Product Code	Product type	Quantity
NS9999	Silent Knight Pill Crusher	1
NS9997	Pill Crusher Deluxe Pouches	1000

The Silent Knight is the quiet, effortless and quick way to crush tablets to a fine powder. Tablets are crushed directly into a pouch, making it clean and reduces the risk of cross contamination.



Infection Control



gelFast Anywear

gelFast Anywear is worn on garments (pant waists, belts, pockets, lanyards etc.) and used by squeezing the sanitising gel into the palm.



Product Code	Product type	Quantity
MX1101	Black Starter Unit	100
MX1102	White Starter Unit	100
MX2100	Gel Cartridge Refill	100

BetterShield

BetterShield provides full-face coverage and is both adjustable to fit over glasses and other safety wear. Does not restrict movement, nor will it fog over.



Product Code	Product type	Quantity per Box
MXFS01	BetterShield Disposable Face Shield	100

If you would like to see our full range of infection control products please contact your local sales person or refer to our website: www.dbhealth.com.au.

Equipment and Instruments



Diagnostic Laryngoscopes

Heine Classic+® Macintosh Fiber Optic Blades

Utilising up to 6500 individual micro-fibers for improved light transmission and longer life.

Compatible with all handles to the ISO 7376 (green standard).

Product Code	Size	Overall Length	Distal Width
HF0022100	Mac 0	82mm	9mm
HF0022101	Mac 1	93mm	11mm
HF0022102	Mac 2	115mm	13mm
HF0022103	Mac 3	135mm	15mm
HF0022104	Mac 4	155mm	15mm
HF0022105	Mac 5	176mm	15mm



Heine XP® Disposable Laryngoscope Blades

With proper use patient to patient cross-infection can be avoided.

Utilise the XHL Xenon Halogen bulb to achieve clarity.

Product Code	Size	Overall Length	Distal Width
HF0022761	XP Mac 0	90mm	9mm
HF0022762	XP Mac 1	114mm	10mm
HF0022763	XP Mac 2	134mm	14mm
HF0022764	XP Mac 4	154mm	14mm



If you would like to see our full range of Laryngoscopes and other diagnostic equipment please contact your local sales person or refer to our website: www.dbhealth.com.au.

Heine F.O. Laryngoscope Handles

Offering fiber optic handles for more light, Heine F.O. laryngoscope handles with Heine XHL® Xenon Halogen bulb technology gives up to 40% more light.

Product Code	Handle Type	Voltage	Batteries
HF0122860	Standard Handle	2.5	2 x C size- not included
HF0722885	Standard Handle	3.5	Liion Rechargeable battery- included
HF0822860	LED Standard Handle	2.5	2 x AA size- not included
HF0822891	LED Standard Handle	3.5	Liion Rechargeable battery- included
HF0122800	Small Handle	2.5	2 x AA size- not included
HF0122812	Short Handle	2.5	2 x AA size- not included
HF0122941	Angled Standard Handle	2.5	2 x C size- not included



HF0122812
(batteries not included)

Heine F.O. Laryngoscope Handles ACCESSORIES

Bulbs, chargers and rechargeable batteries that are sold separately.

Product Code	Product Type	Note	Per Box
HX0799383	Liion Rechargeable Battery		1
HX0188035	2.5v Bulb	Not for use with LED handles	1
HX0288044	3.5v Bulb	Not for use with LED handles	1
HX0299495	NT300 Recharger		1



HX0188035

Trolleys

Maternity Trolleys

Product Code	Product type	Dimensions
OX20-3812	Single Bassinet Trolley	80 x 42 x 106.5cm (frame)
OX20-3810	Premium Mobile Cot	145 x 82 x 145cm
OX20-3840	Twin Bassinet Trolley	80 x 64 x 106.5cm



OX20-3840

Medication Trolleys

Product Code	Product Type	Draws	Dimensions
OX18-4560	Stainless Steel 30 Dosage Medication Trolley	3	60 x 50 x 100cm
OX18-4660	Stainless Steel 50 Dosage Medication Trolley	5	60 x 50 x 100cm
OX55-2040	20-40 Unit Medication Trolley	6	75 x 50 x 100cm

NOTE: Also available in powder coat.
Other sizes available.



OX55-2040

Chart Trolleys

Product Code	Product Type	Dimensions
OX10-1200	24 Chart Trolley	75.2 x 50 x 90cm

NOTE: Also available in sizes 12 and 36.



Chart Holder

Product Code	Product Type	Dimensions
OX11-3036	Bulldog Clip	29 x 23cm
OX11-1000	2 Ring Clip	29 x 23cm

A 3mm thick white polycarbonate board, with curved top- designed to positioned on bed ends.

Dressing Trolley

Product Code	Rails	Draws	Dimensions
OX18-0450	3	1	50 x 50 x 90cm
OX18-0650	3	2	50 x 50 x 90cm

Available in stainless steel or powder coat. With widths in 10cm increments from 50cm to 120cm.



Instrument Trolley

Product Code	Rails	Dimensions
OX18-0250	3	50 x 50 x 90cm

Available in stainless steel or powder coat. With widths in 10cm increments from 50cm to 120cm.



Cylinder Storage Trolley

Product Code	Cylinder Size	Capacity	Dimensions
OX30-0100	C	6	47 x 30 x 75cm or 47 x 30 x 91cm

Other models include: Stationery 12 unit, Mobile 6 unit and Mobile 12 Unit.



Waste Bag Holder

Product Code	Dimensions	Bag Holder Size
OX10-1000	28 x 41 x 81cm	35 x 23.5cm

Yellow Bag for WBH: NS9000Y

Operation Room Trolleys

Product Code	Product Type	Draws	Dimensions
OX18-4410*	Urology Trolley	5	105 x 55 x 112cm
OX18-4300*	General Purpose Equipment Trolley	6	105 x 50 x 91cm
OX10-3612	Mayo Instrument Table Retractable top-Large on 3 wheels	N/A	60 x 95 x 87 to 130cm high



Also available: our range of anaesthetic, plaster, endoscopy, resuscitation, CSSD tables. Mayo tables are available as large or small with 3 wheels (rotatable top) or 4 wheels (non-rotatable).

*Available in Powder coat and Stainless Steel

Emergency Cart

Product Code	Draws	Dimensions
OX55-4471	5	80 x 50 x 94cm

With lift up side shelf, central locking, IV pole, O₂ C-size bracket, backboard and optional rotating tray for defib.



Linen Trolleys

Product Code	Product Type	Dimensions
OX26-1645	Single Soiled Linen Trolley w/ foot operated plastic lid	39 x 59.5 x 87cm
OX26-1655	Double Soiled Linen Trolley w/ foot operated plastic lid	76.5 x 59.5 X 87cm
OX26-1680	Triple Soiled Linen Trolley w/ foot operated plastic lid	113 x 59.5 x 87cm



Note: choice of lid colours available.

NOTE: A range of trolleys are available in Stainless Steel or Powder Coat.

If you would like to see our full range of trolleys, including a more comprehensive list of linen trolleys, IV stands and more, please contact your local sales person or refer to our website: www.dbhealth.com.au.

Electrodes

Disposable pre-gelled Electrodes

Foam ECG electrodes with solid gel and SSC sensor-clip connection.

Product Code	Product type	Size	Quantity per PKT	Order Quantity
F9060	Solid gel Adult	50 x 48mm	50	Carton of 1200
F9060P	Solid gel Pediatric	40mm round	50	Carton of 1200
F2069RM	Solid gel Adult Foam	42 x 55mm	50	Carton of 1200
F2069PRM	Solid gel Pediatric Foam	32 x 36mm	50	Carton of 1200

Foam ECG oval electrodes with solid gel and SSC sensor-clip connection.

Product Code	Product type	Size	Quantity per PKT	Order Quantity
F9079P	Solid gel Pediatric	32 x 36mm	50	Carton of 1250
F9079	Solid gel Adult	43 x 45mm	50	Carton of 1250

Non-woven ECG electrodes with solid gel and SSC sensor-clip connection.

Product Code	Product type	Size	Quantity per PKT	Order Quantity
F2060	Solid gel Adult	55mm round	50	Carton of 1200
F2060P	Solid gel Pediatric	40mm round	50	Carton of 1200

Non-woven ECG oval electrodes with solid gel and SSC sensor-clip connection.

Product Code	Product type	Size	Quantity per PKT	Order Quantity
F2079	Solid gel Adult	43 x 45mm	50	Carton of 1250
F2079P	Solid gel Pediatric	32 x 36mm	50	Carton of 1250

Foam small oval electrodes with solid gel and SSC sensor clip connection.

Product Code	Product type	Size	Quantity per PKT	Order Quantity
PG10S	Solid gel Foam	26mm	30	Carton of 600

Pre-gelled tab electrodes with solid gel, intended for short term monitoring and can be connected with a suitable alligator clip.

Product Code	Product type	Size	Quantity per PKT	Order Quantity
F3001ECG	TAB Electrode	34 x 22mm	100	Carton of 500

Defibrillation Pads

Direct connections to Philips Medical®, Laerdal Medical® and Agilent Philips Medical®.

Product Code	Product type	Quantity per PKT	Order Quantity
F7750	Multifunctional Radiotransparent	2	Carton of 20
F7950P	Multifunctional Defib Pads	2	Carton of 20
F7950W	Multifunctional with Lead Out	2	Carton of 20

Direct connection to Zoll Medical Corp®.

Product Code	Product type	Quantity per PKT	Order Quantity
F7751	Multifunctional Radiotransparent	2	Carton of 20
F7951P	Multifunctional Defib Pads	2	Carton of 20
F7951W	Multifunctional with Lead Out	2	Carton of 20

Direct connections to Medtronic-Physiocontrol®, Osatu Bexen and Cardioline®.

Product Code	Product type	Quantity per PKT	Order Quantity
F7752	Multifunctional Radiotransparent	2	Carton of 20
F7952P	Multifunctional Defib Pads	2	Carton of 20
F7952W	Multifunctional with Lead Out	2	Carton of 20

Direct connections to Agilent, Laerdal Medical® and Hewlett-Packard.

Product Code	Product type	Quantity per PKT	Order Quantity
F7753	Multifunctional Radiotransparent	2	Carton of 20
F7953P	Multifunctional Defib Pads	2	Carton of 20

Direct connections to Schiller and Esaote.

Product Code	Product type	Quantity per PKT	Order Quantity
F7754	Multifunctional Radiotransparent	2	Carton of 20
F7954P	Multifunctional Defib Pads	2	Carton of 20

Direct connections to Nihon Kohden, Welch Allyn® –Innomed and Fiab (transthoraci cardiac stimulator).

Product Code	Product type	Quantity per PKT	Order Quantity
F7755	Multifunctional Radiotransparent	2	Carton of 20
F7955P	Multifunctional Defib Pads	2	Carton of 20

Direct connection to Nihon Kohden.

Product Code	Product type	Quantity per PKT	Order Quantity
F7761	Multifunctional Radiotransparent	2	Carton of 20
F7961P	Multifunctional Defib Pads	2	Carton of 20

Direct connections to GE® and Esaote.

Product Code	Product type	Quantity per PKT	Order Quantity
F7762	Multifunctional Radiotransparent	2	Carton of 20
F7962P	Multifunctional Defib Pads	2	Carton of 20

Direct connection to Schiller.

Product Code	Product type	Quantity per PKT	Order Quantity
F7956P	Multifunctional Defib Pads	2	Carton of 20

Direct connections to Cardiac Science and GE®.

Product Code	Product type	Quantity per PKT	Order Quantity
F7957W	Multifunctional with Lead Out	2	Carton of 20

Electrosurgery

Liquid gel for ultrasound and electrostimulation. Salt free and without stringing effect- these gels are easily removable and resist drying.

Product Code	Product type	Quantity per PKT	Order Quantity
G008	Gel Ultrasound 260ml Clear	1	Carton of 35
G0084	Gel Ultrasound 5L Clear	1	Carton of 4

Connection cables for disposable ground pads.

Product Code	Product type	Quantity per PKT	Order Quantity
F7903	Universal Valleylab type Connection	1	Carton of 1
F7902	Jack Mono Plug 6.3mm Connection	1	Carton of 1
F7902/4	Plug 4mm Connection	1	Carton of 1
F7902/24	2 Plugs 4mm Connection	1	Carton of 1
F7903/F	Valleylab Connection	1	Carton of 1

70mm Sterile disposable electro-surgical electrodes.

Product Code	Product type	Quantity per PKT	Order Quantity
F4044	Disposable Sterile Electrode-Ball	1	Carton of 24
F4046	Disposable Sterile Electrode-Blade	1	Carton of 24
F4048	Disposable Sterile Electrode Needle	1	Carton of 24

70mm Sterile disposable stainless steel- Titanium coated electro-surgical

Product Code	Product type	Quantity per PKT	Order Quantity
F4044T	Disposable Sterile Electrode-Ball Titanium	1	Carton of 24
F4046T	Disposable Sterile Electrode-Blade Titanium	1	Carton of 24
F4048T	Disposable Sterile Electrode Needle Titanium	1	Carton of 24

Sterile disposable curved loop electro-surgical electrodes.

Product Code	Product type	Quantity per PKT	Order Quantity
F4621	Disposable Sterile Electrode Round Loop	1	Carton of 24

Disposable grounding pads with connections.

Product Code	Product type	Quantity per PKT	Order Quantity
F7805W/V	Sterile Abrasive Tip Cleaner Pads with Volleytlab connection cable.	1	Carton of 100
F7805PW/V	Adult Unipolar with Volleytlab connection cable.	1	Carton of 25
F7805W/6.3	Adult Unipolar with 6.3mm bipolar connection.	1	Carton of 25
F7805PW/6.3	Pediatric Unipolar with 6.3mm bipolar connection.	1	Carton of 25

Disposable grounding pads without connections.

Product Code	Product type	Quantity per PKT	Order Quantity
F7805	Adult Unipolar without Cable	1	Carton of 25
F7805P	Pediatric Unipolar	1	Carton of 25
F7805N	Neonatal Unipolar	1	Carton of 25

Disposable electrosurgical pencils.

Product Code	Product type	Quantity per PKT	Order Quantity
F4798	Hand Switch Pencil	1	Carton of 25
F4798/H	Hand Switch Pencil with Holster	1	Carton of 25
F4797/L	Hand Switch Pencil with Long Blade	1	Carton of 50
F4545	Holster - Sterile	1	Carton of 25

Stethoscopes

Spirit® Deluxe Stethoscopes

Product Code	Product type	Head Size	Tubing	Order Quantity
STS601-**	Stainless Steel Dual Head Stethoscope- Adult	47mm Ø	Superflex™ Binaural	Box of 1
STS606-**	Stainless Steel Dual Head Stethoscope- Pediatric	38mm Ø	Superflex™ Binaural	Box of 1
STS607-**	Stainless Steel Dual Head Stethoscope- Infant	32mm Ø	Superflex™ Binaural	Box of 1



Spirit® Majestic Stethoscopes

Product Code	Product type	Head Size	Tubing	Order Quantity
STS603-**	Stethoscope- Adult	47mm Ø	Simplex™ Binaural	Box of 1



If you would like to see our full range of stethoscopes, including a more comprehensive list of colours, please contact your local sales person or refer to our website: www.dbhealth.com.au.

Other



Post Mortem Bags

Product Code	Product Type	Size	Bags per Box
PD55-7500	Post Mortem Bag Adult Curved Zipper	91.4 x 228.6cm	10



Warwick Ware

Lotion Bowls

Product Code	Diameter	Depth	Volume
WALB100	100mm	45mm	280ml
WALB150	150mm	70mm	900ml
WALB200	200mm	90mm	2L
WALB250	250mm	105mm	3L



Wash Bowls

Product Code	Diameter	Depth	Volume
WAWB3112	310mm	105mm	3L
WAWB3513	350mm	120mm	5L



Gallipots

Product Code	Diameter	Depth	Volume
WAGP40	40mm	30mm	35ml
WAGP60	60mm	40mm	94mL
WAGP80	80mm	60mm	280ml



Jugs

Product Code	Volume	Colour
WAJG500	500ml	Natural
WAJG1000	1000ml	Natural
DOJG1000BLUE	1000ml	Light Blue



Instrument Trays

Product Code	Length	Width	Depth
WAIT1813	181mm	134mm	54mm
WAIT1913	195mm	135mm	23mm
WAIT2015	200mm	150mm	30mm
WAIT2715	270mm	150mm	51mm
WAIT3025	300mm	250mm	52mm
WAIT4030	420mm	305mm	75mm



Perforated Instrument Trays

Product Code	Length	Width	Depth
WAPIT2015	200mm	150mm	51mm
WAPIT3025	300mm	250mm	52mm
WAPIT4030	420mm	305mm	75mm
DOTRPD30X20GRN	300mm	205mm	



Compartment Trays

Product Code	Compartments	Length	Width	Depth
WACIT1814	3	181mm	134mm	54mm
WACIT2013	3	205mm	130mm	30mm
WACIT2215*	4	270mm	180mm	41mm
WACIT2718	4	220mm	150mm	24mm



*Lightweight material

Kidney Dishes

Product Code	Length	Width	Depth
WAKD150	150mm	35mm	75mm
WAKD200	200mm	45mm	98mm
WAKD250	250mm	55mm	130mm
WAKD300	300mm	65mm	144mm
MU05-800*	240mm	50mm	120mm



*Disposable

Urinals

Product Code	Type	Size
WAUML1000	Male	1200ml
WAUF1000	Female	1000ml
WASPN	Bed Pan	335 x 245mm



Vomit Bowl

Product Code	Length	Width	Depth	Volume
WAVB300	300mm	144mm	65mm	1.5L





est. 1889

Contact Us

Unit 17
Metro Centre
38-46 South Steet
Rydalmerre NSW 2116

Phone: (02) 9638 2100
Freecall: 1800 451 446
Fax: (02) 9638 2700
Email: customerservice@denyer.com.au

PO BOX 305
Rydalmerre NSW 1701

www.dbhealth.com.au